INSTRUCTIONS FOR USE (IFU)

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MANUFACTURER

Manufacturer	Zima International, Inc. Also DBA Dandy
Address	1320 N. 300 W., Lehi, UT 84043, United States
Contact / Website	www.meetdandy.com
Complaints	complaints@meetdandy.com

DEVICE IDENTITY

Device Identity

Device Family	Sleep-Disordered Breathing Orthosis
GMDN / EMDN Code	47526
Product Code(s)	SLEEPAPD, SLEEPAPH, SLEEPAPPC, SLEEPAPX3
Device Type	Custom-made dental device (CMD) / Patient Matched Medical Device (PMMD)
Sterility	Not supplied sterile
Single-use / Reuse	Reusable
Intended User	Prescriber: Dental professional / dental practice

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AUSTRALIA (ENGLISH)

INSTRUCTIONS FOR USE (IFU)

Australia — Patient Matched Medical Device (PMMD)

GMDN / EMDN	47526
Product Code(s)	SLEEPAPD, SLEEPAPH, SLEEPAPPC, SLEEPAPX3
Manufacturer	Zima International, Inc. Also DBA Dandy 1320 N. 300 W., Lehi, UT 84043 United States
Australian Sponsor	Name: DANDY Labs AUS Pty Ltd Company No.: ABN: 33693955761 Address: Tower One - International Towers, Level 46, 100 Barangaroo Avenue, Sydney, New South Wales 2000, Australia
Intended User	Prescriber — Dental professional / dental practice
Device Type	Patient Matched Medical Device (PMMD)
Sterility	Not supplied sterile
Single-use / Reusable	Reusable

1) Intended Purpose (Intended Use)

A Dandy sleep apnea device is a type of dental device used to treat snoring and mild to moderate obstructive sleep apnea. It is also known as a mandibular advancement device or mandibular repositioning device.

2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: Dandy devices are custom-made to fit your patient's mouth and are designed to be as comfortable and unobtrusive as possible.
- Material: PMMA Block or Nylon

[PMMA Block](#)[Nylon](#)

3) Contraindications

- Sleep apnea devices may not be appropriate for individuals with severe respiratory conditions, such as severe chronic obstructive pulmonary disease (COPD), severe asthma, or other conditions that the use of positive airway pressure could exacerbate.
- Oral devices, which are used to treat sleep apnea by repositioning the jaw and tongue, may not be suitable for individuals with severe dental issues, such as untreated gum disease or extensive

tooth decay. Similarly, individuals with severe jaw problems or temporomandibular joint (TMJ) disorders may require careful evaluation to determine if an oral appliance is appropriate.

4) Warnings / Precautions / Potential Risks

WARNING

It is common for patients to experience some initial discomfort or soreness when the appliance is first placed or adjusted. This discomfort usually subsides as the patient adapts to the device.

CAUTION

N/A

5) Cleaning and Care

Directions for Use

To ensure proper use of the devices, follow these steps:

- Before inserting the sleep apnea device, brush and floss the teeth and rinse the mouth with water to ensure good oral hygiene and prevent any food particles from getting trapped under the device.
- The appliance is crafted to cover the upper and lower teeth and stays secure through adaptable metallic connectors. You may engage the device by means of adjustable lugs or by the Herbst® mechanism, which is a rod and tube assembly that orients the jaws in a predetermined relationship. Place the upper piece of the appliance in the patient's mouth and have the patient bite down gently to position it properly.
- The Herbst® is a one-piece construction held together by two adjustment mechanisms on the buccal or outer area of the upper and lower appliances. The adjustment mechanisms are designed for the doctor to advance the mandible to the desired position incrementally. Dorsal is a two-piece design with separate upper and lower acrylic portions that, when engaged, propel the mandible into a protrusive position via acrylic fins built in the lower acrylic portion.
- The Dorsal is an adjustable design utilizing advancement; once the upper piece is in place, insert the lower piece of the device, and then adjust the metal connectors to hold the lower jaw in a slightly advanced position. Finding the most comfortable and effective position may require some trial and error.
- The patient should breathe through their nose and ensure that the appliance is properly seated and comfortable. It may take some time to adjust to the feeling of the appliance in the mouth, but most patients find that they can sleep comfortably with it in place.
- The patient is advised to use the device according to the guidelines provided by your practice, which may encompass detailed instructions regarding the daily duration of wear and the appropriate times for removal.

Cleaning instructions

- Dandy recommends cleaning your sleep apnea device daily before and after every use. Immediately after removal, use warm water (not to exceed 45°C/113°F) and soap to remove remaining saliva and debris from the device. The manufacturer advises against immersing the device in solutions, which could lead to delamination.

6) Storage

Use a case specifically designed for dental devices to protect them from contamination and damage. The case should have good airflow to prevent moisture buildup.

7) Expected Life / Service

Sleep apnea mouth appliances (oral appliances) typically last 2 to 5 years, but this duration can vary depending on the design and care. Signs indicating the need for replacement include a poor fit, wear, or reduced effectiveness. Therefore, regular dental check-ups are crucial to monitor their condition and ensure they remain functional.

8) Incident / Complaint Reporting

Report suspected serious incidents associated with this device to:

- The Manufacturer and/or Dandy Australia entity using the contacts above, and
- Therapeutic Goods Administration (TGA)

In Australia, serious incidents involving medical devices must be reported to the TGA in accordance with the Therapeutic Goods Act 1989.

INFO

For full contact details of the local Authorized Representative, refer to: Authorized Representatives — end of this document.

9) Disposal

To dispose of a sleep apnea oral appliance, you should primarily consult with your dental specialist or use local e-waste/recycling facilities. Do not dispose of it with regular household trash unless directed by local regulations.

EUROPEAN UNION (ENGLISH)

INSTRUCTIONS FOR USE (IFU)

European Union — Custom-Made Device (CMD)

GMDN / EMDN	Q030102 - ANTISNORING AND OBSTRUCTIVE SLEEP APNOEA DEVICES
Product Code(s)	SLEEPAPD, SLEEPAPH, SLEEPAPPC, SLEEPAPX3
Manufacturer	Zima International, Inc. Also DBA Dandy 1320 N. 300 W., Lehi, UT 84043 United States
EU authorized representative (EURP)	Name: Dandy Labs Europe Company No.: Registration No.: 989857065 Address: 3, boulevard de Sebastopol, 75001 Paris, France
Intended User	Prescriber — Dental professional / dental practice
Device Type	Custom-made dental device (CMD)
Sterility	Not supplied sterile
Single-use / Reusable	Reusable

1) Intended Purpose (Intended Use)

A Dandy sleep apnea device is a type of dental device used to treat snoring and mild to moderate obstructive sleep apnea. It is also known as a mandibular advancement device or mandibular repositioning device.

2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: Dandy devices are custom-made to fit your patient's mouth and are designed to be as comfortable and unobtrusive as possible.
- Material: PMMA Block or Nylon

[PMMA Block](#)[Nylon](#)

3) Contraindications

- Sleep apnea devices may not be appropriate for individuals with severe respiratory conditions, such as severe chronic obstructive pulmonary disease (COPD), severe asthma, or other conditions that the use of positive airway pressure could exacerbate.
- Oral devices, which are used to treat sleep apnea by repositioning the jaw and tongue, may not be suitable for individuals with severe dental issues, such as untreated gum disease or extensive

tooth decay. Similarly, individuals with severe jaw problems or temporomandibular joint (TMJ) disorders may require careful evaluation to determine if an oral appliance is appropriate.

4) Warnings / Precautions / Potential Risks

WARNING

It is common for patients to experience some initial discomfort or soreness when the appliance is first placed or adjusted. This discomfort usually subsides as the patient adapts to the device.

CAUTION

n/a

5) Cleaning and Care

Directions for Use

To ensure proper use of the devices, follow these steps:

- Before inserting the sleep apnea device, brush and floss the teeth and rinse the mouth with water to ensure good oral hygiene and prevent any food particles from getting trapped under the device.
- The appliance is crafted to cover the upper and lower teeth and stays secure through adaptable metallic connectors. You may engage the device by means of adjustable lugs or by the Herbst® mechanism, which is a rod and tube assembly that orients the jaws in a predetermined relationship. Place the upper piece of the appliance in the patient's mouth and have the patient bite down gently to position it properly.
- The Herbst® is a one-piece construction held together by two adjustment mechanisms on the buccal or outer area of the upper and lower appliances. The adjustment mechanisms are designed for the doctor to advance the mandible to the desired position incrementally. Dorsal is a two-piece design with separate upper and lower acrylic portions that, when engaged, propel the mandible into a protrusive position via acrylic fins built in the lower acrylic portion.
- The Dorsal is an adjustable design utilizing advancement; once the upper piece is in place, insert the lower piece of the device, and then adjust the metal connectors to hold the lower jaw in a slightly advanced position. Finding the most comfortable and effective position may require some trial and error.
- The patient should breathe through their nose and ensure that the appliance is properly seated and comfortable. It may take some time to adjust to the feeling of the appliance in the mouth, but most patients find that they can sleep comfortably with it in place.
- The patient is advised to use the device according to the guidelines provided by your practice, which may encompass detailed instructions regarding the daily duration of wear and the appropriate times for removal.

Cleaning instructions

- Dandy recommends cleaning your sleep apnea device daily before and after every use. Immediately after removal, use warm water (not to exceed 45°C/113°F) and soap to remove remaining saliva and debris from the device. The manufacturer advises against immersing the device in solutions, which could lead to delamination.

6) Storage

Use a case specifically designed for dental devices to protect them from contamination and damage. The case should have good airflow to prevent moisture buildup.

7) Expected Life / Service

Sleep apnea mouth appliances (oral appliances) typically last 2 to 5 years, but this duration can vary depending on the design and care. Signs indicating the need for replacement include a poor fit, wear, or reduced effectiveness. Therefore, regular dental check-ups are crucial to monitor their condition and ensure they remain functional.

8) Incident / Complaint Reporting

Report suspected serious incidents associated with this device to:

- The Manufacturer and/or Dandy EU entity using the contacts above, and corresponding local requirements
- The Competent Authority of your Member State, in accordance with Regulation (EU) 2017/745 (MDR).

Competent Authorities:

- France: Agence Nationale de Sécurité du Médicament et des Produits de Santé (ANSM)
- Germany : Bundesinstitut für Arzneimittel und Medizinprodukte (BfArM) (BfArM)
- Italy : Ministero della Salute / Istituto Superiore di Sanità (ISS)
- Spain: Agencia Española de Medicamentos y Productos Sanitarios (AEMPS)
- Netherlands: Inspectie Gezondheidszorg en Jeugd (IGJ)
- Austria: Bundesamt für Sicherheit im Gesundheitswesen (BASG)
- Sweden: Läkemedelsverket
- Denmark: Lægemiddelstyrelsen
- Portugal: Autoridade Nacional do Medicamento e Produtos de Saúde (INFARMED)

INFO

For full contact details of the local Authorized Representative, refer to: Authorized Representatives — end of this document.

9) Disposal

To dispose of a sleep apnea oral appliance, you should primarily consult with your dental specialist or use local e-waste/recycling facilities. Do not dispose of it with regular household trash unless directed by local regulations.

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UNITED KINGDOM (ENGLISH)

INSTRUCTIONS FOR USE (IFU)

United Kingdom — Custom-Made Device (CMD)

GMDN / EMDN	47526
Product Code(s)	SLEEPAPD, SLEEPAPH
Manufacturer	Zima International, Inc. Also DBA Dandy 1320 N. 300 W., Lehi, UT 84043 United States
UK Responsible Person (UKRP)	Name: DANDY LABS GB, LTD Company No.: Company Number: 16873608 Address: 5 New Street Square, London EC4A 3TW, United Kingdom
Intended User	Prescriber — Dental professional / dental practice
Device Type	Custom-made dental device (CMD)
Sterility	Not supplied sterile
Single-use / Reusable	Reusable

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- The Manufacturer and/or UK Responsible Person (UKRP) using the contacts above, and
- The UK Medicines and Healthcare products Regulatory Agency (MHRA) — www.gov.uk/mhra

In accordance with The Medical Devices Regulations 2002 (as amended).

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The table below lists the local Authorized Representative or Sponsor for each country / jurisdiction where this device is marketed. For incident reporting, contact the representative for your country.

La tabla siguiente enumera el Representante Autorizado o Patrocinador local para cada país. | Le tableau ci-dessous liste le Représentant Autorisé local pour chaque pays. | Die folgende Tabelle enthält die bevollmächtigten Vertreter je Land.

Jurisdiction	Authorised Rep / Sponsor	Registered Address / Contact Email	Reg. Number	Regulatory Body
Australia (AUS)	DANDY Labs AUS Pty Ltd	Tower One, Level 46, 100 Barangaroo Ave, Sydney NSW 2000 Compliance@meetdandy.com	ABN: 33693955761	TGA (Therapeutic Goods Administration)
EU — English	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, France Compliance@meetdandy.com	Reg: 989857065	Competent Authority / EUDAMED
EU — España	Dandy Labs Europa	3, boulevard de Sebastopol, 75001 París, Francia Compliance@meetdandy.com	Reg: 989857065	AEMPS / Autoridad Competente
EU — France	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, France Compliance@meetdandy.com	Reg: 989857065	ANSM / Autorité Compétente
EU — Deutschland	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, Frankreich Compliance@meetdandy.com	Reg: 989857065	BfArM / Zuständige Behörde
EU — Italia	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Parigi, Francia Compliance@meetdandy.com	Reg: 989857065	Min. della Salute / ISS
United Kingdom (UK)	DANDY LABS GB, LTD	5 New Street Square, London EC4A 3TW Compliance@meetdandy.com	Co. No.: 16873608	MHRA

NOTE

Contact details are subject to change. Always verify current information on the official product labelling or the manufacturer's website before reporting an incident. / Les coordonnées sont susceptibles de changer. / Angaben können sich ändern. / I recapiti sono soggetti a modifiche.

**Manufacturer / Hersteller /
Fabbricante / Fabricant /
Fabricante**

Zima International, Inc. Also DBA Dandy
1320 N. 300 W., Lehi, UT 84043, United States
compliance@meetdandy.com

This is a representation of an electronic record that was signed electronically and this page is the manifestation of the electronic signatures.

Signatory Table

Action Name	User Name	Title	Signature Date
Approve	Kasha Lim	Supplier Quality Engineer	18-Jun-2026 06:41

* Datetime fields are displayed according to the system time zone: (GMT-04:00) Eastern Daylight Time (America/New_York)