INSTRUCTIONS FOR USE (IFU)

Dandy / Zima International, Inc. | AUS · NZ · UK · EU/EN · EU/ES · EU/FR · EU/DE · EU/IT · EU/DA ·
EU/NL · EU/PL · EU/SV · EU/NB · EU/PT

MANUFACTURER

Manufacturer	Zima International, Inc. Also DBA Dandy
Address	1320 N. 300 W., Lehi, UT 84043, United States
Contact / Website	www.meetdandy.com
Complaints	complaints@meetdandy.com

DEVICE IDENTITY

Device Identity

Device Family	Denture Reliner - Hard
GMDN / EMDN Code	17609
Product Code(s)	ACRFDENTRLN, DENTFLRLN, MILLDENTRLN
Device Type	Custom-made dental device (CMD) / Patient Matched Medical Device (PMMD)
Sterility	Not supplied sterile
Single-use / Reuse	Reusable
Intended User	Prescriber: Dental professional / dental practice

TABLE OF CONTENT

AUSTRALIA (ENGLISH)	4
1) Intended Purpose (Intended Use)	4
2) Device Description / Key Specifications	4
3) Contraindications	4
4) Warnings / Precautions / Potential Risks	5
5) Cleaning and Care	5
6) Storage	5
7) Expected Life / Service	5
8) Incident / Complaint Reporting	6
9) Disposal	6
NEW ZEALAND (ENGLISH)	7
1) Intended Purpose (Intended Use)	7
2) Device Description / Key Specifications	7
3) Contraindications	7
4) Warnings / Precautions / Potential Risks	8
5) Cleaning and Care	8
6) Storage	8
7) Expected Life / Service	8
8) Incident / Complaint Reporting	8
9) Disposal	9
UNITED KINGDOM (ENGLISH)	10
1) Intended Purpose (Intended Use)	10
2) Device Description / Key Specifications	10
3) Contraindications	10
4) Warnings / Precautions / Potential Risks	11
5) Cleaning and Care	11
6) Storage	11
7) Expected Life / Service	11
8) Incident / Complaint Reporting	12
9) Disposal	12
EUROPEAN UNION (ENGLISH)	13
1) Intended Purpose (Intended Use)	13
2) Device Description / Key Specifications	13
3) Contraindications	13
4) Warnings / Precautions / Potential Risks	14
5) Cleaning and Care	14
6) Storage	14
7) Expected Life / Service	14
8) Incident / Complaint Reporting	14
9) Disposal	15
UNIÓN EUROPEA (ESPAÑOL)	16
1) Finalidad prevista (uso previsto)	16
2) Descripción del dispositivo / Especificaciones clave	16
3) Contraindicaciones	16

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AUSTRALIA (ENGLISH)

INSTRUCTIONS FOR USE (IFU)

Australia — Patient Matched Medical Device (PMMD)

GMDN / EMDN	17609
Product Code(s)	ACRFDENTRLN, DENTFLRLN, MILLDENTRLN
Manufacturer	Zima International, Inc. Also DBA Dandy 1320 N. 300 W., Lehi, UT 84043 United States
Australian Sponsor	Name: DANDY Labs AUS Pty Ltd Company No.: ABN: 33693955761 Address: Tower One - International Towers, Level 46, 100 Barangaroo Avenue, Sydney, New South Wales 2000, Australia
Intended User	Prescriber — Dental professional / dental practice
Device Type	Patient Matched Medical Device (PMMD)
Sterility	Not supplied sterile
Single-use / Reusable	Reusable

1) Intended Purpose (Intended Use)

This custom-made denture reliner is intended to improve denture fit/retention when a denture no longer fits as intended and requires relining, as prescribed by a licensed dental professional.

2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: Hard denture reline (full arch) (as applicable)
- Material: Sofreliner (Silicone Dioxide)

[Liner Material Safety Data Sheets \(Paste\)](#)[Liner Material Safety Data Sheets \(Primer\)](#)

3) Contraindications

Do not use in patients with:

- Known hypersensitivity to the device material
- Untreated oral lesions

4) Warnings / Precautions / Potential Risks

⚠ **WARNING**

Potential risks associated with denture lining materials may include:

- Gingival irritation, mucosal soreness, ulceration
- Poor adaptation may cause instability or trauma
- Allergic response
- Separation of material

If the patient experiences pain, irritation, or other adverse effects, discontinue use and contact the treating dental professional.

⚠ **CAUTION**

NA

5) Cleaning and Care

Daily cleaning:

- Brush gently using a soft-bristled denture brush or soft toothbrush.
- Avoid abrasive products (e.g., toothpaste, abrasive cleaners, strong soaps).
- Use mild, non-abrasive liquid dish soap as directed by the dental professional.
- Rinse thoroughly under cool or lukewarm running water.
- Avoid high heat; do not use hot or boiling water.

Soaking:

- Store/soak as directed by the dental professional to prevent drying.
- Avoid prolonged soaking in effervescent cleaners if this may damage the material.

Do not use:

- Avoid cleaning/immersing the device in sodium hypochlorite (chlorine bleach), hydrogen peroxide, mouthwash, or alcohol products, as these may damage or discolour the material.

6) Storage

When not in use, store the device as directed by the dental professional. Avoid excessive heat, pressure, or conditions that may deform or damage the device.

7) Expected Life / Service

With appropriate care and maintenance, the device is intended for ongoing use as directed by the dental professional. If the device becomes damaged or fit is no longer acceptable, the patient should contact their treating dental professional for evaluation.

8) Incident / Complaint Reporting

Report suspected serious incidents associated with this device to:

- The Manufacturer and/or Dandy Australia entity using the contacts above, and
- Therapeutic Goods Administration (TGA)

In Australia, serious incidents involving medical devices must be reported to the TGA in accordance with the Therapeutic Goods Act 1989.

INFO

For full contact details of the local Authorized Representative, refer to: Authorized Representatives — end of this document.

9) Disposal

Dispose of the device in accordance with local regulations and the dental practice's procedures. If the device is contaminated with bodily fluids, treat it as clinical/infectious waste per local requirements.

NEW ZEALAND (ENGLISH)

INSTRUCTIONS FOR USE (IFU)

NEW ZEALAND — Patient Matched Medical Device (PMMD)

GMDN / EMDN	17609
Product Code(s)	ACRFDENTRLN, DENTFLRLN, MILLDENTRLN
Manufacturer	Zima International, Inc. Also DBA Dandy 1320 N. 300 W., Lehi, UT 84043 United States
New Zealand Sponsor	Name: DANDY Labs NZ Limited Address: Quigg Partners, Level 7, 36 Brandon Street, Wellington, NZ 6011
Intended User	Prescriber — Dental professional / dental practice
Device Type	Patient Matched Medical Device (PMMD)
Sterility	Not supplied sterile
Single-use / Reusable	Reusable

1) Intended Purpose (Intended Use)

This custom-made denture reliner is intended to improve denture fit/retention when a denture no longer fits as intended and requires relining, as prescribed by a licensed dental professional.

2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: Hard denture reline (full arch) (as applicable)
- Material: Sofreliner (Silicone Dioxide)

[Liner Material Safety Data Sheets \(Paste\)](#)[Liner Material Safety Data Sheets \(Primer\)](#)

3) Contraindications

Do not use in patients with:

- Known hypersensitivity to the device material
- Untreated oral lesions

4) Warnings / Precautions / Potential Risks

⚠ **WARNING**

Potential risks associated with denture lining materials may include:

- Gingival irritation, mucosal soreness, ulceration
- Poor adaptation may cause instability or trauma
- Allergic response
- Separation of material

If the patient experiences pain, irritation, or other adverse effects, discontinue use and contact the treating dental professional.

⚠ **CAUTION**

NA

5) Cleaning and Care

Daily cleaning:

- Brush gently using a soft-bristled denture brush or soft toothbrush.
- Avoid abrasive products (e.g., toothpaste, abrasive cleaners, strong soaps).
- Use mild, non-abrasive liquid dish soap as directed by the dental professional.
- Rinse thoroughly under cool or lukewarm running water.
- Avoid high heat; do not use hot or boiling water.

Soaking:

- Store/soak as directed by the dental professional to prevent drying.
- Avoid prolonged soaking in effervescent cleaners if this may damage the material.

Do not use:

- Avoid cleaning/immersing the device in sodium hypochlorite (chlorine bleach), hydrogen peroxide, mouthwash, or alcohol products, as these may damage or discolour the material.

6) Storage

When not in use, store the device as directed by the dental professional. Avoid excessive heat, pressure, or conditions that may deform or damage the device.

7) Expected Life / Service

With appropriate care and maintenance, the device is intended for ongoing use as directed by the dental professional. If the device becomes damaged or fit is no longer acceptable, the patient should contact their treating dental professional for evaluation.

8) Incident / Complaint Reporting

Report suspected serious incidents associated with this device to:

- The Manufacturer and/or Dandy New Zealand entity using the contacts above, and
- Medsafe

In New Zealand, serious incidents and quality issues involving medical devices should be reported to Medsafe in accordance with New Zealand's medical device adverse event reporting requirements.

INFO

For full contact details of the local Authorized Representative, refer to: Authorized Representatives — end of this document.

9) Disposal

Dispose of the device in accordance with local regulations and the dental practice's procedures. If the device is contaminated with bodily fluids, treat it as clinical/infectious waste per local requirements.

UNITED KINGDOM (ENGLISH)

INSTRUCTIONS FOR USE (IFU)

United Kingdom — Custom-Made Device (CMD)

GMDN / EMDN	17609
Product Code(s)	ACRFDENTRLN, DENTFLRLN, MILLDENTRLN
Manufacturer	Zima International, Inc. Also DBA Dandy 1320 N. 300 W., Lehi, UT 84043 United States
UK Responsible Person (UKRP)	Name: DANDY LABS GB, LTD Company No.: Company Number: 16873608 Address: 5 New Street Square, London EC4A 3TW, United Kingdom
Intended User	Prescriber — Dental professional / dental practice
Device Type	Custom-made dental device (CMD)
Sterility	Not supplied sterile (if applicable)
Single-use / Reusable	Reusable

1) Intended Purpose (Intended Use)

This custom-made denture reliner is intended to improve denture fit/retention when a denture no longer fits as intended and requires relining, as prescribed by a licensed dental professional.

2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: Hard denture reline (full arch) (as applicable)
- Material: Sofreliner (Silicone Dioxide)

[Liner Material Safety Data Sheets \(Paste\)](#)[Liner Material Safety Data Sheets \(Primer\)](#)

3) Contraindications

Do not use in patients with:

- Known hypersensitivity to the device material
- Untreated oral lesions

4) Warnings / Precautions / Potential Risks

⚠ **WARNING**

Potential risks associated with denture lining materials may include:

- Gingival irritation, mucosal soreness, ulceration
- Poor adaptation may cause instability or trauma
- Allergic response
- Separation of material

If the patient experiences pain, irritation, or other adverse effects, discontinue use and contact the treating dental professional.

⚠ **CAUTION**

NA

5) Cleaning and Care

Daily cleaning:

- Brush gently using a soft-bristled denture brush or soft toothbrush.
- Avoid abrasive products (e.g., toothpaste, abrasive cleaners, strong soaps).
- Use mild, non-abrasive liquid dish soap as directed by the dental professional.
- Rinse thoroughly under cool or lukewarm running water.
- Avoid high heat; do not use hot or boiling water.

Soaking:

- Store/soak as directed by the dental professional to prevent drying.
- Avoid prolonged soaking in effervescent cleaners if this may damage the material.

Do not use:

- Avoid cleaning/immersing the device in sodium hypochlorite (chlorine bleach), hydrogen peroxide, mouthwash, or alcohol products, as these may damage or discolour the material.

6) Storage

When not in use, store the device as directed by the dental professional. Avoid excessive heat, pressure, or conditions that may deform or damage the device.

7) Expected Life / Service

With appropriate care and maintenance, the device is intended for ongoing use as directed by the dental professional. If the device becomes damaged or fit is no longer acceptable, the patient should contact their treating dental professional for evaluation.

8) Incident / Complaint Reporting

Report suspected serious incidents associated with this device to:

- The Manufacturer and/or UK Responsible Person (UKRP) using the contacts above, and
- The UK Medicines and Healthcare products Regulatory Agency (MHRA) — www.gov.uk/mhra

In accordance with The Medical Devices Regulations 2002 (as amended).

INFO

For full contact details of the local Authorized Representative, refer to: Authorized Representatives — end of this document.

9) Disposal

Dispose of the device in accordance with local regulations and the dental practice's procedures. If the device is contaminated with bodily fluids, treat it as clinical/infectious waste per local requirements.

EUROPEAN UNION (ENGLISH)

INSTRUCTIONS FOR USE (IFU)

European Union — Custom-Made Device (CMD)

GMDN / EMDN	Q01020601: Mobile Dental Prostheses
Product Code(s)	ACRFDENTRLN, DENTFLRLN, MILLDENTRLN
Manufacturer	Zima International, Inc. Also DBA Dandy 1320 N. 300 W., Lehi, UT 84043 United States
EU authorized representative (EURP)	Name: Dandy Labs Europe Company No.: Registration No.: 989857065 Address: 3, boulevard de Sebastopol, 75001 Paris, France
Intended User	Prescriber — Dental professional / dental practice
Device Type	Custom-made dental device (CMD)
Sterility	Not supplied sterile
Single-use / Reusable	Reusable

1) Intended Purpose (Intended Use)

This custom-made denture reliner is intended to improve denture fit/retention when a denture no longer fits as intended and requires relining, as prescribed by a licensed dental professional.

2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: Hard denture reline (full arch) (as applicable)
- Material: Sofreliner (Silicone Dioxide)

[Liner Material Safety Data Sheets \(Paste\)](#)[Liner Material Safety Data Sheets \(Primer\)](#)

3) Contraindications

Do not use in patients with:

- Known hypersensitivity to the device material
- Untreated oral lesions

4) Warnings / Precautions / Potential Risks

⚠ **WARNING**

Potential risks associated with denture lining materials may include:

- Gingival irritation, mucosal soreness, ulceration
- Poor adaptation may cause instability or trauma
- Allergic response
- Separation of material

If the patient experiences pain, irritation, or other adverse effects, discontinue use and contact the treating dental professional.

⚠ **CAUTION**

NA

5) Cleaning and Care

Daily cleaning:

- Brush gently using a soft-bristled denture brush or soft toothbrush.
- Avoid abrasive products (e.g., toothpaste, abrasive cleaners, strong soaps).
- Use mild, non-abrasive liquid dish soap as directed by the dental professional.
- Rinse thoroughly under cool or lukewarm running water.
- Avoid high heat; do not use hot or boiling water.

Soaking:

- Store/soak as directed by the dental professional to prevent drying.
- Avoid prolonged soaking in effervescent cleaners if this may damage the material.

Do not use:

- Avoid cleaning/immersing the device in sodium hypochlorite (chlorine bleach), hydrogen peroxide, mouthwash, or alcohol products, as these may damage or discolour the material.

6) Storage

When not in use, store the device as directed by the dental professional. Avoid excessive heat, pressure, or conditions that may deform or damage the device.

7) Expected Life / Service

With appropriate care and maintenance, the device is intended for ongoing use as directed by the dental professional. If the device becomes damaged or fit is no longer acceptable, the patient should contact their treating dental professional for evaluation.

8) Incident / Complaint Reporting

Report suspected serious incidents associated with this device to:

- The Manufacturer and/or Dandy EU entity using the contacts above, and corresponding local requirements
- The Competent Authority of your Member State, in accordance with Regulation (EU) 2017/745 (MDR).

Competent Authorities:

- France: Agence Nationale de Sécurité du Médicament et des Produits de Santé (ANSM)
- Spain : Spanish Agency for Medicines and Health Products (AEMPS)
- Germany : Bundesinstitut für Arzneimittel und Medizinprodukte (BfArM) (BfArM)
- Italy : Ministero della Salute / Istituto Superiore di Sanità (ISS)
- Spain: Agencia Española de Medicamentos y Productos Sanitarios (AEMPS)
- Netherlands: Inspectie Gezondheidszorg en Jeugd (IGJ)
- Austria: Bundesamt für Sicherheit im Gesundheitswesen (BASG)
- Sweden: Läkemedelsverket
- Denmark: Lægemiddelstyrelsen
- Portugal: Autoridade Nacional do Medicamento e Produtos de Saúde (INFARMED)
- Poland: Urząd Rejestracji Produktów Leczniczych, Wyrobów Medycznych i Produktów Biobójczych (URPL)
- Norway: Direktoratet for medisinske produkter (DMP)

INFO

For the full contact information of the local Authorized Representative, see: Authorized Representatives — at the end of the document.

9) Disposal

Dispose of the device in accordance with local regulations and the dental practice's procedures. If the device is contaminated with bodily fluids, treat it as clinical/infectious waste per local requirements.

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The table below lists the local Authorized Representative or Sponsor for each country / jurisdiction where this device is marketed. For incident reporting, contact the representative for your country.

Le tableau ci-dessous liste le Représentant Autorisé local pour chaque pays. | Die folgende Tabelle enthält die bevollmächtigten Vertreter je Land. | La tabella seguente elenca il rappresentante autorizzato o lo sponsor locale per ciascun paese. | La tabla siguiente enumera el Representante Autorizado o Patrocinador local para cada país. | Følgende tabel viser den autoriserede repræsentant eller lokale sponsor for hvert land. | De volgende tabel geeft een overzicht van de gemachtigde vertegenwoordiger of lokale sponsor voor elk land. | W poniższej tabeli wymieniono upoważnionych przedstawicieli lub lokalnych sponsorów dla każdego kraju. | Följande tabell listar den auktoriserade representanten eller den lokala sponsorn för varje land. | A tabela seguinte lista o Representante Autorizado ou Patrocinador local para cada país. | Tabellen nedenfor viser den autoriserte representanten eller den lokale sponsoren for hvert land.

Jurisdiction	Authorised Rep / Sponsor	Registered Address / Contact Email	Reg. Number	Regulatory Body
Australia (AUS)	DANDY Labs AUS Pty Ltd	Tower One, Level 46, 100 Barangaroo Ave, Sydney NSW 2000 Compliance@meetdandy.com	ABN: 33693955761	TGA (Therapeutic Goods Administration)
EU - English	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, France Compliance@meetdandy.com	Reg: 989857065	Competent Authority / EUDAMED
EU - España	Dandy Labs Europa	3, boulevard de Sebastopol, 75001 París, Francia Compliance@meetdandy.com	Reg: 989857065	AEMPS / Autoridad Competente
EU - France	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, France Compliance@meetdandy.com	Reg: 989857065	ANSM / Autorité Compétente
EU - Deutschland	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, Frankreich Compliance@meetdandy.com	Reg: 989857065	BfArM / Zuständige Behörde
EU - Italia	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Parigi, Francia Compliance@meetdandy.com	Reg: 989857065	Min. della Salute / ISS
EU - Denmark	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Parigi, Francia Compliance@meetdandy.com	Reg: 989857065	Lægemiddelstyrelsen
EU - Nederland	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Parigi, Francia Compliance@meetdandy.com	Reg: 989857065	Inspectie Gezondheidszorg en Jeugd (IGJ)
EU - Polska	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Parigi, Francia Compliance@meetdandy.com	Reg: 989857065	Urząd Rejestracji Produktów Leczniczych (UR PL)
EU - Sverige	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Parigi, Francia	Reg: 989857065	Läkemedelsverket

		Compliance@meetdandy.com		
EU - Portugal	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Parigi, Francia Compliance@meetdandy.com	Reg: 989857065	Autoridade Nacional do Medicamento e Produtos de Saúde (INFARMED)
EU - Norge	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Parigi, Francia Compliance@meetdandy.com	Reg: 989857065	Direktoratet for medisinske produkter (DMP)
United Kingdom (UK)	DANDY LABS GB, LTD	5 New Street Square, London EC4A 3TW Compliance@meetdandy.com	Co. No.: 16873608	MHRA
New Zealand (NZ)	DANDY Labs NZ Limited	Quigg Partners, Level 7, 36 Brandon Street, Wellington, NZ 6011	NZBN: 942905375844 6	Medsafe

NOTE

Contact details are subject to change. Always verify current information on the official product labelling or the manufacturer's website before reporting an incident. / Les coordonnées sont susceptibles de changer. / Angaben können sich ändern. / I recapiti sono soggetti a modifiche.

**Manufacturer / Hersteller /
Fabbricante / Fabricant /
Fabricante**

Zima International, Inc. Also DBA Dandy
1320 N. 300 W., Lehi, UT 84043, United States
compliance@meetdandy.com

This is a representation of an electronic record that was signed electronically and this page is the manifestation of the electronic signatures.

Signatory Table

Action Name	User Name	Title	Signature Date
Approve	Prasen Modi	Quality & Regulatory Affairs Specialist	07-Aug-2026 17:20

* Datetime fields are displayed according to the system time zone: (GMT-04:00) Eastern Daylight Time (America/New_York)