INSTRUCTIONS FOR USE (IFU)

Dandy / Zima International, Inc. | AUS · NZ · UK · EU/EN · EU/ES · EU/FR · EU/DE · EU/IT · EU/DA ·
EU/NL · EU/PL · EU/SV · EU/NB · EU/PT

MANUFACTURER

Manufacturer	Zima International, Inc. Also DBA Dandy
Address	22 Cortlandt St, 30th Floor, New York, NY, 10007, US
Contact / Website	www.meetdandy.com
Complaints	complaints@meetdandy.com

DEVICE IDENTITY

Device Identity

Device Family	Dental Occlusal Splint
GMDN / EMDN Code	43025 / Q0199
Product Code(s)	TMJNTI, TMJAD, TMJDBS, TMJGA
Device Type	Custom-made dental device (CMD) / Patient Matched Medical Device (PMMD)
Sterility	Not supplied sterile
Single-use / Reuse	Single-use
Intended User	Prescriber: Dental professional / dental practice

TABLE OF CONTENT

AUSTRALIA (ENGLISH)	6
1) Intended Purpose (Intended Use)	6
2) Device Description / Key Specifications	6
3) Contraindications	6
4) Warnings / Precautions / Potential Risks	6
5) Cleaning and Care	7
6) Storage	7
7) Expected Life / Service	7
8) Incident / Complaint Reporting	7
9) Disposal	8
NEW ZEALAND (ENGLISH)	9
1) Intended Purpose (Intended Use)	9
2) Device Description / Key Specifications	9
3) Contraindications	9
4) Warnings / Precautions / Potential Risks	9
5) Cleaning and Care	10
6) Storage	10
7) Expected Life / Service	10
8) Incident / Complaint Reporting	10
9) Disposal	11
UNITED KINGDOM (ENGLISH)	12
1) Intended Purpose (Intended Use)	12
2) Device Description / Key Specifications	12
3) Contraindications	12
4) Warnings / Precautions / Potential Risks	13
5) Cleaning and Care	13
6) Storage	13
7) Expected Life / Service	13
8) Incident / Complaint Reporting	13
9) Disposal	14
EUROPEAN UNION (ENGLISH)	15
1) Intended Purpose (Intended Use)	15
2) Device Description / Key Specifications	15
3) Contraindications	15
4) Warnings / Precautions / Potential Risks	16
5) Cleaning and Care	16
6) Storage	16
7) Expected Life / Service	16

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AUSTRALIA (ENGLISH)

INSTRUCTIONS FOR USE (IFU)

Australia — Patient Matched Medical Device (PMMD)

GMDN / EMDN	43025 / Q0199
Product Code(s)	TMJNTI, TMJAD, TMJDBS, TMJGA
Manufacturer	Zima International, Inc. Also DBA Dandy 22 Cortlandt St, 30th Floor, New York, NY, 10007 US
Australian Sponsor	Name: DANDY Labs AUS Pty Ltd Company No.: ABN: 33693955761 Address: Tower One - International Towers, Level 46, 100 Barangaroo Avenue, Sydney, New South Wales 2000, Australia
Intended User	Prescriber — Dental professional / dental practice
Device Type	Patient Matched Medical Device (PMMD)
Sterility	Not supplied sterile (if applicable)
Single-use / Reusable	Single-use

1) Intended Purpose (Intended Use)

Dandy TMJ Appliances are intended to reduce pain and discomfort in the TMJ and sometimes help re-train muscles to avoid these movements.

2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: The devices are worn while sleeping to support patients' TMJ in the position prescribed by the dentists. The appliances are removable by the patient.
- Material: Erkodent Disc

[Night Guard Material Safety Data Sheets](#)

3) Contraindications

For the select subset of patients who have surgery to treat TMD or any TMJ issues, any further treatment with TMJ appliances or splints should be prescribed by the surgeon.

4) Warnings / Precautions / Potential Risks

⚠ WARNING

- Tooth movement or changes in dental occlusion
- Gingival or dental soreness
- Pain or soreness to the temporomandibular joint
- Obstruction of oral breathing
- Excessive salivation

⚠ CAUTION

Dentists should consider the medical history of their patients. Including history of arthritis, surgery, impact trauma, any previous TMJ treatment, respiratory disorders, or other relevant health problems, and refer the patient to the appropriate healthcare provider before prescribing the device.

5) Cleaning and Care

Appliance Fit & Adjustments

- The dental staff or dentist should remove the appliance from its packaging and rinse it under water. Then, place the device in the patient's mouth, checking for the proper fit and function of the appliance. The prescribing dentist should always perform the appliance adjustments. Any adjustments in the fit of the device should easily be made chairside to ensure that the patient is pleased with the fit and function of their new appliance.

Operating Instructions

CLEANING

- Brush with a soft-bristle toothbrush before each use.

FIT ADJUSTMENT

- The prescribing dentist or physician should make all adjustments.

REMOVAL INSTRUCTIONS

- To remove the appliance properly, use both index fingers to pull down the front area of the upper appliance.

For Dawson B-splint,

- Do this for the upper first; second, for the lower, pull up in the front area.
- To avoid breakage, do not pull down the back of the appliance.

6) Storage

Dry and store in the provided case when not in use..

7) Expected Life / Service

Expect a 2-3 year product life with the appliance.

Appliance replacement signs are cracks in the acrylic and failure of the device to maintain retention in the patient's mouth.

8) Incident / Complaint Reporting

Report suspected serious incidents associated with this device to:

- The Manufacturer and/or Dandy Australia entity using the contacts above, and
- Therapeutic Goods Administration (TGA)

In Australia, serious incidents involving medical devices must be reported to the TGA in accordance with the Therapeutic Goods Act 1989.

INFO

For full contact details of the local Authorized Representative, refer to: Authorized Representatives — end of this document.

9) Disposal

Disposing of a TMJ appliance depends on its material and your location, but generally involves treating it as a small electronic/medical item, focusing on local rules for special waste.

If the item requires refabrication, Dandy requests that the item be returned so the hardware can be retrieved.

NEW ZEALAND (ENGLISH)

INSTRUCTIONS FOR USE (IFU)

New Zealand— Patient Matched Medical Device (PMMD)

GMDN / EMDN	43025 / Q0199
Product Code(s)	TMJNTI, TMJAD, TMJDBS, TMJGA
Manufacturer	Zima International, Inc. Also DBA Dandy 22 Cortlandt St, 30th Floor, New York, NY, 10007 US
New Zealand Sponsor	Name: DANDY Labs NZ Limited Address: Quigg Partners, Level 7, 36 Brandon Street, Wellington, NZ 6011
Intended User	Prescriber — Dental professional / dental practice
Device Type	Patient Matched Medical Device (PMMD)
Sterility	Not supplied sterile (if applicable)
Single-use / Reusable	Single-use

1) Intended Purpose (Intended Use)

Dandy TMJ Appliances are intended to reduce pain and discomfort in the TMJ and sometimes help re-train muscles to avoid these movements.

2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: The devices are worn while sleeping to support patients' TMJ in the position prescribed by the dentists. The appliances are removable by the patient.
- Material: Erkodent Disc

[Night Guard Material Safety Data Sheets](#)

3) Contraindications

For the select subset of patients who have surgery to treat TMD or any TMJ issues, any further treatment with TMJ appliances or splints should be prescribed by the surgeon.

4) Warnings / Precautions / Potential Risks

⚠ WARNING

- Tooth movement or changes in dental occlusion
- Gingival or dental soreness
- Pain or soreness to the temporomandibular joint
- Obstruction of oral breathing
- Excessive salivation

⚠ CAUTION

Dentists should consider the medical history of their patients. Including history of arthritis, surgery, impact trauma, any previous TMJ treatment, respiratory disorders, or other relevant health problems, and refer the patient to the appropriate healthcare provider before prescribing the device.

5) Cleaning and Care

Appliance Fit & Adjustments

- The dental staff or dentist should remove the appliance from its packaging and rinse it under water. Then, place the device in the patient's mouth, checking for the proper fit and function of the appliance. The prescribing dentist should always perform the appliance adjustments. Any adjustments in the fit of the device should easily be made chairside to ensure that the patient is pleased with the fit and function of their new appliance.

Operating Instructions

CLEANING

- Brush with a soft-bristle toothbrush before each use.

FIT ADJUSTMENT

- The prescribing dentist or physician should make all adjustments.

REMOVAL INSTRUCTIONS

- To remove the appliance properly, use both index fingers to pull down the front area of the upper appliance.

For Dawson B-splint,

- Do this for the upper first; second, for the lower, pull up in the front area.
- To avoid breakage, do not pull down the back of the appliance.

6) Storage

Dry and store in the provided case when not in use..

7) Expected Life / Service

Expect a 2-3 year product life with the appliance.

Appliance replacement signs are cracks in the acrylic and failure of the device to maintain retention in the patient's mouth.

8) Incident / Complaint Reporting

Report suspected serious incidents associated with this device to:

- The Manufacturer and/or Dandy New Zealand entity using the contacts above, and
- Medsafe

In New Zealand, serious incidents and quality issues involving medical devices should be reported to Medsafe in accordance with New Zealand's medical device adverse event reporting requirements.

INFO

For full contact details of the local Authorized Representative, refer to: Authorized Representatives — end of this document.

9) Disposal

Disposing of a TMJ appliance depends on its material and your location, but generally involves treating it as a small electronic/medical item, focusing on local rules for special waste.

If the item requires refabrication, Dandy requests that the item be returned so the hardware can be retrieved.

UNITED KINGDOM (ENGLISH)

INSTRUCTIONS FOR USE (IFU)

United Kingdom — Custom-Made Device (CMD)

GMDN / EMDN	43025 / Q0199
Product Code(s)	TMJNTI, TMJAD, TMJDBS, TMJGA
Manufacturer	Zima International, Inc. Also DBA Dandy 22 Cortlandt St, 30th Floor, New York, NY, 10007 US
UK Responsible Person (UKRP)	Name: DANDY LABS GB, LTD Company No.: Company Number: 16873608 Address: 5 New Street Square, London EC4A 3TW, United Kingdom
Intended User	Prescriber — Dental professional / dental practice
Device Type	Custom-made dental device (CMD)
Sterility	Not supplied sterile
Single-use / Reusable	Single-use

1) Intended Purpose (Intended Use)

Dandy TMJ Appliances are intended to reduce pain and discomfort in the TMJ and sometimes help re-train muscles to avoid these movements.

2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: The devices are worn while sleeping to support patients' TMJ in the position prescribed by the dentists. The appliances are removable by the patient.
- Material: Erkodent Disc

[Night Guard Material Safety Data Sheets](#)

3) Contraindications

For the select subset of patients who have surgery to treat TMD or any TMJ issues, any further treatment with TMJ appliances or splints should be prescribed by the surgeon.

4) Warnings / Precautions / Potential Risks

⚠ WARNING

- Tooth movement or changes in dental occlusion
- Gingival or dental soreness
- Pain or soreness to the temporomandibular joint
- Obstruction of oral breathing
- Excessive salivation

⚠ CAUTION

Dentists should consider the medical history of their patients. Including history of arthritis, surgery, impact trauma, any previous TMJ treatment, respiratory disorders, or other relevant health problems, and refer the patient to the appropriate healthcare provider before prescribing the device.

5) Cleaning and Care

Appliance Fit & Adjustments

- The dental staff or dentist should remove the appliance from its packaging and rinse it under water. Then, place the device in the patient's mouth, checking for the proper fit and function of the appliance. The prescribing dentist should always perform the appliance adjustments. Any adjustments in the fit of the device should easily be made chairside to ensure that the patient is pleased with the fit and function of their new appliance.

Operating Instructions

CLEANING

- Brush with a soft-bristle toothbrush before each use.

FIT ADJUSTMENT

- The prescribing dentist or physician should make all adjustments.

REMOVAL INSTRUCTIONS

- To remove the appliance properly, use both index fingers to pull down the front area of the upper appliance.

For Dawson B-splint,

- Do this for the upper first; second, for the lower, pull up in the front area.
- To avoid breakage, do not pull down the back of the appliance.

6) Storage

Dry and store in the provided case when not in use. .

7) Expected Life / Service

Expect a 2-3 year product life with the appliance.

Appliance replacement signs are cracks in the acrylic and failure of the device to maintain retention in the patient's mouth.

8) Incident / Complaint Reporting

Report suspected serious incidents associated with this device to:

- The Manufacturer and/or UK Responsible Person (UKRP) using the contacts above, and
- The UK Medicines and Healthcare products Regulatory Agency (MHRA) — www.gov.uk/mhra

In accordance with The Medical Devices Regulations 2002 (as amended).

INFO

For full contact details of the local Authorized Representative, refer to: Authorized Representatives — end of this document.

9) Disposal

- Disposing of a TMJ appliance depends on its material and your location, but generally involves treating it as a small electronic/medical item, focusing on local rules for special waste.
- If the item requires refabrication, Dandy requests that the item be returned so the hardware can be retrieved.

EUROPEAN UNION (ENGLISH)

INSTRUCTIONS FOR USE (IFU)

European Union — Custom-Made Device (CMD)

GMDN / EMDN	43025 / Q0199: ODONTOLOGY DEVICES - OTHER
Product Code(s)	TMJNTI, TMJAD, TMJDBS, TMJGA
Manufacturer	Zima International, Inc. Also DBA Dandy 22 Cortlandt St, 30th Floor, New York, NY, 10007 US
EU authorized representative (EURP)	Name: Dandy Labs Europe Company No.: Registration No.: 989857065 Address: 3, boulevard de Sebastopol, 75001 Paris, France
Intended User	Prescriber — Dental professional / dental practice
Device Type	Custom-made dental device (CMD)
Sterility	Not supplied sterile
Single-use / Reusable	Single-use

1) Intended Purpose (Intended Use)

Dandy TMJ Appliances are intended to reduce pain and discomfort in the TMJ and sometimes help re-train muscles to avoid these movements.

2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: The devices are worn while sleeping to support patients' TMJ in the position prescribed by the dentists. The appliances are removable by the patient.
- Material: Erkodent Disc

[Night Guard Material Safety Data Sheets](#)

3) Contraindications

For the select subset of patients who have surgery to treat TMD or any TMJ issues, any further treatment with TMJ appliances or splints should be prescribed by the surgeon.

4) Warnings / Precautions / Potential Risks

⚠ WARNING

- Tooth movement or changes in dental occlusion
- Gingival or dental soreness
- Pain or soreness to the temporomandibular joint
- Obstruction of oral breathing
- Excessive salivation

⚠ CAUTION

Dentists should consider the medical history of their patients. Including history of arthritis, surgery, impact trauma, any previous TMJ treatment, respiratory disorders, or other relevant health problems, and refer the patient to the appropriate healthcare provider before prescribing the device.

5) Cleaning and Care

Appliance Fit & Adjustments

- The dental staff or dentist should remove the appliance from its packaging and rinse it under water. Then, place the device in the patient's mouth, checking for the proper fit and function of the appliance. The prescribing dentist should always perform the appliance adjustments. Any adjustments in the fit of the device should easily be made chairside to ensure that the patient is pleased with the fit and function of their new appliance.

Operating Instructions

CLEANING

- Brush with a soft-bristle toothbrush before each use.

FIT ADJUSTMENT

- The prescribing dentist or physician should make all adjustments.

REMOVAL INSTRUCTIONS

- To remove the appliance properly, use both index fingers to pull down the front area of the upper appliance.
- For Dawson B-splint,
- Do this for the upper first; second, for the lower, pull up in the front area.
- To avoid breakage, do not pull down the back of the appliance.

6) Storage

Dry and store in the provided case when not in use..

7) Expected Life / Service

Expect a 2-3 year product life with the appliance.

Appliance replacement signs are cracks in the acrylic and failure of the device to maintain retention in the patient's mouth.

8) Incident / Complaint Reporting

Report suspected serious incidents associated with this device to:

- The Manufacturer and/or Dandy EU entity using the contacts above, and corresponding local requirements
- The Competent Authority of your Member State, in accordance with Regulation (EU) 2017/745 (MDR).

Competent Authorities:

- France: Agence Nationale de Sécurité du Médicament et des Produits de Santé (ANSM)
- Spain : Spanish Agency for Medicines and Health Products (AEMPS)
- Germany : Bundesinstitut für Arzneimittel und Medizinprodukte (BfArM) (BfArM)
- Italy : Ministero della Salute / Istituto Superiore di Sanità (ISS)
- Netherlands: Inspectie Gezondheidszorg en Jeugd (IGJ)
- Austria: Bundesamt für Sicherheit im Gesundheitswesen (BASG)
- Sweden: Läkemedelsverket
- Denmark: Lægemiddelstyrelsen
- Portugal: Autoridade Nacional do Medicamento e Produtos de Saúde (INFARMED)
- Ireland: Health Products Regulatory Authority (HPRA)
- Poland: Urząd Rejestracji Produktów Leczniczych, Wyrobów Medycznych i Produktów Biobójczych (URPL)
- Norway: Direktoratet for medisinske produkter (DMP)

INFO

For the full contact information of the local Authorized Representative, see: Authorized Representatives — at the end of the document.

9) Disposal

- Disposing of a TMJ appliance depends on its material and your location, but generally involves treating it as a small electronic/medical item, focusing on local rules for special waste.
- If the item requires refabrication, Dandy requests that the item be returned so the hardware can be retrieved.

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The table below lists the local Authorized Representative or Sponsor for each country / jurisdiction where this device is marketed. For incident reporting, contact the representative for your country.

Le tableau ci-dessous liste le Représentant Autorisé local pour chaque pays. | Die folgende Tabelle enthält die bevollmächtigten Vertreter je Land. | La tabella seguente elenca il rappresentante autorizzato o lo sponsor locale per ciascun paese. | La tabla siguiente enumera el Representante Autorizado o Patrocinador local para cada país. | Følgende tabel viser den autoriserede repræsentant eller lokale sponsor for hvert land. | De volgende tabel geeft een overzicht van de gemachtigde vertegenwoordiger of lokale sponsor voor elk land. | W poniższej tabeli wymieniono upoważnionych przedstawicieli lub lokalnych sponsorów dla każdego kraju. | Följande tabell listar den auktoriserade representanten eller den lokala sponsorn för varje land. | A tabela seguinte lista o Representante Autorizado ou Patrocinador local para cada país. | Tabellen nedenfor viser den autoriserte representanten eller den lokale sponsoren for hvert land.

Jurisdiction	Authorised Rep / Sponsor	Registered Address / Contact Email	Reg. Number	Regulatory Body
Australia (AUS)	DANDY Labs AUS Pty Ltd	Tower One, Level 46, 100 Barangaroo Ave, Sydney NSW 2000 Compliance@meetdandy.com	ABN: 33693955761	TGA (Therapeutic Goods Administration)
EU - English	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, France Compliance@meetdandy.com	Reg: 989857065	Competent Authority / EUDAMED
EU - España	Dandy Labs Europa	3, boulevard de Sebastopol, 75001 Paris, Francia Compliance@meetdandy.com	Reg: 989857065	AEMPS / Autoridad Competente
EU - France	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, France Compliance@meetdandy.com	Reg: 989857065	ANSM / Autorité Compétente
EU - Deutschland	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, Frankreich Compliance@meetdandy.com	Reg: 989857065	BfArM / Zuständige Behörde
EU - Italia	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Parigi, Francia Compliance@meetdandy.com	Reg: 989857065	Min. della Salute / ISS
EU - Austria	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, France Compliance@meetdandy.com	Reg: 989857065	Bundesamt für Sicherheit im Gesundheitswesen (BASG)
EU - Ireland	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, France Compliance@meetdandy.com	Reg: 989857065	Health Products Regulatory Authority (HPRA)
EU - Denmark	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, France Compliance@meetdandy.com	Reg: 989857065	Lægemiddelstyrelsen

EU - Nederland	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, France Compliance@meetdandy.com	Reg: 989857065	Inspectie Gezondheidszorg en Jeugd (IGJ)
EU - Polska	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, France Compliance@meetdandy.com	Reg: 989857065	Urząd Rejestracji Produktów Leczniczych (URPL)
EU - Sverige	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, France Compliance@meetdandy.com	Reg: 989857065	Läkemedelsverket
EU - Portugal	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, France Compliance@meetdandy.com	Reg: 989857065	Autoridade Nacional do Medicamento e Produtos de Saúde (INFARMED)
EU - Norge	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, France Compliance@meetdandy.com	Reg: 989857065	Direktoratet for medisinske produkter (DMP)
United Kingdom (UK)	DANDY LABS GB, LTD	5 New Street Square, London EC4A 3TW Compliance@meetdandy.com	Co. No.: 16873608	MHRA
New Zealand (NZ)	DANDY Labs NZ Limited	Quigg Partners, Level 7, 36 Brandon Street, Wellington, NZ 6011	NZBN: 9429053758446	Medsafe

NOTE

Contact details are subject to change. Always verify current information on the official product labelling or the manufacturer's website before reporting an incident. / Les coordonnées sont susceptibles de changer. / Angaben können sich ändern. / I recapiti sono soggetti a modifiche.

**Manufacturer / Hersteller /
Fabbricante / Fabricant /
Fabricante**

Zima International, Inc. Also DBA Dandy
22 Cortlandt St, 30th Floor, New York, NY, 10007, US
compliance@meetdandy.com

This is a representation of an electronic record that was signed electronically and this page is the manifestation of the electronic signatures.

Signatory Table

Action Name	User Name	Title	Signature Date
Approve	Garvin Tran	Manager, Operations - QMS	28-Aug-2026 20:34

* Datetime fields are displayed according to the system time zone: (GMT-04:00) Eastern Daylight Time (America/New_York)