

# dandy

## Chairside Software

### Instructions for Use

Issued on: 2026-09

CE 3022

|                                                         |           |
|---------------------------------------------------------|-----------|
| <b>1. Device Description</b>                            | <b>2</b>  |
| 1.1. Device details                                     | 2         |
| 1.2. Intended Purpose                                   | 2         |
| 1.3. Intended Patient population                        | 2         |
| 1.4. Contraindications                                  | 2         |
| 1.5. Intended User                                      | 2         |
| 1.6. Intended Use Environment                           | 2         |
| 1.7. How Chairside works                                | 3         |
| 1.7.1. General description                              | 3         |
| 1.7.2. Clinical Benefits                                | 3         |
| 1.8. AI/ML—functions                                    | 3         |
| 1.8.1. General description                              | 3         |
| 1.8.2. AI/ML models                                     | 4         |
| 1.8.3. How it is used                                   | 4         |
| 1.8.4. Clinical Benefits                                | 4         |
| 1.8.5. AI/ML performance                                | 5         |
| 1.8.6. Quantitative performance results                 | 5         |
| 1.8.7. Interoperability                                 | 5         |
| <b>2. Warnings &amp; Precautions</b>                    | <b>6</b>  |
| 2.1. Explanation                                        | 6         |
| 2.2. List of Warnings and cautions                      | 6         |
| <b>3. Residual Risks</b>                                | <b>7</b>  |
| <b>4. Data Privacy &amp; Cybersecurity</b>              | <b>7</b>  |
| 4.1. Data privacy                                       | 7         |
| 4.2. Cybersecurity                                      | 8         |
| 4.2.1. Security Patches & Updates                       | 9         |
| <b>5. System &amp; Network Requirements</b>             | <b>9</b>  |
| 5.1. System Requirements                                | 9         |
| 5.2. Network Requirements                               | 9         |
| <b>6. Installation &amp; Operation</b>                  | <b>9</b>  |
| 6.1. Installation                                       | 9         |
| 6.2. Operation Workflow                                 | 9         |
| 6.3. AI Scan Review Instructions:                       | 9         |
| 6.4. Software decommissioning                           | 10        |
| <b>7. Accessing the Electronic IFU (eIFU):</b>          | <b>10</b> |
| <b>8. Troubleshooting &amp; Support</b>                 | <b>10</b> |
| 8.1. Troubleshooting                                    | 10        |
| 8.2. Support                                            | 10        |
| <b>9. Serious Incident Reporting</b>                    | <b>10</b> |
| <b>10. Manufacturer &amp; Authorized Representative</b> | <b>11</b> |
| 10.1. Legal Manufacturer                                | 11        |
| 10.2. EU Authorized Representative (EC REP)             | 11        |
| 10.3. EU Importer's Name and Address                    | 11        |
| <b>11. Product Label</b>                                | <b>11</b> |
| <b>12. Symbols Glossary</b>                             | <b>11</b> |

## 1. Device Description

### 1.1. Device details

The current software version can be found within the application by navigating to the 'About' screen.

|     |                                     |
|-----|-------------------------------------|
| #   | Version 1.113.0                     |
| UDI | 081020609CHAIRSIDENZ (Basic UDI-DI) |

### 1.2. Intended Purpose

Dandy Chairside supports trained dental professionals in determining whether a digital impression acquired for a restorative or prosthetic case is technically complete and of sufficient quality to support fabrication of a custom-made dental restoration. Operating on scan data received from a compatible intraoral scanner, the software analyzes the scan for defined technical quality characteristics and presents the resulting findings to the user. The user reviews those findings together with the scan and decides whether to proceed with the scan or to acquire a new one. The software supports, and does not replace, the dental professional's judgment on scan acceptability. It is intended for use by trained dental professionals in dental clinics and dental laboratories. It does not diagnose oral conditions, determine treatment, or design the restoration.

Chairside is not intended for the diagnosis of oral conditions or identification of the required treatment.

### 1.3. Intended Patient population

The intended patient population consists of patients undergoing intraoral digital scanning as part of restorative or prosthetic dental procedures.

This includes patients receiving dental restorations such as crowns, bridges, implants, inlays, onlays, or other prosthetic devices that require a digital impression or intraoral scan as part of the clinical workflow.

The Chairside device does not directly interact with the patient and does not analyze the patient's health status or disease.

### 1.4. Contraindications

There are no device-specific contraindications.

### 1.5. Intended User

Licensed dentists, Dental hygienists, Dental assistants, Dental laboratory staff.

Users are expected to have appropriate training in dental workflows and the use of intraoral scanning systems. Chairside is not intended for use by untrained personnel.

### 1.6. Intended Use Environment

Chairside is intended for use in Dental clinics and Dental laboratories.

## 1.7. How Chairside works

### 1.7.1. General description

Dandy Chairside is a software medical device intended to support digital dental workflows by enabling the creation, management, and submission of dental cases based on intraoral scan data acquired from compatible third-party scanning systems. The software is designed for use in dental clinics and laboratories to facilitate coordination between clinicians and dental laboratories during restorative and prosthetic procedures.

Chairside operates as a desktop application deployed on managed Windows-based systems. It interfaces with external intraoral scanners through a standardized integration layer, allowing users to initiate and guide scanning workflows, visualize and manipulate three-dimensional (3D) scan data, and prepare digital impressions for submission.

Chairside includes an optional software feature, AI Scan Review, which provides automated, algorithm-based analysis of intraoral scan data to assess technical quality and completeness. This feature identifies potential scan issues, such as incomplete data capture, gaps, or unclear margin regions, and presents advisory indicators and optional recommendations for rescan or refinement. AI Scan Review may also generate a preliminary margin line as a visual aid; however, this output is intended solely as a starting point for user review and refinement.

The device operates within a distributed system architecture consisting of a local application (Chairside), a version management component (Chairside Launcher), and cloud-based backend services (Orthly Web). The backend services provide centralized data management, case processing, and execution of machine learning algorithms associated with AI Scan Review. Secure data transmission is used to communicate between the local application and backend services.

### 1.7.2. Clinical Benefits

Dandy Chairside does not provide diagnostic or therapeutic benefit directly; the software improves patient management by enabling a digital workflow that allows clinicians to review and modify the scan in real-time.

Description of the Clinical Benefits for the AI/ML functions is provided under 1.8.5.

**Note:** Dandy Chairside is a workflow optimization tool. It does not provide direct diagnostic or therapeutic benefits.

## 1.8. AI/ML-functions

### 1.8.1. General description

AI Scan Review is an optional, advisory feature of Dandy Chairside that automatically analyses intraoral scan data to help the clinician assess the technical quality and completeness of a scan before a case is submitted. It runs on Dandy's managed backend using a fixed ("locked") set of models and algorithms that do not learn or change during use; any change is delivered only as a new, validated software version. AI Scan Review is non-diagnostic: it does not interpret clinical images, diagnose disease, detect caries or pathology, or influence the care pathway. Its outputs are non-binding aids, colour-coded overlays, highlighted regions and a findings list, and the licensed clinician remains the sole decision-maker, reviewing every scan and remaining free to accept, edit or reject the software's findings.

### 1.8.2. AI/ML models

AI Scan Review combines machine-learning models with deterministic geometric algorithms across the following functions:

- **Jaw segmentation (supporting model).** A 3D convolutional neural network classifies the scan surface into tooth, gingiva and preparation regions. It is the shared upstream input on which the clearance, undercut and margin-line functions depend.
- **Unclear margins / margin-line prediction.** A machine-learning model predicts the preparation margin and fits an editable margin-line spline, presented as a visual aid the clinician can accept or redraw.
- **Occlusal clearance detection.** A geometric algorithm compares the prepared surface against the opposing jaw relative to the minimum material thickness and flags teeth with insufficient clearance via a colour-coded overlay.
- **Path-of-insertion (POI) undercut detection.** A geometric algorithm computes candidate insertion axes and alerts the user when no valid path clears the margin.
- **Gaps analysis (denture scans).** A geometric algorithm aligns the scan to a template and highlights critical versus non-critical regions of missing data.

### 1.8.3. How it is used

Within the normal workflow the clinician loads the completed scan, runs AI Scan Review, inspects the highlighted low-quality or incomplete regions and the findings list, re-scans the areas indicated where needed, and then approves or rejects the scan for submission. The margin-line prediction can be reviewed and edited directly on the 3D model.

|                                                                                     |         |                                                                                                                                                                                                                                                                                                     |
|-------------------------------------------------------------------------------------|---------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | WARNING | Because its detection thresholds are deliberately tuned to limit false positives, the function does not flag every genuine scan-quality issue; an absence of a flag must not be interpreted as confirmation that a scan is adequate, and the clinician must still verify scan quality independently |
|  | WARNING | The predicted margin line is an advisory aid and may require correction. AI Scan Review is non-diagnostic and does not detect caries or pathology, it must not be used for any diagnostic purpose.                                                                                                  |

### 1.8.4. AI/ML-input requirements

To obtain results which are representative for clinical use, the following AI/ML data input requirements are specified:

- Intra-oral scans obtained with a supported intra-oral scanner specified in section 1.8.8
- Intra-oral scans which fulfil the respective quality criteria:
  - Data should be complete
  - Free from clinical artifacts and noise
  - Sufficient coverage and accurate occlusal registration
  - A standard CAD format (e.g. OBJ, STL or PLY)

### 1.8.5. Clinical Benefits

The following clinical benefits are defined when using the Dandy AI/ML-functions:

The AI/ML functions benefit users by providing the clinician with accurate information related to the quality of the scan.

|                                                                                   |         |                                                                                                                                                                   |
|-----------------------------------------------------------------------------------|---------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | Caution | AI Scan Review is intended to support the clinician in assessing scan quality prior to submission. It does not provide direct diagnostic or therapeutic benefits. |
|-----------------------------------------------------------------------------------|---------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------|

### 1.8.6. AI/ML performance

The below table visualizes the performance of the several AI/ML-models and heuristics models embedded within the Dandy Chairside software. Performance characteristics may change over time as specified in section 1.8.6, performance will exceed the metrics displayed below.

| Metric                     | Performance specification |
|----------------------------|---------------------------|
| Accuracy                   | ≥ 90%                     |
| Sensitivity / Recall (TPR) | ≥ 50%                     |
| Specificity (TNR)          | ≥ 90%                     |
| Precision (PPV)            | ≥ 85%                     |

**Note:** Separate performance specifications apply to ML models: Jaw segmentation MeanIoU ≥ 80% (0.80); Margin-line prediction F1 ≥ 70%.

### 1.8.7. Quantitative performance results

For each model embedded within the Dandy Chairside software, performance has been tested on an independent test dataset. For each of the models the test dataset has different sizes and characteristics. The performance results are shown in the table below. Each of the results exceeds the intended performance specification set out in section 1.8.5.

| Function (config)                                                      | Metric             | Result       | 95% CI (Wilson)              |
|------------------------------------------------------------------------|--------------------|--------------|------------------------------|
| <b>Occlusal Clearance</b><br>(0.4 mm leeway;<br>n=152 pos / n=100 neg) | Accuracy           | 94.0%        | [90.3%, 96.3%] (n=252)       |
|                                                                        | Precision (PPV)    | 97.3%        | [92.2%, 99.1%] (n=96)        |
|                                                                        | Sensitivity (TPR)  | 61.0%        | [53.1%, 68.4%] (n=152)       |
|                                                                        | Specificity (TNR)  | 97.5%        | [92.3%, 99.2%] (n=100)       |
| <b>POI Undercut</b> (0.05 threshold; n~100 pos / n~100 neg)            | Accuracy           | 94.6%        | [88.7%, 97.5%] (n=110)       |
|                                                                        | Precision (PPV)    | 91.5%        | [81.7%, 96.3%] (n=59)        |
|                                                                        | Sensitivity (TPR)  | 54.0%        | [44.3%, 63.4%] (n=100)       |
|                                                                        | Specificity (TNR)  | 95.0%        | [88.8%, 97.8%] (n=100)       |
| <b>Jaw segmentation</b><br>(supporting ML, n= > 1,000)                 | Tooth-only MeanIoU | 0.92         | n/a                          |
|                                                                        | Prep-only MeanIoU  | 0.83         | n/a                          |
| <b>Margin-line prediction</b><br>(ML; n=1,196)                         | F1 score @ 70 µm   | 0.71         | Tolerance used for the claim |
| <b>Gaps Analysis</b><br>(supporting ML)                                | Accuracy           | <b>98.2%</b> | [97.2%, 98.8%] (n=1099)      |
|                                                                        | Precision (PPV)    | <b>87.5%</b> | [71.9%, 95.0%] (n=32)        |

| Function (config) | Metric            | Result       | 95% CI (Wilson)         |
|-------------------|-------------------|--------------|-------------------------|
|                   | Sensitivity (TPR) | <b>63.6%</b> | [48.9%, 76.2%] (n=44)   |
|                   | Specificity (TNR) | <b>99.6%</b> | [99.0%, 99.9%] (n=1055) |

### 1.8.8. Interoperability

Dandy Chairside Software integrates with intraoral scanners to support digital workflows. The software bundle delivered to customers is an Electron-based application provided by Dandy.

Chairside supports two scanner-specific application variants:

- Trios Chairside (integrating with 3Shape Unite Scanner)
  - 3Shape Trios 3 Scanner
  - 3Shape Trios 5 Scanner
- Vision Chairside (integrating with AlliedStar ScanPro)
  - Dandy Vision Scanner

When there is no scanner connected, or a scanner is connected that Chairside does not support, the application clearly communicates that there is no scanner connected.

|                                                                                   |                |                                                                                              |
|-----------------------------------------------------------------------------------|----------------|----------------------------------------------------------------------------------------------|
|  | <b>WARNING</b> | Dandy Chairside should not be used with unsupported scanner manufacturer types and variants. |
|-----------------------------------------------------------------------------------|----------------|----------------------------------------------------------------------------------------------|

## 2. Warnings & Precautions

Before using this software, you must have read and understood all instructions to ensure safe and correct use. Be aware of the following warnings and cautions when using the software.

### 2.1. Explanation

|                                                                                     |                |                                                                                                                               |
|-------------------------------------------------------------------------------------|----------------|-------------------------------------------------------------------------------------------------------------------------------|
|  | <b>WARNING</b> | A warning informs the user of potential risks and details the expected behaviour to avoid injury, user harm or device damage. |
|  | <b>Caution</b> | A caution indicates a situation in which an operator requires action or awareness in order to avoid undesirable consequences. |

### 2.2. List of Warnings and cautions

|                                                                                     |                |                                                                                                                                                                                                                                                                                                      |
|-------------------------------------------------------------------------------------|----------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | <b>WARNING</b> | The AI Scan results are provided for information purposes only, results generated are not intended for direct diagnosis or for direct therapeutic decision making. The results require review and correction (as needed) by qualified dental professionals.                                          |
|  | <b>WARNING</b> | The system does not detect caries or pathology.                                                                                                                                                                                                                                                      |
|  | <b>WARNING</b> | Because its detection thresholds are deliberately tuned to limit false positives, the function does not flag every genuine scan-quality issue; an absence of a flag must not be interpreted as confirmation that a scan is adequate, and the clinician must still verify scan quality independently. |
|  | <b>WARNING</b> | The predicted margin line is an advisory aid and may require correction. AI Scan Review is non-diagnostic and does not detect caries or pathology, it must not be used for any diagnostic purpose.                                                                                                   |

|                                                                                     |         |                                                                                                                                                                                                                                                                      |
|-------------------------------------------------------------------------------------|---------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | WARNING | Dandy Chairside should not be used with unsupported scanner manufacturer types and variants.                                                                                                                                                                         |
|    | WARNING | The software is solely intended to inform qualified dental professionals; the software is not intended as a replacement for clinical diagnosis or therapeutic decision making.                                                                                       |
|    | Caution | Network interruptions may delay case submission.                                                                                                                                                                                                                     |
|    | Caution | Verify patient identity and case information.                                                                                                                                                                                                                        |
|    | Caution | AI Scan Review is intended to support the clinician in assessing scan quality prior to submission. It does not provide direct diagnostic or therapeutic benefits.                                                                                                    |
|    | Caution | Provision of the software on the preinstalled computerised system by Dandy should only be used for the intended purpose as described within this Instructions for Use.                                                                                               |
|    | Caution | It is strictly prohibited to adjust any of the system's security settings; adjustment of such settings places patient safety and data confidentiality, integrity and availability at risk.                                                                           |
|   | Caution | It is your responsibility to ensure that you connect the system to the internet over secure connections only; using insecure connections places patient safety and data confidentiality, integrity and availability at risk.                                         |
|  | Caution | In the event of any suspected loss of confidentiality, availability or integrity caused by the use of the Dandy software, you are responsible for immediately informing the Dandy organization through the contact information provided under the section 'support'. |
|  | Caution | In the event of any suspicious or unexpected software behaviour you are responsible for immediately informing the Dandy organization through the contact information provided under the section 'support'.                                                           |
| Important note                                                                      |         | Re-scan if data corruption or upload failure occurs.                                                                                                                                                                                                                 |

### 3. Residual Risks

It is the responsibility of the user to take note of the warnings and cautions available within these instructions for use.

|                                                                                     |         |                                                                                                                                                                                |
|-------------------------------------------------------------------------------------|---------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | Warning | The software is solely intended to inform qualified dental professionals, the software is not intended as a replacement for clinical diagnosis or therapeutic decision making. |
|-------------------------------------------------------------------------------------|---------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

There are no known undesirable side-effects.

## 4. Data Privacy & Cybersecurity

### 4.1. Data privacy

The software is operated remotely by the Dandy organisation. As data processing will include processing personal health information (PHI) to achieve its intended purpose, Dandy qualifies as a Data Processor of Personal Data under regulatory jurisdictions, such as the EU and the UK. Data

processing is executed under assignment by the Data Controller through a Data Processing Agreement.

|                                                                                   |         |                                                                                                                                                                                                                                                                      |
|-----------------------------------------------------------------------------------|---------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | Caution | In the event of any suspected loss of confidentiality, availability or integrity caused by the use of the Dandy software, you are responsible for immediately informing the Dandy organization through the contact information provided under the section 'support'. |
|-----------------------------------------------------------------------------------|---------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

As the use of the software you bear the responsibility to secure the systems that you use to avoid breach of confidentiality. You can take the following actions to further strengthen secure use of the Dandy software:

- Never document your access credentials and make use of a password management system
- Do not share your access credentials with any other person, instead create or request the creation of additional user accounts
- Always terminate active user sessions when leaving a mobile device with access to the software unattended

Where individual patient data is deleted from the Dandy Chairside software, the data is deleted within a month from any remaining back-ups within Dandy. Once permanently deleted data can no longer be restored or reconstructed otherwise.

## 4.2. Cybersecurity

In addition to the measures described above, the Dandy organization is committed to securing the software to the latest standards, the following security measures have been implemented:

- PHI when stored locally is encrypted at rest via OS disk encryption (e.g., BitLocker)
- The computer system runs a anti-virus and malware software (e.g., Microsoft Defender)
- Data transmitted to cloud services uses TLS 1.2+ encryption
- Cloud databases and storage are encrypted at rest
- Users authenticate with individual accounts; backend services use JWT tokens
- Role-Based Access Control (RBAC) limits access to authorized functions
- PHI is not shared across clinics; each practice can only access its own data
- AI Scan Review workflow transmits anonymized scan data (no PHI)
- Cloud infrastructure hosted on Google Cloud Platform (GCP) with IAM, audit logging, and Cloud Armor protections

|                                                                                     |         |                                                                                                                                                                                                                              |
|-------------------------------------------------------------------------------------|---------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | Caution | Provision of the software on the preinstalled computerised system by Dandy should only be used for the intended purpose as described within this Instructions for Use.                                                       |
|  | Caution | It is strictly prohibited to adjust any of the system's security settings, adjustment of such settings places patient safety and data confidentiality, integrity and availability at risk.                                   |
|  | Caution | In the event of any suspicious or unexpected software behaviour you are responsible for immediately informing the Dandy organization through the contact information provided under the section 'support'.                   |
|  | Caution | It is your responsibility to ensure that you connect the system to the internet over secure connections only, using insecure connections places patient safety and data confidentiality, integrity and availability at risk. |

#### 4.2.1. Security Patches & Updates

The organization will, from time to time implement security and non-security related updates to ensure on-going performance and security of the software system. Users will be notified by Dandy upon pending updates which may affect availability of the Dandy Chairside software.

#### 4.2.2. Reporting vulnerabilities

If there are suspected vulnerabilities identified within the Dandy Chairside software, we urge you to report such vulnerabilities to [bugbounty@meetdandy.com](mailto:bugbounty@meetdandy.com) or [vulnerabilities@meetdandy.com](mailto:vulnerabilities@meetdandy.com).

## 5. System & Network Requirements

### 5.1. System Requirements

Chairside runs on a managed workstation distributed and controlled via Microsoft Intune MDM. Minimum workstation requirements:

- Windows 11pro (64-bit)
- Intel i7 CPU or better
- 32 GB RAM
- 477GB Storage HDD
- Administrative privileges for installation (required for MDM enrollment)

### 5.2. Network Requirements

- Chairside communicates with Dandy backend services over secure outbound connections
- Stable internet connection
- Recommended minimum bandwidth: 20 Mbps download / 5 Mbps upload (scan uploads)
- Network access to TCP port 443 outbound (HTTPS)
- TLS 1.2+ is required for all connections to the Dandy cloud infrastructure
- Chairside does not require inbound firewall exceptions.

## 6. Installation & Operation

### 6.1. Installation

The full installation of the Dandy Chairside software is performed by Dandy. There are no actions required on the end of the user to start using the software. Once the user starts the use of the software, the user is provided with the relevant user credentials and required to change their password in order to start using the software.

### 6.2. Operation Workflow

- Launch the application and select the patient.
- Acquire an intraoral scan.
- Review scan quality.
- Submit the case to the lab.

### 6.3. AI Scan Review Instructions:

- Load completed scan.
- Run AI Scan Review.
- Inspect highlighted low-quality regions.
- Re-scan as needed.
- Approve or reject scan.

## 6.4. Software decommissioning

The software is provided pre-installed on a Dandy owned computer system. Upon replacement or termination of the cooperation with Dandy, the computer system is to be returned to Dandy.

Dandy is responsible for the full in-house decommissioning of the software from the computer system and will ensure the full and irreversible deletion of any customer and or PHI data remaining on the computer system, within the period specified by the Data Controller within the Data Processing Agreement.

## 7. Accessing the Electronic IFU (eIFU):

- A copy of this IFU is available on this website: <http://dandy.dental/chairside-product-ifu>
- A printed copy of this IFU can be requested and will be provided free of charge within 7 calendar days.

## 8. Troubleshooting & Support

### 8.1. Troubleshooting

| Event                  | Resolution                                                                                             |
|------------------------|--------------------------------------------------------------------------------------------------------|
| Failure to submit case | <ul style="list-style-type: none"> <li>• Check internet connection</li> </ul>                          |
| Unable to log in       | <ul style="list-style-type: none"> <li>• Verify credentials</li> <li>• Request new password</li> </ul> |
| Missing scan data      | <ul style="list-style-type: none"> <li>• Confirm scanner import completed</li> </ul>                   |

### 8.2. Support

In the event of any detected software failures or bugs, or in the event of any suspected security incident, data breach, or potential patient harm, the user is responsible for immediately contacting Dandy to request for further support. Within the software a direct link to the support channels of Dandy is available (phone, email and support form) to reach out to the Dandy team.

Alternatively, you can contact Dandy directly via the details provided in Section 10 of this instructions for use.

## 9. Serious Incident Reporting

The user is responsible for reporting any potential concerns with regards to the performance, security or privacy of the software. Such concerns must be reported to Dandy and their authorised representative at the contact details provided in section 10 (specifically section 10.1 & 10.2).

In addition, in the event of a serious incident (for example, use of the software has (potentially) led to serious injury that has occurred in relation to the software must be reported to the competent authority of the member state in which the user and/or patient is established.

The following website includes an overview with the relevant contact details for all competent authorities within the European Union:

[https://health.ec.europa.eu/document/download/900ad9e7-fc79-4617-93be-f6712a3e3306\\_en](https://health.ec.europa.eu/document/download/900ad9e7-fc79-4617-93be-f6712a3e3306_en)

## 10. Manufacturer & Authorized Representative

### 10.1. Legal Manufacturer

|                                                                                   |                                                                                                                                                                                                                                                            |
|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | <p>Zima International, Inc. d/b/a Dandy<br/>22 Cortlandt St, 30th Floor,<br/>New York,<br/>NY 10007<br/>USA</p> <p><b>Direct contact details:</b><br/>Phone +1-866-310-8665<br/>Email <a href="mailto:Doctors@meetdandy.com">Doctors@meetdandy.com</a></p> |
|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

### 10.2. EU Authorized Representative (EC REP)

|                                                                                   |                                                                                    |
|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------|
|  | <p>Dandy Labs Europe<br/>3, boulevard de Sébastopol<br/>75001 Paris<br/>France</p> |
|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------|

### 10.3. EU Importer's Name and Address

|                                                                                     |                                                                                    |
|-------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|
|  | <p>Dandy Labs Europe<br/>3, boulevard de Sébastopol<br/>75001 Paris<br/>France</p> |
|-------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|

## 11. Product Label

An example of the product label which can be located within the software is presented below.

|                                                                                                                                                                                                                                                                                 |                                                                                                                                                                                  |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p><b>Dandy Chairside</b> </p> <p>Software version: 1.0.0</p>                                                                                                                                | <p><b>CE 3022</b></p>                                                                                                                                                            |
| <p> <b>Zima International, Inc. d/b/a Dandy</b><br/>22 Cortlandt St, 30th Floor<br/>New York, NY 10007, USA</p>                                                                              | <p> <b>Dandy Labs Europe</b><br/>3 boulevard de Sébastopol<br/>75001 Paris, France</p>        |
| <p> <b>Date of manufacture:</b><br/><b>2026-04</b></p>                                                                                                                                       | <p> <b>Importer: Dandy Labs Europe</b><br/>3 boulevard de Sébastopol, 75001 Paris, France</p> |
| <p> <b>UDI-DI:</b> (01) 00850012345008<br/><b>UDI-PI:</b> (8012) 1.0.0</p>                                                                                                                   |                                                                                                                                                                                  |
| <p><b>(EN) Chairside dental workflow software — Class IIa medical device software</b><br/>Software for managing digital dental workflows and coordinating intraoral scan data acquired by third-party scanning systems. Includes AI Scan Review (advisory, non-diagnostic).</p> |                                                                                                                                                                                  |
| <p> <b>Consult instructions for use (eIFU):</b><br/><a href="http://dandy.dental/chairside-product-ifu">http://dandy.dental/chairside-product-ifu</a></p>                                    | <p> <b>Caution</b><br/>Consult the IFU for warnings &amp; precautions.</p>                    |

Note: This representation is shown for clarification purposes only. The applicable device label is available in the "About" box within the software application.

## 12. Symbols Glossary

| Symbol                                                                              | Title                                                                   | Description                                                                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | Conformité Européenne                                                   | The software complies with the European legislation with regards to medical devices (MDR 2017/745). The Notified Body number is displayed adjacent to the CE marking (3022).                                                       |
|    | Consult instructions for use or consult electronic instructions for use | Indicates the need for the user to consult the instructions for use                                                                                                                                                                |
|    | Legal manufacturer                                                      | Indicates the medical device manufacturer                                                                                                                                                                                          |
|    | Medical Device                                                          | Indicates the item is a medical device                                                                                                                                                                                             |
|   | Authorized representative in the European Community/ European Union     | Indicates the authorized representative in the European Community/ European Union                                                                                                                                                  |
|  | Date of manufacture                                                     | Indicates the date when the medical device was manufactured. For Dandy this refers to the software release date.                                                                                                                   |
|  | Model number                                                            | Indicates the model number or type number of a product. For Dandy this refers to the software version number.                                                                                                                      |
|  | Importer                                                                | Indicates the entity importing the medical device into the locale                                                                                                                                                                  |
|  | Caution                                                                 | Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences |
|  | Unique device identifier                                                | Indicates a carrier that contains unique device identifier information                                                                                                                                                             |



This is a representation of an electronic record that was signed electronically and this page is the manifestation of the electronic signatures.

Signatory Table

| Action Name                  | User Name   | Title                                   | Signature Date    |
|------------------------------|-------------|-----------------------------------------|-------------------|
| Send for Review (Written By) | Prasen Modi | Quality & Regulatory Affairs Specialist | 02-Sep-2026 14:14 |
| Review                       | Garvin Tran | Manager, Operations - QMS               | 03-Sep-2026 15:22 |
| Approve                      | Garvin Tran | Manager, Operations - QMS               | 03-Sep-2026 15:25 |

\* Datetime fields are displayed according to the system time zone: (GMT-04:00) Eastern Daylight Time (America/New\_York)