INSTRUCTIONS FOR USE (IFU)

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MANUFACTURER

Manufacturer	Zima International, Inc. Also DBA Dandy
Address	1320 N. 300 W., Lehi, UT 84043, United States
Contact / Website	www.meetdandy.com
Complaints	complaints@meetdandy.com

DEVICE IDENTITY

Device Identity

Device Family	Partial Denture
GMDN / EMDN Code	17845
Product Code(s)	ACRPARDENT, ACTLPARDENT, CCPARDENT, DURPARDENT, TITPARDENT, TRYPARDENT, VALPARDENT, PDWAXRIM, PDHYBRID
Device Type	Custom-made dental device (CMD) / Patient Matched Medical Device (PMMD)
Sterility	Not supplied sterile
Single-use / Reuse	Reusable
Intended User	Prescriber: Dental professional / dental practice

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AUSTRALIA (ENGLISH)**INSTRUCTIONS FOR USE (IFU)**

Australia — Patient Matched Medical Device (PMMD)

GMDN / EMDN	17845
Product Code(s)	ACRPARDENT, ACTLPARDENT, CCPARDENT, DURPARDENT, TITPARDENT, TRYPARDENT, VALPARDENT, PDWAXRIM, PDHYBRID
Manufacturer	Zima International, Inc. Also DBA Dandy 1320 N. 300 W., Lehi, UT 84043 United States
Australian Sponsor	Name: DANDY Labs AUS Pty Ltd Company No.: ABN: 33693955761 Address: Tower One - International Towers, Level 46, 100 Barangaroo Avenue, Sydney, New South Wales 2000, Australia
Intended User	Prescriber — Dental professional / dental practice
Device Type	Patient Matched Medical Device (PMMD)
Sterility	Not supplied sterile
Single-use / Reusable	Reusable

1) Intended Purpose (Intended Use)

Partial dentures are a type of removable dental prosthetic device used to replace one or more missing teeth. They are typically made of either a metal, acrylic, or a flexible framework with replacement teeth that are attached to an acrylic base, which is colored to match the gums

2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: Dandy acrylic partial dentures are custom-made to fit your mouth and are designed to be as comfortable and unobtrusive as possible.
- Material: Chrome Cobalt, Titanium, Acrylic, TCS, Valplast, Duraflex, Acetal, PMMA

[Partial Dentures Material Safety Data Sheets](#)

3) Contraindications

Partial dentures should not be used if you have any allergies to acrylic, plastic, or certain metals. Partial dentures should not be used in patients who have severe gum disease.

4) Warnings / Precautions / Potential Risks

WARNING

Risks include sore spots, discomfort, and possible damage to remaining natural teeth.

CAUTION

NA

5) Cleaning and Care

Directions for Use

- To ensure proper use of acrylic partial dentures, follow these steps:
- Ensure that the partial denture is clean and free of any debris.
- Rinse the patient's mouth with warm water before inserting the partial denture.
- Gently place the partial denture in the mouth, and use your fingers to guide it into place.
- Avoid using excessive force or pushing the partial denture too far.
- Instruct the patient to gently bite down on the partial denture to ensure proper fit and stability.
- Use carbon articulation paper to adjust any heavy occlusion points.
- Provide specific instructions to the patient on how long to wear it each day and when to remove it.

Cleaning instructions

- Dandy recommends cleaning your partial dentures daily before and after every use. Immediately after removal, use warm water (not to exceed 45°C/113°F) and soap to remove remaining saliva and debris from your partial dentures.

Avoid

- Soaking your partial dentures in any liquid for longer than the manufacturer of the cleaning solution recommends.
- Soaking your partial dentures in hard water. Over time, exposure to tap water containing elevated chlorine levels could tarnish the metal.
- Do not immerse your partial dentures in water exceeding 45°C/113°F to prevent warping.
- Use only approved cleaners for meta/acrylic/flexible partial dentures and avoid abrasive toothpastes, as they can damage the device's appearance.
- Soaking your partial dentures in a bleach-based cleaning solution this may tarnish them.

Do

- Dispose of the cleaning solution after every use. Rinse your partial dentures well before use.

6) Storage

Store the partial dentures in a clean, dry place when not in use. Avoid exposing partial dentures to excessive heat or pressure, as this may cause warping or other damage.

7) Expected Life / Service

Partial dentures are designed to last several years with proper care and maintenance. Patients should contact the dental practice for evaluation and repair if dentures become damaged or lose fit.

8) Incident / Complaint Reporting

Report suspected serious incidents associated with this device to:

- The Manufacturer and/or Dandy Australia entity using the contacts above, and Document State: Effective (Simple)

- Therapeutic Goods Administration (TGA)

In Australia, serious incidents involving medical devices must be reported to the TGA in accordance with the Therapeutic Goods Act 1989.

INFO

For full contact details of the local Authorized Representative, refer to: Authorized Representatives — end of this document.

9) Disposal

- Partial denture disposal usually involves handling it as regulated biohazardous or general waste, often via specialized handlers for incineration or secure landfill.
- Always avoid putting them in regular trash without checking local rules.
- Biohazard/Infectious Waste : Treat used or unclaimed dentures as biohazardous due to saliva/blood; use red bags and authorized carriers.
- Incineration : A common method for final destruction.
- Confidentiality: Redact or destroy patient info on lab boxes and bases before disposal.

EUROPEAN UNION (ENGLISH)**INSTRUCTIONS FOR USE (IFU)**

European Union — Custom-Made Device (CMD)

GMDN / EMDN	Q01020601 - Mobile Dental Prostheses
Product Code(s)	ACRPARDENT, ACTLPARDENT, CCPARDENT, DURPARDENT, TITPARDENT, TRYPARDENT, VALPARDENT, PDWAXRIM, PDHYBRID
Manufacturer	Zima International, Inc. Also DBA Dandy 1320 N. 300 W., Lehi, UT 84043 United States
EU authorized representative (EURP)	Name: Dandy Labs Europe Company No.: Registration No.: 989857065 Address: 3, boulevard de Sebastopol, 75001 Paris, France
Intended User	Prescriber — Dental professional / dental practice
Device Type	Custom-made dental device (CMD)
Sterility	Not supplied sterile
Single-use / Reusable	Reusable

1) Intended Purpose (Intended Use)

Partial dentures are a type of removable dental prosthetic device used to replace one or more missing teeth. They are typically made of either a metal, acrylic, or a flexible framework with replacement teeth that are attached to an acrylic base, which is colored to match the gums

2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: Dandy acrylic partial dentures are custom-made to fit your mouth and are designed to be as comfortable and unobtrusive as possible.
- Material: Chrome Cobalt, Titanium, Acrylic, TCS, Valplast, Duraflex, Acetal, PMMA

[Partial Dentures Material Safety Data Sheets](#)

3) Contraindications

Partial dentures should not be used if you have any allergies to acrylic, plastic, or certain metals. Partial dentures should not be used in patients who have severe gum disease.

4) Warnings / Precautions / Potential Risks

WARNING

Risks include sore spots, discomfort, and possible damage to remaining natural teeth.

CAUTION

NA

5) Cleaning and Care

Directions for Use

- To ensure proper use of acrylic partial dentures, follow these steps:
- Ensure that the partial denture is clean and free of any debris.
- Rinse the patient's mouth with warm water before inserting the partial denture.
- Gently place the partial denture in the mouth, and use your fingers to guide it into place.
- Avoid using excessive force or pushing the partial denture too far.
- Instruct the patient to gently bite down on the partial denture to ensure proper fit and stability.
- Use carbon articulation paper to adjust any heavy occlusion points.
- Provide specific instructions to the patient on how long to wear it each day and when to remove it.

Cleaning instructions

- Dandy recommends cleaning your partial dentures daily before and after every use. Immediately after removal, use warm water (not to exceed 45°C/113°F) and soap to remove remaining saliva and debris from your partial dentures.

Avoid

- Soaking your partial dentures in any liquid for longer than the manufacturer of the cleaning solution recommends.
- Soaking your partial dentures in hard water. Over time, exposure to tap water containing elevated chlorine levels could tarnish the metal.
- Do not immerse your partial dentures in water exceeding 45°C/113°F to prevent warping.
- Use only approved cleaners for meta/acrylic/flexible partial dentures and avoid abrasive toothpastes, as they can damage the device's appearance.
- Soaking your partial dentures in a bleach-based cleaning solution this may tarnish them.

Do

- Dispose of the cleaning solution after every use. Rinse your partial dentures well before use.

6) Storage

Store the partial dentures in a clean, dry place when not in use. Avoid exposing partial dentures to excessive heat or pressure, as this may cause warping or other damage.

7) Expected Life / Service

Partial dentures are designed to last several years with proper care and maintenance. Patients should contact the dental practice for evaluation and repair if dentures become damaged or lose fit.

8) Incident / Complaint Reporting

Report suspected serious incidents associated with this device to:

- The Manufacturer and/or Dandy EU entity using the contacts above, and corresponding local requirements

- The Competent Authority of your Member State, in accordance with Regulation (EU) 2017/745 (MDR).

Competent Authorities:

- France: Agence Nationale de Sécurité du Médicament et des Produits de Santé (ANSM)
- Spain : Spanish Agency for Medicines and Health Products (AEMPS)
- Germany : Bundesinstitut für Arzneimittel und Medizinprodukte (BfArM) (BfArM)
- Italy : Ministero della Salute / Istituto Superiore di Sanità (ISS)
- Spain: Agencia Española de Medicamentos y Productos Sanitarios (AEMPS) for Spain
- Netherlands: Inspectie Gezondheidszorg en Jeugd (IGJ)
- Austria: Bundesamt für Sicherheit im Gesundheitswesen (BASG)
- Sweden: Läkemedelsverket
- Denmark: Lægemiddelstyrelsen
- Portugal: Autoridade Nacional do Medicamento e Produtos de Saúde (INFARMED)

INFO

For the full contact information of the local Authorized Representative, see: Authorized Representatives — at the end of the document.

9) Disposal

- Partial denture disposal usually involves handling it as regulated biohazardous or general waste, often via specialized handlers for incineration or secure landfill.
- Always avoid putting them in regular trash without checking local rules.
- Biohazard/Infectious Waste : Treat used or unclaimed dentures as biohazardous due to saliva/blood; use red bags and authorized carriers.
- Incineration : A common method for final destruction.
- Confidentiality: Redact or destroy patient info on lab boxes and bases before disposal.

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UNITED KINGDOM (ENGLISH)**INSTRUCTIONS FOR USE (IFU)**

United Kingdom — Custom-Made Device (CMD)

GMDN / EMDN	17845
Product Code(s)	ACRPARDENT, ACTLPARDENT, CCPARDENT, DURPARDENT, TITPARDENT, TRYPARDENT, VALPARDENT, PDWAXRIM, PDHYBRID
Manufacturer	Zima International, Inc. Also DBA Dandy 1320 N. 300 W., Lehi, UT 84043 United States
UK Responsible Person (UKRP)	Name: DANDY LABS GB, LTD Company No.: Company Number: 16873608 Address: 5 New Street Square, London EC4A 3TW, United Kingdom
Intended User	Prescriber — Dental professional / dental practice
Device Type	Custom-made dental device (CMD)
Sterility	Not supplied sterile
Single-use / Reusable	Reusable

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- The Manufacturer and/or UK Responsible Person (UKRP) using the contacts above, and

Document State: Effective (Simple)

Page 23 of 26

Confidential

FRM-0002002, revision 2

- The UK Medicines and Healthcare products Regulatory Agency (MHRA) — www.gov.uk/mhra

In accordance with The Medical Devices Regulations 2002 (as amended).

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The table below lists the local Authorized Representative or Sponsor for each country / jurisdiction where this device is marketed. For incident reporting, contact the representative for your country.

La tabla siguiente enumera el Representante Autorizado o Patrocinador local para cada país. | Le tableau ci-dessous liste le Représentant Autorisé local pour chaque pays. | Die folgende Tabelle enthält die bevollmächtigten Vertreter je Land.

Jurisdiction	Authorised Rep / Sponsor	Registered Address / Contact Email	Reg. Number	Regulatory Body
Australia (AUS)	DANDY Labs AUS Pty Ltd	Tower One, Level 46, 100 Barangaroo Ave, Sydney NSW 2000 Compliance@meetdandy.com	ABN: 33693955761	TGA (Therapeutic Goods Administration)
EU — English	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, France Compliance@meetdandy.com	Reg: 989857065	Competent Authority / EUDAMED
EU — España	Dandy Labs Europa	3, boulevard de Sebastopol, 75001 Paris, Francia Compliance@meetdandy.com	Reg: 989857065	AEMPS / Autoridad Competente
EU — France	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, France Compliance@meetdandy.com	Reg: 989857065	ANSM / Autorité Compétente
EU — Deutschland	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, Frankreich Compliance@meetdandy.com	Reg: 989857065	BfArM / Zuständige Behörde
EU — Italia	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Parigi, Francia Compliance@meetdandy.com	Reg: 989857065	Min. della Salute / ISS
United Kingdom (UK)	DANDY LABS GB, LTD	5 New Street Square, London EC4A 3TW Compliance@meetdandy.com	Co. No.: 16873608	MHRA

NOTE

Contact details are subject to change. Always verify current information on the official product labelling or the manufacturer's website before reporting an incident. / Les coordonnées sont susceptibles de changer. / Angaben können sich ändern. / I recapiti sono soggetti a modifiche.

**Manufacturer / Hersteller /
Fabbicante / Fabricant /
Fabricante**

Zima International, Inc. Also DBA Dandy
1320 N. 300 W., Lehi, UT 84043, United States
compliance@meetdandy.com

This is a representation of an electronic record that was signed electronically and this page is the manifestation of the electronic signatures.

Signatory Table

Action Name	User Name	Title	Signature Date
Approve	Kasha Lim	Supplier Quality Engineer	29-Jun-2026 10:38

* Datetime fields are displayed according to the system time zone: (GMT-04:00) Eastern Daylight Time (America/New_York)