# INSTRUCTIONS FOR USE (IFU)

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## MANUFACTURER

|                   |                                               |
|-------------------|-----------------------------------------------|
| Manufacturer      | Zima International, Inc. Also DBA Dandy       |
| Address           | 1320 N. 300 W., Lehi, UT 84043, United States |
| Contact / Website | www.meetdandy.com                             |
| Complaints        | complaints@meetdandy.com                      |

## DEVICE IDENTITY

### Device Identity

|                    |                                                                         |
|--------------------|-------------------------------------------------------------------------|
| Device Family      | Surgical Guide, STD Abutment Insertion Guide                            |
| GMDN / EMDN Code   | 57949                                                                   |
| Product Code(s)    | SAIG                                                                    |
| Device Type        | Custom-made dental device (CMD) / Patient Matched Medical Device (PMMD) |
| Sterility          | Not supplied sterile                                                    |
| Single-use / Reuse | Single-use                                                              |
| Intended User      | Prescriber: Dental professional / dental practice                       |

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**AUSTRALIA (ENGLISH)****INSTRUCTIONS FOR USE (IFU)**

Australia — Patient Matched Medical Device (PMMD)

|                              |                                                                                                                                                                                                              |
|------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>GMDN / EMDN</b>           | 57949                                                                                                                                                                                                        |
| <b>Product Code(s)</b>       | SAIG                                                                                                                                                                                                         |
| <b>Manufacturer</b>          | <b>Zima International, Inc. Also DBA Dandy</b><br>1320 N. 300 W., Lehi, UT 84043<br>United States                                                                                                            |
| <b>Australian Sponsor</b>    | <b>Name: DANDY Labs AUS Pty Ltd</b><br><b>Company No.:</b> ABN: 33693955761<br><b>Address:</b> Tower One - International Towers, Level 46, 100<br>Barangaroo Avenue, Sydney, New South Wales 2000, Australia |
| <b>Intended User</b>         | Prescriber — Dental professional / dental practice                                                                                                                                                           |
| <b>Device Type</b>           | Patient Matched Medical Device (PMMD)                                                                                                                                                                        |
| <b>Sterility</b>             | Not supplied sterile                                                                                                                                                                                         |
| <b>Single-use / Reusable</b> | Single use                                                                                                                                                                                                   |

**1) Intended Purpose (Intended Use)**

The Insertion Guide is intended for use with single or multiple abutments during placement. It is intended for short-term use and is then discarded.

**2) Device Description / Key Specifications**

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: Abutment placement on a dental implant or on multiple implants in the correct orientation, performed by a dental professional.
- Material: Acrylic - monomer and powder

[Lucitone 199](#)[Monomer](#)[Powder](#)**3) Contraindications**

Allergy to or hypersensitivity to the used materials

## 4) Warnings / Precautions / Potential Risks

### **WARNING**

Surgery in the oral cavity and oral rehabilitation include general risks.

**These include (but are not limited to):**

- Small components accidentally dropped in the patient's mouth can be swallowed or aspirated.
- Allergy or hypersensitivity to the chemical ingredients of the material used.

### **CAUTION**

- The Abutment Insertion Guide should not be exposed to thermal treatment such as steam or dry heat sterilization.
- This product is intended for single use and may not be reused. Reusing the product introduces:
- The risk of infection and complications due to a loss of accuracy of fit and the precision of the components.
- Performance failure of the device.
- Infection, inflammation, or allergy related to the device.
- Swallowing or inhalation of the device or a part thereof.

## 5) Cleaning and Care

Do not autoclave or use dry heat sterilization as these may affect the mechanical properties. Chemical disinfection is recommended before use.

E.g.. Cidex OPA Solution - 0.55% ortho-phthaldehyde

Use

- 1) Check the osseointegration of the implant/implants to be restored
- 2) Try in the Insertion Guide
- 3) Before placing the abutment/abutments, try the guide in the patient's mouth. Ensure a tight fit and a secure fit.
- 4) Placement of the abutment in the guide
- 5) Place the abutment/abutments in the insertion guide and ensure proper fit and orientation (timing). There should be a slight resistance while inserting the abutment /abutments into the guide, and it should hold them securely.
- 6) Abutment insertion
- 7) Make sure to hold the guide gently.
- 8) Place the guide on the patient's teeth and seat it in place.
- 9) Each abutment should adapt to the implant's orientation.
- 10) Insert each abutment/abutment into the appropriate implant.
- 11) If resistance is felt, ensure abutment/abutments are aligned with the implant.
- 12) Abutment Screw Tightening
- 13) While holding the guide in place, insert the Abutment Screw (should be located in the model) and tighten.
- 14) Take a radiograph to confirm seating.
- 15) Final tightening
- 16) Tighten the abutment screw using a torque wrench according to the implant manufacturer's torque specifications.
- 17) If the soft tissue creates hard resistance and blanches during seating, it is recommended that you proceed slowly before final torquing of the screw.
- 18) Cover the screw head
  - a) Use your preferred method.
  - b) Dandy recommends:
    - i) Either a material similar to plumber's tape or a cotton pellet.
- 19) Then cover with a material similar to Cavit.

## 20) Cementation

- a) Cement the final restoration to the abutment using your preferred method (make sure that the cement method is consistent with the restoration type).

## 6) Storage

The products should be stored in their packaging in a dry place at a normal temperature (18–25°C/64–77°F). Use the sterilized components within the stated time period from the sterile bag.

## 7) Expected Life / Service

The dental abutment placement guide itself is a single-use or short-term tool for precise implant placement

## 8) Incident / Complaint Reporting

**Report suspected serious incidents associated with this device to:**

- The Manufacturer and/or Dandy Australia entity using the contacts above, and
- Therapeutic Goods Administration (TGA)

In Australia, serious incidents involving medical devices must be reported to the TGA in accordance with the Therapeutic Goods Act 1989.

### **INFO**

*For full contact details of the local Authorized Representative, refer to: Authorized Representatives — end of this document.*

## 9) Disposal

Comply with the current applicable national waste disposal regulations.

## EUROPEAN UNION (ENGLISH)

### INSTRUCTIONS FOR USE (IFU)

European Union — Custom-Made Device (CMD)

|                                            |                                                                                                                                                      |
|--------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>GMDN / EMDN</b>                         | Q010399- Surgical Dental Devices - Other                                                                                                             |
| <b>Product Code(s)</b>                     | SAIG                                                                                                                                                 |
| <b>Manufacturer</b>                        | <b>Zima International, Inc. Also DBA Dandy</b><br>1320 N. 300 W., Lehi, UT 84043<br>United States                                                    |
| <b>EU authorized representative (EURP)</b> | <b>Name:</b> Dandy Labs Europe<br><b>Company No.:</b> Registration No.: 989857065<br><b>Address:</b> 3, boulevard de Sebastopol, 75001 Paris, France |
| <b>Intended User</b>                       | Prescriber — Dental professional / dental practice                                                                                                   |
| <b>Device Type</b>                         | Custom-made dental device (CMD)                                                                                                                      |
| <b>Sterility</b>                           | Not supplied sterile                                                                                                                                 |
| <b>Single-use / Reusable</b>               | Single use                                                                                                                                           |

#### 1) Intended Purpose (Intended Use)

The Insertion Guide is intended for use with single or multiple abutments during placement. It is intended for short-term use and is then discarded.

#### 2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: Abutment placement on a dental implant or on multiple implants in the correct orientation, performed by a dental professional.
- Material: Keyguide

[KeyGuide](#)

#### 3) Contraindications

Allergy to or hypersensitivity to the used materials

#### 4) Warnings / Precautions / Potential Risks

**⚠ WARNING**

Surgery in the oral cavity and oral rehabilitation include general risks.

**These include (but are not limited to):**

- Small components accidentally dropped in the patient's mouth can be swallowed or aspirated.
- Allergy or hypersensitivity to the chemical ingredients of the material used.

**⚠ CAUTION**

- The Abutment Insertion Guide should not be exposed to thermal treatment such as steam or dry heat sterilization.
- This product is intended for single use and may not be reused. Reusing the product introduces:
- The risk of infection and complications due to a loss of accuracy of fit and the precision of the components.
- Performance failure of the device.
- Infection, inflammation, or allergy related to the device.
- Swallowing or inhalation of the device or a part thereof.

## 5) Cleaning and Care

Do not autoclave or use dry heat sterilization as these may affect the mechanical properties. Chemical disinfection is recommended before use.

E.g., Cidex OPA Solution - 0.55% ortho-phthalaldehyde

Use

- 21) Check the osseointegration of the implant/implants to be restored
- 22) Try in the Insertion Guide
- 23) Before placing the abutment/abutments, try the guide in the patient's mouth. Ensure a tight fit and a secure fit.
- 24) Placement of the abutment in the guide
- 25) Place the abutment/abutments in the insertion guide and ensure proper fit and orientation (timing). There should be a slight resistance while inserting the abutment /abutments into the guide, and it should hold them securely.
- 26) Abutment insertion
- 27) Make sure to hold the guide gently.
- 28) Place the guide on the patient's teeth and seat it in place.
- 29) Each abutment should adapt to the implant's orientation.
- 30) Insert each abutment/abutment into the appropriate implant.
- 31) If resistance is felt, ensure abutment/abutments are aligned with the implant.
- 32) Abutment Screw Tightening
- 33) While holding the guide in place, insert the Abutment Screw (should be located in the model) and tighten.
- 34) Take a radiograph to confirm seating.
- 35) Final tightening
- 36) Tighten the abutment screw using a torque wrench according to the implant manufacturer's torque specifications.
- 37) If the soft tissue creates hard resistance and blanches during seating, it is recommended that you proceed slowly before final torquing of the screw.
- 38) Cover the screw head
  - a) Use your preferred method.
  - b) Dandy recommends:
    - i) Either a material similar to plumber's tape or a cotton pellet.
- 39) Then cover with a material similar to Cavit.
- 40) Cementation
  - a) Cement the final restoration to the abutment using your preferred method (make sure that the cement method is consistent with the restoration type).

## 6) Storage

The products should be stored in their packaging in a dry place at a normal temperature (18–25°C/64–77°F). Use the sterilized components within the stated time period from the sterile bag.

## 7) Expected Life / Service

The dental abutment placement guide itself is a single-use or short-term tool for precise implant placement

## 8) Incident / Complaint Reporting

### Report suspected serious incidents associated with this device to:

- The Manufacturer and/or Dandy EU entity using the contacts above, and corresponding local requirements
- The Competent Authority of your Member State, in accordance with Regulation (EU) 2017/745 (MDR).

### Competent Authorities:

- France: Agence Nationale de Sécurité du Médicament et des Produits de Santé (ANSM)
- Spain : Spanish Agency for Medicines and Health Products (AEMPS)
- Germany : Bundesinstitut für Arzneimittel und Medizinprodukte (BfArM) (BfArM)
- Italy : Ministero della Salute / Istituto Superiore di Sanità (ISS)
- Spain: Agencia Española de Medicamentos y Productos Sanitarios (AEMPS) for Spain
- Netherlands: Inspectie Gezondheidszorg en Jeugd (IGJ)
- Austria: Bundesamt für Sicherheit im Gesundheitswesen (BASG)
- Sweden: Läkemedelsverket
- Denmark: Lægemiddelstyrelsen
- Portugal: Autoridade Nacional do Medicamento e Produtos de Saúde (INFARMED)

### INFO

For the full contact information of the local Authorized Representative, see: Authorized Representatives — at the end of the document.

## 9) Disposal

Comply with the current applicable national waste disposal regulations.

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**UNITED KINGDOM (ENGLISH)****INSTRUCTIONS FOR USE (IFU)**

United Kingdom — Custom-Made Device (CMD)

|                                     |                                                                                                                                                         |
|-------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>GMDN / EMDN</b>                  | 57949                                                                                                                                                   |
| <b>Product Code(s)</b>              | SAIG                                                                                                                                                    |
| <b>Manufacturer</b>                 | <b>Zima International, Inc. Also DBA Dandy</b><br>1320 N. 300 W., Lehi, UT 84043<br>United States                                                       |
| <b>UK Responsible Person (UKRP)</b> | <b>Name:</b> DANDY LABS GB, LTD<br><b>Company No.:</b> Company Number: 16873608<br><b>Address:</b> 5 New Street Square, London EC4A 3TW, United Kingdom |
| <b>Intended User</b>                | Prescriber — Dental professional / dental practice                                                                                                      |
| <b>Device Type</b>                  | Custom-made dental device (CMD)                                                                                                                         |
| <b>Sterility</b>                    | Not supplied sterile                                                                                                                                    |
| <b>Single-use / Reusable</b>        | Reusable                                                                                                                                                |

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Use

- 41) Check the osseointegration of the implant/implants to be restored
- 42) Try in the Insertion Guide
- 43) Before placing the abutment/abutments, try the guide in the patient's mouth. Ensure a tight fit and a secure fit.
- 44) Placement of the abutment in the guide
- 45) Place the abutment/abutments in the insertion guide and ensure proper fit and orientation (timing). There should be a slight resistance while inserting the abutment /abutments into the guide, and it should hold them securely.
- 46) Abutment insertion
- 47) Make sure to hold the guide gently.
- 48) Place the guide on the patient's teeth and seat it in place.
- 49) Each abutment should adapt to the implant's orientation.
- 50) Insert each abutment/abutment into the appropriate implant.
- 51) If resistance is felt, ensure abutment/abutments are aligned with the implant.
- 52) Abutment Screw Tightening
- 53) While holding the guide in place, insert the Abutment Screw (should be located in the model) and tighten.
- 54) Take a radiograph to confirm seating.
- 55) Final tightening
- 56) Tighten the abutment screw using a torque wrench according to the implant manufacturer's torque specifications.
- 57) If the soft tissue creates hard resistance and blanches during seating, it is recommended that you proceed slowly before final torquing of the screw.
- 58) Cover the screw head
  - a) Use your preferred method.
  - b) Dandy recommends:
    - i) Either a material similar to plumber's tape or a cotton pellet.
- 59) Then cover with a material similar to Cavit.

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## 8) Incident / Complaint Reporting

Report suspected serious incidents associated with this device to:

- The Manufacturer and/or UK Responsible Person (UKRP) using the contacts above, and
- The UK Medicines and Healthcare products Regulatory Agency (MHRA) — [www.gov.uk/mhra](http://www.gov.uk/mhra)

In accordance with The Medical Devices Regulations 2002 (as amended).

**INFO**

For full contact details of the local Authorized Representative, refer to: Authorized Representatives — end of this document.

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The table below lists the local Authorized Representative or Sponsor for each country / jurisdiction where this device is marketed. For incident reporting, contact the representative for your country.

*La tabla siguiente enumera el Representante Autorizado o Patrocinador local para cada país. | Le tableau ci-dessous liste le Représentant Autorisé local pour chaque pays. | Die folgende Tabelle enthält die bevollmächtigten Vertreter je Land.*

| Jurisdiction        | Authorised Rep / Sponsor | Registered Address / Contact Email                                                   | Reg. Number       | Regulatory Body                        |
|---------------------|--------------------------|--------------------------------------------------------------------------------------|-------------------|----------------------------------------|
| Australia (AUS)     | DANDY Labs AUS Pty Ltd   | Tower One, Level 46, 100 Barangaroo Ave, Sydney NSW 2000<br>Compliance@meetdandy.com | ABN: 33693955761  | TGA (Therapeutic Goods Administration) |
| EU — English        | Dandy Labs Europe        | 3, boulevard de Sebastopol, 75001 Paris, France<br>Compliance@meetdandy.com          | Reg: 989857065    | Competent Authority / EUDAMED          |
| EU — España         | Dandy Labs Europa        | 3, boulevard de Sebastopol, 75001 París, Francia<br>Compliance@meetdandy.com         | Reg: 989857065    | AEMPS / Autoridad Competente           |
| EU — France         | Dandy Labs Europe        | 3, boulevard de Sebastopol, 75001 Paris, France<br>Compliance@meetdandy.com          | Reg: 989857065    | ANSM / Autorité Compétente             |
| EU — Deutschland    | Dandy Labs Europe        | 3, boulevard de Sebastopol, 75001 Paris, Frankreich<br>Compliance@meetdandy.com      | Reg: 989857065    | BfArM / Zuständige Behörde             |
| EU — Italia         | Dandy Labs Europe        | 3, boulevard de Sebastopol, 75001 Parigi, Francia<br>Compliance@meetdandy.com        | Reg: 989857065    | Min. della Salute / ISS                |
| United Kingdom (UK) | DANDY LABS GB, LTD       | 5 New Street Square, London EC4A 3TW<br>Compliance@meetdandy.com                     | Co. No.: 16873608 | MHRA                                   |

### NOTE

Contact details are subject to change. Always verify current information on the official product labelling or the manufacturer's website before reporting an incident. / Les coordonnées sont susceptibles de changer. / Angaben können sich ändern. / I recapiti sono soggetti a modifiche.

**Manufacturer / Hersteller /  
Fabbicante / Fabricant /  
Fabricante**

**Zima International, Inc. Also DBA Dandy**  
1320 N. 300 W., Lehi, UT 84043, United States  
compliance@meetdandy.com

This is a representation of an electronic record that was signed electronically and this page is the manifestation of the electronic signatures.

Signatory Table

| Action Name | User Name | Title                     | Signature Date    |
|-------------|-----------|---------------------------|-------------------|
| Approve     | Kasha Lim | Supplier Quality Engineer | 12-Jun-2026 12:50 |

\* Datetime fields are displayed according to the system time zone: (GMT-04:00) Eastern Daylight Time (America/New\_York)