INSTRUCTIONS FOR USE (IFU)

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MANUFACTURER

Manufacturer	Zima International, Inc. Also DBA Dandy
Address	1320 N. 300 W., Lehi, UT 84043, United States
Contact / Website	www.meetdandy.com
Complaints	complaints@meetdandy.com

DEVICE IDENTITY

Device Identity

Device Family	MouthGuard
GMDN / EMDN Code	38620
Product Code(s)	NIGHTGRDS, THERMONGH, THERMONGS, THERMONGHS, THERMONGA, THERMONGNTI, THERMOSG, THERMOHFA, THERMOHOS, THERMONGS, THERMOHSOS, THERMOHSFA, THERMONGA, THERMONGNTI, THERMOSFA
Device Type	Custom-made dental device (CMD) / Patient Matched Medical Device (PMMD)
Sterility	Not supplied sterile
Single-use / Reuse	Reusable
Intended User	Prescriber: Dental professional / dental practice

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AUSTRALIA (ENGLISH)**INSTRUCTIONS FOR USE (IFU)**

Australia — Patient Matched Medical Device (PMMD)

GMDN / EMDN	38620
Product Code(s)	NIGHTGRDS, THERMONGH, THERMONGS, THERMONGHS, THERMONGA, THERMONGNTI, THERMOSG, THERMOHFA, THERMOHOS, THERMONGS, THERMOHSOS, THERMOHSFA, THERMONGA, THERMONGNTI, THERMOSFA
Manufacturer	Zima International, Inc. Also DBA Dandy 1320 N. 300 W., Lehi, UT 84043 United States
Australian Sponsor	Name: DANDY Labs AUS Pty Ltd Company No.: ABN: 33693955761 Address: Tower One - International Towers, Level 46, 100 Barangaroo Avenue, Sydney, New South Wales 2000, Australia
Intended User	Prescriber — Dental professional / dental practice
Device Type	Patient Matched Medical Device (PMMD)
Sterility	Not supplied sterile
Single-use / Reusable	Reusable

1) Intended Purpose (Intended Use)

The device is made of biocompatible materials and is custom-made for each patient. The nightguard fits comfortably over the upper or lower teeth and provides protection against damage caused by bruxism.

2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: Hard, Soft, and Dual-Laminated Night Guards (Thermoformed) are dental devices that protect teeth from nighttime grinding and clenching.
- Material: Hard, Soft, Hard/Soft, Astron, and 3D Printed

[Erkoflex Thermoforming Discs - Soft](#)

[Erkoflex Thermoforming Discs - Hard](#)

[Erkoflex Thermoforming Discs - Hard/Soft](#)

[Clearsplint - Powder](#)

[Clearsplint - Liquid](#)

[Keysplint](#)

3) Contraindications

Patients with loose or damaged teeth or those with temporomandibular joint disorder (TMD) should not use the device without first consulting a dental professional.

4) Warnings / Precautions / Potential Risks

⚠ WARNING

The potential risks include mouth or tooth soreness, which should resolve after a few days of use.

⚠ CAUTION

NA

5) Cleaning and Care

Directions for Use

- Wash your hands before handling the nightguard.
- Place the nightguard over the upper or lower teeth as the dentist instructs.
- Bite down gently to ensure a comfortable fit.
- Remove the nightguard in the morning and rinse with warm water.

Cleaning instructions

- Dandy recommends cleaning your device daily before and after every use. Immediately after removal,
- Use warm water (not to exceed 45°C/113°F) and dishwashing soap to remove remaining saliva and debris from the device.

Only use soft, antibacterial dish soap (liquid)

- Do not use flavored or scented soaps
- Do not use soaps with heavy dyes and colors
- Rinse the device thoroughly with cool water.
- Do not use hot water, as it can warp the device.
- Do not use mouthwash, as it can discolor the device.

Clean the device daily.

USING CLEANING AGENTS

- For periodic debris removal and maintenance, use a soft bristle brush to clean the inside of your appliance for debris removal.
- The product should be cleaned with antibacterial cleaning solutions like Polident and Efferdent.
- These cleaners are designed for use in brief sessions of 10–15 minutes.
- Always make sure to follow the manufacturer's instructions for proper use.

Do not use:

- Avoid immersing the product in, or cleaning it with, any of the following substances: sodium hypochlorite (chlorine bleach), hydrogen peroxide, mouthwash, or alcohol products. Using these can lead to discoloration and damage to the material.

6) Storage

The device should be stored in a dry, clean container when not in use.

7) Expected Life / Service

.A mouthguard's life expectancy varies by type: 1-5 years for durable custom-fit guards, depending heavily on usage, material, and care, with heavy grinding or poor cleaning shortening their lifespan significantly and making them less protective. Regular checks for cracks, fit, odor, or damage are key indicators for replacement, ensuring consistent protection against sports impacts or teeth grinding (bruxism).

8) Incident / Complaint Reporting

Report suspected serious incidents associated with this device to:

- The Manufacturer and/or Dandy Australia entity using the contacts above, and
- Therapeutic Goods Administration (TGA)

In Australia, serious incidents involving medical devices must be reported to the TGA in accordance with the Therapeutic Goods Act 1989.

INFO

For full contact details of the local Authorized Representative, refer to: Authorized Representatives — end of this document.

9) Disposal

Specialized services handle biohazardous waste, but for home users, standard household disposal is key. This involves keeping it out of landfills whenever possible, and definitely out of sewers and recycling facilities.

EUROPEAN UNION (ENGLISH)**INSTRUCTIONS FOR USE (IFU)**

European Union — Custom-Made Device (CMD)

GMDN / EMDN	Q0199 - ODONTOLOGY DEVICES - OTHER
Product Code(s)	NIGHTGRDS, THERMONGH, THERMONGS, THERMONGHS, THERMONGA, THERMONGNTI, THERMOSG, THERMOHFA, THERMOHOS, THERMONGS, THERMOHSOS, THERMOHSFA, THERMONGA, THERMONGNTI, THERMOSFA
Manufacturer	Zima International, Inc. Also DBA Dandy 1320 N. 300 W., Lehi, UT 84043 United States
EU authorized representative (EURP)	Name: Dandy Labs Europe Company No.: Registration No.: 989857065 Address: 3, boulevard de Sebastopol, 75001 Paris, France
Intended User	Prescriber — Dental professional / dental practice
Device Type	Custom-made dental device (CMD)
Sterility	Not supplied sterile
Single-use / Reusable	Reusable

1) Intended Purpose (Intended Use)

The device is made of biocompatible materials and is custom-made for each patient. The nightguard fits comfortably over the upper or lower teeth and provides protection against damage caused by bruxism.

2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: Hard, Soft, and Dual-Laminated Night Guards (Thermoformed) are dental devices that protect teeth from nighttime grinding and clenching.
- Material: Hard, Soft, Hard/Soft, Astron, and 3D Printed

[Erkoflex Thermoforming Discs - Soft](#)

[Erkoflex Thermoforming Discs - Hard](#)

[Erkoflex Thermoforming Discs - Hard/Soft](#)

[Clearsplint - Powder](#)

[Clearsplint - Liquid](#)

[Keysplint](#)

3) Contraindications

Patients with loose or damaged teeth or those with temporomandibular joint disorder (TMD) should not use the device without first consulting a dental professional.

4) Warnings / Precautions / Potential Risks

⚠ WARNING

The potential risks include mouth or tooth soreness, which should resolve after a few days of use.

⚠ CAUTION

[Insert precaution text]

5) Cleaning and Care

Directions for Use

- Wash your hands before handling the nightguard.
- Place the nightguard over the upper or lower teeth as the dentist instructs.
- Bite down gently to ensure a comfortable fit.
- Remove the nightguard in the morning and rinse with warm water.

Cleaning instructions

- Dandy recommends cleaning your device daily before and after every use. Immediately after removal,
- Use warm water (not to exceed 45°C/113°F) and dishwashing soap to remove remaining saliva and debris from the device.

Only use soft, antibacterial dish soap (liquid)

- Do not use flavored or scented soaps
- Do not use soaps with heavy dyes and colors
- Rinse the device thoroughly with cool water.
- Do not use hot water, as it can warp the device.
- Do not use mouthwash, as it can discolor the device.

Clean the device daily.

USING CLEANING AGENTS

- For periodic debris removal and maintenance, use a soft bristle brush to clean the inside of your appliance for debris removal.
- The product should be cleaned with antibacterial cleaning solutions like Polident and Efferdent.
- These cleaners are designed for use in brief sessions of 10–15 minutes.
- Always make sure to follow the manufacturer's instructions for proper use.

Do not use:

- Avoid immersing the product in, or cleaning it with, any of the following substances: sodium hypochlorite (chlorine bleach), hydrogen peroxide, mouthwash, or alcohol products. Using these can lead to discoloration and damage to the material.

6) Storage

The device should be stored in a dry, clean container when not in use.

7) Expected Life / Service

A mouthguard's life expectancy varies by type: 1-5 years for durable custom-fit guards, depending heavily on usage, material, and care, with heavy grinding or poor cleaning shortening their lifespan significantly and making them less protective. Regular checks for cracks, fit, odor, or damage are key

indicators for replacement, ensuring consistent protection against sports impacts or teeth grinding (bruxism).

8) Incident / Complaint Reporting

Report suspected serious incidents associated with this device to:

- The Manufacturer and/or Dandy EU entity using the contacts above, and corresponding local requirements
- The Competent Authority of your Member State, in accordance with Regulation (EU) 2017/745 (MDR).

Competent Authorities:

- France: Agence Nationale de Sécurité du Médicament et des Produits de Santé (ANSM)
- Spain : Spanish Agency for Medicines and Health Products (AEMPS)
- Germany : Bundesinstitut für Arzneimittel und Medizinprodukte (BfArM) (BfArM)
- Italy : Ministero della Salute / Istituto Superiore di Sanità (ISS)
- Spain: Agencia Española de Medicamentos y Productos Sanitarios (AEMPS) for Spain
- Netherlands: Inspectie Gezondheidszorg en Jeugd (IGJ)
- Austria: Bundesamt für Sicherheit im Gesundheitswesen (BASG)
- Sweden: Läkemedelsverket
- Denmark: Lægemiddelstyrelsen
- Portugal: Autoridade Nacional do Medicamento e Produtos de Saúde (INFARMED)

INFO

For full contact details of the local Authorized Representative, refer to: Authorized Representatives — end of this document.

9) Disposal

Specialized services handle biohazardous waste, but for home users, standard household disposal is key. This involves keeping it out of landfills whenever possible, and definitely out of sewers and recycling facilities.

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UNITED KINGDOM (ENGLISH)**INSTRUCTIONS FOR USE (IFU)**

United Kingdom — Custom-Made Device (CMD)

GMDN / EMDN	38620
Product Code(s)	NIGHTGRDS, THERMONGH, THERMONGS, THERMONGHS, THERMONGA, THERMONGNTI, THERMOSG, THERMOHFA, THERMOHOS, THERMONGS, THERMOHSOS, THERMOHSFA, THERMONGA, THERMONGNTI, THERMOSFA
Manufacturer	Zima International, Inc. Also DBA Dandy 1320 N. 300 W., Lehi, UT 84043 United States
UK Responsible Person (UKRP)	Name: DANDY LABS GB, LTD Company No.: Company Number: 16873608 Address: 5 New Street Square, London EC4A 3TW, United Kingdom
Intended User	Prescriber — Dental professional / dental practice
Device Type	Custom-made dental device (CMD)
Sterility	Not supplied sterile
Single-use / Reusable	Reusable

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7) Expected Life / Service

A mouthguard's life expectancy varies by type: 1-5 years for durable custom-fit guards, depending heavily on usage, material, and care, with heavy grinding or poor cleaning shortening their lifespan significantly and making them less protective. Regular checks for cracks, fit, odor, or damage are key indicators for replacement, ensuring consistent protection against sports impacts or teeth grinding (bruxism).

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- The Manufacturer and/or UK Responsible Person (UKRP) using the contacts above, and
- The UK Medicines and Healthcare products Regulatory Agency (MHRA) — www.gov.uk/mhra

In accordance with The Medical Devices Regulations 2002 (as amended).

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The table below lists the local Authorized Representative or Sponsor for each country / jurisdiction where this device is marketed. For incident reporting, contact the representative for your country.

La tabla siguiente enumera el Representante Autorizado o Patrocinador local para cada país. | Le tableau ci-dessous liste le Représentant Autorisé local pour chaque pays. | Die folgende Tabelle enthält die bevollmächtigten Vertreter je Land.

Jurisdiction	Authorised Rep / Sponsor	Registered Address / Contact Email	Reg. Number	Regulatory Body
Australia (AUS)	DANDY Labs AUS Pty Ltd	Tower One, Level 46, 100 Barangaroo Ave, Sydney NSW 2000 Compliance@meetdandy.com	ABN: 33693955761	TGA (Therapeutic Goods Administration)
EU — English	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, France Compliance@meetdandy.com	Reg: 989857065	Competent Authority / EUDAMED
EU — España	Dandy Labs Europa	3, boulevard de Sebastopol, 75001 París, Francia Compliance@meetdandy.com	Reg: 989857065	AEMPS / Autoridad Competente
EU — France	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, France Compliance@meetdandy.com	Reg: 989857065	ANSM / Autorité Compétente
EU — Deutschland	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, Frankreich Compliance@meetdandy.com	Reg: 989857065	BfArM / Zuständige Behörde
EU — Italia	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Parigi, Francia Compliance@meetdandy.com	Reg: 989857065	Min. della Salute / ISS
United Kingdom (UK)	DANDY LABS GB, LTD	5 New Street Square, London EC4A 3TW Compliance@meetdandy.com	Co. No.: 16873608	MHRA

NOTE

Contact details are subject to change. Always verify current information on the official product labelling or the manufacturer's website before reporting an incident. / Les coordonnées sont susceptibles de changer. / Angaben können sich ändern. / I recapiti sono soggetti a modifiche.

**Manufacturer / Hersteller /
Fabbicante / Fabricant /
Fabricante**

Zima International, Inc. Also DBA Dandy
1320 N. 300 W., Lehi, UT 84043, United States
compliance@meetdandy.com

This is a representation of an electronic record that was signed electronically and this page is the manifestation of the electronic signatures.

Signatory Table

Action Name	User Name	Title	Signature Date
Approve	Kasha Lim	Supplier Quality Engineer	12-Jun-2026 01:55

* Datetime fields are displayed according to the system time zone: (GMT-04:00) Eastern Daylight Time (America/New_York)