INSTRUCTIONS FOR USE (IFU)

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MANUFACTURER

Manufacturer	Zima International, Inc. Also DBA Dandy
Address	1320 N. 300 W., Lehi, UT 84043, United States
Contact / Website	www.meetdandy.com
Complaints	complaints@meetdandy.com

DEVICE IDENTITY

Device Identity

Device Family	Implant Abutment
GMDN / EMDN Code	63349
Product Code(s)	TICLOCDES, TICLOCMED, TICLOCMEG, TICLOCNOB, TICLOCSTR, TICLOCZIM
Device Type	Custom-made dental device (CMD) / Patient Matched Medical Device (PMMD)
Sterility	Not supplied sterile
Single-use / Reuse	Single-use
Intended User	Prescriber: Dental professional / dental practice

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AUSTRALIA (ENGLISH)**INSTRUCTIONS FOR USE (IFU)**

Australia — Patient Matched Medical Device (PMMD)

GMDN / EMDN	63349
Product Code(s)	TICLOCDES, TICLOCMED, TICLOCMEG, TICLOCNOB, TICLOCSTR, TICLOCZIM
Manufacturer	Zima International, Inc. Also DBA Dandy 1320 N. 300 W., Lehi, UT 84043 United States
Australian Sponsor	Name: DANDY Labs AUS Pty Ltd Company No.: ABN: 33693955761 Address: Tower One - International Towers, Level 46, 100 Barangaroo Avenue, Sydney, New South Wales 2000, Australia
Intended User	Prescriber — Dental professional / dental practice
Device Type	Patient Matched Medical Device (PMMD)
Sterility	Not supplied sterile
Single-use / Reusable	Single use

1) Intended Purpose (Intended Use)

These dental abutments are connected directly to an endosseous implant for use as an aid in prosthetic rehabilitation. The abutment is designed with a scalloped cuff that profiles the natural contours of soft tissue and is custom-made to meet your patient's specific needs. Made from titanium, titanium alloy, or Zirconia. The abutment gets restored with a crown, bridge, or is used to anchor a restorative device.

2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: Dental abutments are designed to be used in the Maxilla or mandible to support a prosthesis on endosseous implants, thereby restoring chewing function for the patient.
- Material: Titanium/Zirconia

[Zirconia](#)**3) Contraindications****Do not use in patients:**

- Those who are medically unfit for dental implant procedures.
- Those who are allergic or have hypersensitivity to pure titanium or titanium alloy.

- Those where adequate numbers of implants could not be placed to achieve full functional support for a prosthesis.

4) Warnings / Precautions / Potential Risks

THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING

WARNING

- For safe and effective use of dental implants, it is strongly suggested that specialised training be undertaken, including hands-on training to learn proper technique, biomechanical requirements, and radiographic evaluations.
- Responsibility for proper patient selection, adequate training, experience in implant placement, and providing appropriate information for informed consent rests with the practitioner. Improper technique can result in implant failure, damage to nerves/vessels, and/or loss of supporting bone.

CAUTION

- New and experienced Implant users should do training before using a new system or attempting a new treatment method.
- Take special care when treating patients who have local or systemic factors that could affect the healing of the bone and soft tissue. (e.g., poor oral hygiene, uncontrolled diabetes, are on steroid therapy, smokers, infection in the nearby bone, and patients who had oro-facial radiotherapy.)
- Remove the temporary restoration if applicable.
- Verify and re-tighten the screw if applicable. Seal the screw access hole with silicone and cement the final crown using conventional procedures.
- Make sure that all excess cement is removed from the margin.
- Overtightening of the abutment may lead to screw fracture)
- Cement the final crown or framework using conventional procedures after sealing the access hole, ensuring there's no excess cement on the margin.
- For screw retained restorations
- The final restoration (crown and abutment) is attached to the implant
- Verify the seating
- Tighten the screw to the required torque. (Do not exceed the recommended tightening torque for the abutment screw.

5) Cleaning and Care

These abutments are not supplied sterile and are intended for single use

6) Storage

The product must be stored in a dry, well-ventilated area at room temperature and protected from direct sunlight.

Incorrect storage may influence device characteristics

7) Expected Life / Service

A dental abutment typically lasts 10 to 20+ years, often aligning with the crown's lifespan, but can last longer with excellent care, it usually outlasts the crown (10-15 yrs) but is less durable than the implant

Document State: Effective (Simple)

Page 5 of 30

post. Longevity depends on material, chewing forces (molars wear faster), and habits like grinding (bruxism). Still, good hygiene and regular checkups are key to maximizing its life, with some studies showing great success with all-ceramic options.

8) Incident / Complaint Reporting

Report suspected serious incidents associated with this device to:

- The Manufacturer and/or Dandy Australia entity using the contacts above, and
- Therapeutic Goods Administration (TGA)

In Australia, serious incidents involving medical devices must be reported to the TGA in accordance with the Therapeutic Goods Act 1989.

INFO

For full contact details of the local Authorized Representative, refer to: Authorized Representatives — end of this document.

9) Disposal

Must be disposed of as biological or medical waste if they have been used in a patient's mouth
Unused abutments (Titanium or Zirconia) may be disposed of as general waste

Guidelines:

Used Abutments:

Any dental restoration or component that has been in the patient's mouth is considered biological waste and must be handled and disposed of according to strict medical waste management regulations.

Unused Abutments:

If the custom abutment is unused and has not been sterilized or come into contact with bodily fluids, it may not fall under the same strict biohazard rules.

Material:

Titanium/Zirconia: The materials themselves (titanium alloy, zirconia) are not inherently hazardous materials in their solid form, but their classification as waste is determined by their potential biohazard contamination during use or the preparation process.

Recycling:

While some materials are recyclable in general construction/demolition waste streams, dental practices must use certified medical waste streams to avoid regulatory non-compliance.

Follow Local Regulations: Disposal rules are governed by federal, state, and local health and environmental regulations. Ensure your practice is compliant with the specific requirements for your location.

EUROPEAN UNION (ENGLISH)**INSTRUCTIONS FOR USE (IFU)**

European Union — Custom-Made Device (CMD)

GMDN / EMDN	P01020180: Dental Implants - Accessories
Product Code(s)	TICLOCDES, TICLOCMED, TICLOCMEG, TICLOCNOB, TICLOCSTR, TICLOCZIM
Manufacturer	Zima International, Inc. Also DBA Dandy 1320 N. 300 W., Lehi, UT 84043 United States
EU authorized representative (EURP)	Name: Dandy Labs Europe Company No.: Registration No.: 989857065 Address: 3, boulevard de Sebastopol, 75001 Paris, France
Intended User	Prescriber — Dental professional / dental practice
Device Type	Custom-made dental device (CMD)
Sterility	Not supplied sterile
Single-use / Reusable	Single use

1) Intended Purpose (Intended Use)

These dental abutments are connected directly to an endosseous implant for use as an aid in prosthetic rehabilitation. The abutment is designed with a scalloped cuff that profiles the natural contours of soft tissue and is custom-made to meet your patient's specific needs. Made from titanium, titanium alloy, or Zirconia. The abutment gets restored with a crown, bridge, or is used to anchor a restorative device.

2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: Dental abutments are designed to be used in the Maxilla or mandible to support a prosthesis on endosseous implants, thereby restoring chewing function for the patient.
- Material: Titanium/Zirconia

[Zirconia](#)**3) Contraindications****Do not use in patients:**

- Those who are medically unfit for dental implant procedures.
- Those who are allergic or have hypersensitivity to pure titanium or titanium alloy.

Document State: Effective (Simple)

Page 7 of 30

Confidential

FRM-0002002, revision 2

- Those where adequate numbers of implants could not be placed to achieve full functional support for a prosthesis.

4) Warnings / Precautions / Potential Risks

THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING

WARNING

- For safe and effective use of dental implants, it is strongly suggested that specialised training be undertaken, including hands-on training to learn proper technique, biomechanical requirements, and radiographic evaluations.
- Responsibility for proper patient selection, adequate training, experience in implant placement, and providing appropriate information for informed consent rests with the practitioner. Improper technique can result in implant failure, damage to nerves/vessels, and/or loss of supporting bone.

CAUTION

- New and experienced Implant users should do training before using a new system or attempting a new treatment method.
- Take special care when treating patients who have local or systemic factors that could affect the healing of the bone and soft tissue. (e.g., poor oral hygiene, uncontrolled diabetes, are on steroid therapy, smokers, infection in the nearby bone, and patients who had oro -facial radiotherapy.)
- Remove the temporary restoration if applicable.
- Verify and re-tighten the screw if applicable. Seal the screw access hole with silicone and cement the final "closed" crown using conventional procedures.
- Make sure that all excess cement is removed from the margin.
- Overtightening of the abutment may lead to screw fracture)
- Cement the final crown or framework using conventional procedures after sealing the access hole, ensuring there's no excess cement on the margin.
- For screw retained restorations
- The final restoration (crown and abutment) is attached to the implant
- Verify the seating
- Tighten the screw to the required torque. (Do not exceed the recommended tightening torque for the abutment screw.

5) Cleaning and Care

These abutments are not supplied sterile and are intended for single use

6) Storage

The product must be stored in a dry, well-ventilated area at room temperature and protected from direct sunlight.

Incorrect storage may influence device characteristics

7) Expected Life / Service

A dental abutment typically lasts 10 to 20+ years, often aligning with the crown's lifespan, but can last longer with excellent care, it usually outlasts the crown (10-15 yrs) but is less durable than the implant

Document State: Effective (Simple)

Page 8 of 30

post. Longevity depends on material, chewing forces (molars wear faster), and habits like grinding (bruxism). Still, good hygiene and regular checkups are key to maximizing its life, with some studies showing great success with all-ceramic options.

8) Incident / Complaint Reporting

Report suspected serious incidents associated with this device to:

- The Manufacturer and/or Dandy EU entity using the contacts above, and corresponding local requirements
- The Competent Authority of your Member State, in accordance with Regulation (EU) 2017/745 (MDR).

Competent Authorities:

- France: Agence Nationale de Sécurité du Médicament et des Produits de Santé (ANSM)
- Spain : Spanish Agency for Medicines and Health Products (AEMPS)
- Germany : Bundesinstitut für Arzneimittel und Medizinprodukte (BfArM) (BfArM)
- Italy : Ministero della Salute / Istituto Superiore di Sanità (ISS)
- Spain: Agencia Española de Medicamentos y Productos Sanitarios (AEMPS) for Spain
- Netherlands: Inspectie Gezondheidszorg en Jeugd (IGJ)
- Austria: Bundesamt für Sicherheit im Gesundheitswesen (BASG)
- Sweden: Läkemedelsverket
- Denmark: Lægemiddelstyrelsen
- Portugal: Autoridade Nacional do Medicamento e Produtos de Saúde (INFARMED)

INFO

For the full contact information of the local Authorized Representative, see: Authorized Representatives — at the end of the document.

9) Disposal

Must be disposed of as biological or medical waste if they have been used in a patient's mouth

Unused abutments (Titanium or Zirconia) may be disposed of as general waste

Guidelines:

Used Abutments:

Any dental restoration or component that has been in the patient's mouth is considered biological waste and must be handled and disposed of according to strict medical waste management regulations.

Unused Abutments:

If the custom abutment is unused and has not been sterilized or come into contact with bodily fluids, it may not fall under the same strict biohazard rules.

Material:

Titanium/Zirconia: The materials themselves (titanium alloy, zirconia) are not inherently hazardous materials in their solid form, but their classification as waste is determined by their potential biohazard contamination during use or the preparation process.

Recycling:

While some materials are recyclable in general construction/demolition waste streams, dental practices must use certified medical waste streams to avoid regulatory non-compliance.

Follow Local Regulations: Disposal rules are governed by federal, state, and local health and environmental regulations. Ensure your practice is compliant with the specific requirements for your location.

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UNITED KINGDOM (ENGLISH)**INSTRUCTIONS FOR USE (IFU)**

United Kingdom — Custom-Made Device (CMD)

GMDN / EMDN	63349
Product Code(s)	TICLOCDES, TICLOCMED, TICLOCMEG, TICLOCNOB, TICLOCSTR, TICLOCZIM
Manufacturer	Zima International, Inc. Also DBA Dandy 1320 N. 300 W., Lehi, UT 84043 United States
UK Responsible Person (UKRP)	Name: DANDY LABS GB, LTD Company No.: Company Number: 16873608 Address: 5 New Street Square, London EC4A 3TW, United Kingdom
Intended User	Prescriber — Dental professional / dental practice
Device Type	Custom-made dental device (CMD)
Sterility	Not supplied sterile
Single-use / Reusable	Single use

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CAUTION

- New and experienced Implant users should do training before using a new system or attempting a new treatment method.
- Take special care when treating patients who have local or systemic factors that could affect the healing of the bone and soft tissue. (e.g., poor oral hygiene, uncontrolled diabetes, are on steroid therapy, smokers, infection in the nearby bone, and patients who had oro-facial radiotherapy.)
- Remove the temporary restoration if applicable.
- Verify and re-tighten the screw if applicable. Seal the screw access hole with silicone and cement the final "closed" crown using conventional procedures.
- Make sure that all excess cement is removed from the margin.
- Overtightening of the abutment may lead to screw fracture)
- Cement the final crown or framework using conventional procedures after sealing the access hole, ensuring there's no excess cement on the margin.
- For screw retained restorations
- The final restoration (crown and abutment) is attached to the implant
- Verify the seating
- Tighten the screw to the required torque. (Do not exceed the recommended tightening torque for the abutment screw.

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7) Expected Life / Service

A dental abutment typically lasts 10 to 20+ years, often aligning with the crown's lifespan, but can last longer with excellent care; it usually outlasts the crown (10-15 yrs) but is less durable than the implant post. Longevity depends on material, chewing forces (molars wear faster), and habits like grinding (bruxism). Still, good hygiene and regular checkups are key to maximizing its life, with some studies showing great success with all-ceramic options.

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- The Manufacturer and/or UK Responsible Person (UKRP) using the contacts above, and
- The UK Medicines and Healthcare products Regulatory Agency (MHRA) — www.gov.uk/mhra

In accordance with The Medical Devices Regulations 2002 (as amended).

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The table below lists the local Authorized Representative or Sponsor for each country / jurisdiction where this device is marketed. For incident reporting, contact the representative for your country.

La tabla siguiente enumera el Representante Autorizado o Patrocinador local para cada país. | Le tableau ci-dessous liste le Représentant Autorisé local pour chaque pays. | Die folgende Tabelle enthält die bevollmächtigten Vertreter je Land.

Jurisdiction	Authorised Rep / Sponsor	Registered Address / Contact Email	Reg. Number	Regulatory Body
Australia (AUS)	DANDY Labs AUS Pty Ltd	Tower One, Level 46, 100 Barangaroo Ave, Sydney NSW 2000 Compliance@meetdandy.com	ABN: 33693955761	TGA (Therapeutic Goods Administration)
EU — English	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, France Compliance@meetdandy.com	Reg: 989857065	Competent Authority / EUDAMED
EU — España	Dandy Labs Europa	3, boulevard de Sebastopol, 75001 París, Francia Compliance@meetdandy.com	Reg: 989857065	AEMPS / Autoridad Competente
EU — France	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, France Compliance@meetdandy.com	Reg: 989857065	ANSM / Autorité Compétente
EU — Deutschland	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, Frankreich Compliance@meetdandy.com	Reg: 989857065	BfArM / Zuständige Behörde
EU — Italia	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Parigi, Francia Compliance@meetdandy.com	Reg: 989857065	Min. della Salute / ISS
United Kingdom (UK)	DANDY LABS GB, LTD	5 New Street Square, London EC4A 3TW Compliance@meetdandy.com	Co. No.: 16873608	MHRA

NOTE

Contact details are subject to change. Always verify current information on the official product labelling or the manufacturer's website before reporting an incident. / Les coordonnées sont susceptibles de changer. / Angaben können sich ändern. / I recapiti sono soggetti a modifiche.

**Manufacturer / Hersteller /
Fabbricante / Fabricant /
Fabricante**

Zima International, Inc. Also DBA Dandy
1320 N. 300 W., Lehi, UT 84043, United States
compliance@meetdandy.com

This is a representation of an electronic record that was signed electronically and this page is the manifestation of the electronic signatures.

Signatory Table

Action Name	User Name	Title	Signature Date
Approve	Kasha Lim	Supplier Quality Engineer	12-Jun-2026 01:47

* Datetime fields are displayed according to the system time zone: (GMT-04:00) Eastern Daylight Time (America/New_York)