INSTRUCTIONS FOR USE (IFU)

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MANUFACTURER

Manufacturer	Zima International, Inc. Also DBA Dandy
Address	1320 N. 300 W., Lehi, UT 84043, United States
Contact / Website	www.meetdandy.com
Complaints	complaints@meetdandy.com

DEVICE IDENTITY

Device Identity

Device Family	Dental Guided Surgery Reference Plate
GMDN / EMDN Code	57949
Product Code(s)	SSURGG
Device Type	Custom-made dental device (CMD) / Patient Matched Medical Device (PMMD)
Sterility	Not supplied sterile
Single-use / Reuse	Single-use
Intended User	Prescriber: Dental professional / dental practice

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AUSTRALIA (ENGLISH)**INSTRUCTIONS FOR USE (IFU)**

Australia — Patient Matched Medical Device (PMMD)

GMDN / EMDN	57949
Product Code(s)	SSURGG
Manufacturer	Zima International, Inc. Also DBA Dandy 1320 N. 300 W., Lehi, UT 84043 United States
Australian Sponsor	Name: DANDY Labs AUS Pty Ltd Company No.: ABN: 33693955761 Address: Tower One - International Towers, Level 46, 100 Barangaroo Avenue, Sydney, New South Wales 2000, Australia
Intended User	Prescriber — Dental professional / dental practice
Device Type	Patient Matched Medical Device (PMMD)
Sterility	Not supplied sterile
Single-use / Reusable	Single-use

1) Intended Purpose (Intended Use)

Dandy's surgical guide is a disposable medical device customized (3D-printed) for a single patient. The 3D-printed medical device requires steam sterilization before use in surgery.

2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: It is indicated for intraoral use in guided dental implant surgery, performed by a dental professional.
- Material: Keyguide

[KeyGuide](#)

3) Contraindications

Dental surgical guide contraindications aren't about the guide itself, but the underlying implant procedure; they involve uncontrolled systemic diseases (diabetes, immune issues, clotting disorders, cancer), heavy smoking, insufficient bone, certain medications (bisphosphonates, blood thinners), pregnancy, uncontrolled bruxism, or conditions preventing good oral hygiene, all impacting healing and implant success, with older age (not old) and fully grown jaws being key factors for success

4) Warnings / Precautions / Potential Risks

WARNING

Potential inaccuracies may arise from poor scans or planning, guide fit issues (such as fractures/poor fit), technology dependence, limited surgical adjustments, and general implant risks, including infection, nerve damage, or sinus perforation. However, guides aim to reduce common surgical errors by improving precision. Key risks stem from the process (planning, imaging) and the procedure itself (implant failure, healing issues).

CAUTION

NA

5) Cleaning and Care

CLEANING

- Thoroughly remove any visible soil remaining on the surface. Alcohol can be used to wipe the surface of the surgical guide to ensure a clean part results from this process. No contamination
- shall be found by the naked eye or felt by running fingers on the surface of the part. Special attention
- should be given to edges, hollows, and other geometries.

DRYING

- Dry the part, especially in the hollows and geometries. Compressed air can be used for air drying.

INSPECTION

- Once dried, inspect the medical device for damage or cracks. Replace the part if needed.

DIRECTIONS FOR USE

- To ensure the proper use of the device, dental professionals should follow these steps:
- Wash your hands thoroughly before handling the device.
- Carefully insert the device into the patient's teeth, starting from the back of the patient's mouth and working forward.
- Use your fingers to gently press the device onto the patient's teeth to ensure a snug fit.

6) Storage

The products should be stored in their packaging in a dry place at a normal temperature (18–25°C/64–77°F). Use the sterilized components within the stated time period from the sterile bag.

7) Expected Life / Service

Dental surgical guide itself is a single-use or short-term tool for precise implant placement

8) Incident / Complaint Reporting

Report suspected serious incidents associated with this device to:

- The Manufacturer and/or Dandy Australia entity using the contacts above, and
- Therapeutic Goods Administration (TGA)

In Australia, serious incidents involving medical devices must be reported to the TGA in accordance with the Therapeutic Goods Act 1989.

INFO

For full contact details of the local Authorized Representative, refer to: Authorized Representatives — end of this document.

9) Disposal

Dental surgical guide disposal depends on contamination: if blood/saliva-soaked (Regulated Medical Waste), use red biohazard bags/containers for licensed pickup; if clean, it is generally considered regular clinic waste (check local rules). For sharps (drills, wires), use puncture-proof sharps containers. Always follow federal, state, and local regulations and use a licensed medical waste hauler for hazardous/infectious waste.

EUROPEAN UNION (ENGLISH)**INSTRUCTIONS FOR USE (IFU)**

European Union — Custom-Made Device (CMD)

GMDN / EMDN	Q010399- Surgical Dental Devices - Other
Product Code(s)	SSURGG
Manufacturer	Zima International, Inc. Also DBA Dandy 1320 N. 300 W., Lehi, UT 84043 United States
EU authorized representative (EURP)	Name: Dandy Labs Europe Company No.: Registration No.: 989857065 Address: 3, boulevard de Sebastopol, 75001 Paris, France
Intended User	Prescriber — Dental professional / dental practice
Device Type	Custom-made dental device (CMD)
Sterility	Not supplied sterile
Single-use / Reusable	Single-use

1) Intended Purpose (Intended Use)

Dandy's surgical guide is a disposable medical device customized (3D-printed) for a single patient. The 3D-printed medical device requires steam sterilization before use in surgery.

2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: It is indicated for intraoral use in guided dental implant surgery, performed by a dental professional.
- Material: Keyguide

[KeyGuide](#)**3) Contraindications**

Dental surgical guide contraindications aren't about the guide itself, but the underlying implant procedure; they involve uncontrolled systemic diseases (diabetes, immune issues, clotting disorders, cancer), heavy smoking, insufficient bone, certain medications (bisphosphonates, blood thinners), pregnancy, uncontrolled bruxism, or conditions preventing good oral hygiene, all impacting healing and implant success, with older age (not old) and fully grown jaws being key factors for success

4) Warnings / Precautions / Potential Risks

WARNING

Potential inaccuracies may arise from poor scans or planning, guide fit issues (such as fractures/poor fit), technology dependence, limited surgical adjustments, and general implant risks, including infection, nerve damage, or sinus perforation. However, guides aim to reduce common surgical errors by improving precision. Key risks stem from the process (planning, imaging) and the procedure itself (implant failure, healing issues).

CAUTION

NA

5) Cleaning and Care

CLEANING

- Thoroughly remove any visible soil remaining on the surface. Alcohol can be used to wipe the surface of the surgical guide to ensure a clean part results from this process. No contamination
- shall be found by the naked eye or felt by running fingers on the surface of the part. Special attention
- should be given to edges, hollows, and other geometries.

DRYING

- Dry the part, especially in the hollows and geometries. Compressed air can be used for air drying.

INSPECTION

- Once dried, inspect the medical device for damage or cracks. Replace the part if needed.

DIRECTIONS FOR USE

- To ensure the proper use of the device, dental professionals should follow these steps:
- Wash your hands thoroughly before handling the device.
- Carefully insert the device into the patient's teeth, starting from the back of the patient's mouth and working forward.
- Use your fingers to gently press the device onto the patient's teeth to ensure a snug fit.

6) Storage

The products should be stored in their packaging in a dry place at a normal temperature (18–25°C/64–77°F). Use the sterilized components within the stated time period from the sterile bag.

7) Expected Life / Service

Dental surgical guide itself is a single-use or short-term tool for precise implant placement

8) Incident / Complaint Reporting

Report suspected serious incidents associated with this device to:

- The Manufacturer and/or Dandy EU entity using the contacts above, and corresponding local requirements
- The Competent Authority of your Member State, in accordance with Regulation (EU) 2017/745 (MDR).

Competent Authorities:

- France: Agence Nationale de Sécurité du Médicament et des Produits de Santé (ANSM)
- Spain : Spanish Agency for Medicines and Health Products (AEMPS)
- Germany : Bundesinstitut für Arzneimittel und Medizinprodukte (BfArM) (BfArM)
- Italy : Ministero della Salute / Istituto Superiore di Sanità (ISS)
- Spain: Agencia Española de Medicamentos y Productos Sanitarios (AEMPS) for Spain
- Netherlands: Inspectie Gezondheidszorg en Jeugd (IGJ)
- Austria: Bundesamt für Sicherheit im Gesundheitswesen (BASG)
- Sweden: Läkemedelsverket
- Denmark: Lægemiddelstyrelsen
- Portugal: Autoridade Nacional do Medicamento e Produtos de Saúde (INFARMED)

INFO

For the full contact information of the local Authorized Representative, see: Authorized Representatives — at the end of the document.

9) Disposal

Dental surgical guide disposal depends on contamination: if blood/saliva-soaked (Regulated Medical Waste), use red biohazard bags/containers for licensed pickup; if clean, it is generally considered regular clinic waste (check local rules). For sharps (drills, wires), use puncture-proof sharps containers. Always follow federal, state, and local regulations and use a licensed medical waste hauler for hazardous/infectious waste.

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UNITED KINGDOM (ENGLISH)**INSTRUCTIONS FOR USE (IFU)**

United Kingdom — Custom-Made Device (CMD)

GMDN / EMDN	57949
Product Code(s)	SSURGG
Manufacturer	Zima International, Inc. Also DBA Dandy 1320 N. 300 W., Lehi, UT 84043 United States
UK Responsible Person (UKRP)	Name: DANDY LABS GB, LTD Company No.: Company Number: 16873608 Address: 5 New Street Square, London EC4A 3TW, United Kingdom
Intended User	Prescriber — Dental professional / dental practice
Device Type	Custom-made dental device (CMD)
Sterility	Not supplied sterile (if applicable)
Single-use / Reusable	Single use

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5) Cleaning and Care

CLEANING

- Thoroughly remove any visible soil remaining on the surface. Alcohol can be used to wipe the surface of the surgical guide to ensure a clean part results from this process. No contamination shall be found by the naked eye or felt by running fingers on the surface of the part. Special attention should be given to edges, hollows, and other geometries.

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- The Manufacturer and/or UK Responsible Person (UKRP) using the contacts above, and
- The UK Medicines and Healthcare products Regulatory Agency (MHRA) — www.gov.uk/mhra

In accordance with The Medical Devices Regulations 2002 (as amended).

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The table below lists the local Authorized Representative or Sponsor for each country / jurisdiction where this device is marketed. For incident reporting, contact the representative for your country.

La tabla siguiente enumera el Representante Autorizado o Patrocinador local para cada país. | Le tableau ci-dessous liste le Représentant Autorisé local pour chaque pays. | Die folgende Tabelle enthält die bevollmächtigten Vertreter je Land.

Jurisdiction	Authorised Rep / Sponsor	Registered Address / Contact Email	Reg. Number	Regulatory Body
Australia (AUS)	DANDY Labs AUS Pty Ltd	Tower One, Level 46, 100 Barangaroo Ave, Sydney NSW 2000 Compliance@meetdandy.com	ABN: 33693955761	TGA (Therapeutic Goods Administration)
EU — English	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, France Compliance@meetdandy.com	Reg: 989857065	Competent Authority / EUDAMED
EU — España	Dandy Labs Europa	3, boulevard de Sebastopol, 75001 París, Francia Compliance@meetdandy.com	Reg: 989857065	AEMPS / Autoridad Competente
EU — France	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, France Compliance@meetdandy.com	Reg: 989857065	ANSM / Autorité Compétente
EU — Deutschland	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, Frankreich Compliance@meetdandy.com	Reg: 989857065	BfArM / Zuständige Behörde
EU — Italia	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Parigi, Francia Compliance@meetdandy.com	Reg: 989857065	Min. della Salute / ISS
United Kingdom (UK)	DANDY LABS GB, LTD	5 New Street Square, London EC4A 3TW Compliance@meetdandy.com	Co. No.: 16873608	MHRA

NOTE

Contact details are subject to change. Always verify current information on the official product labelling or the manufacturer's website before reporting an incident. / Les coordonnées sont susceptibles de changer. / Angaben können sich ändern. / I recapiti sono soggetti a modifiche.

**Manufacturer / Hersteller /
Fabbicante / Fabricant /
Fabricante**

Zima International, Inc. Also DBA Dandy
1320 N. 300 W., Lehi, UT 84043, United States
compliance@meetdandy.com

This is a representation of an electronic record that was signed electronically and this page is the manifestation of the electronic signatures.

Signatory Table

Action Name	User Name	Title	Signature Date
Approve	Kasha Lim	Supplier Quality Engineer	12-Jun-2026 01:36

* Datetime fields are displayed according to the system time zone: (GMT-04:00) Eastern Daylight Time (America/New_York)