# INSTRUCTIONS FOR USE (IFU)

Dandy / Zima International, Inc. | AUS · NZ · UK · EU/EN · EU/ES · EU/FR · EU/DE · EU/IT · EU/DA ·  
EU/NL · EU/PL · EU/SV · EU/NB · EU/PT

## MANUFACTURER

|                          |                                               |
|--------------------------|-----------------------------------------------|
| <b>Manufacturer</b>      | Zima International, Inc. Also DBA Dandy       |
| <b>Address</b>           | 1320 N. 300 W., Lehi, UT 84043, United States |
| <b>Contact / Website</b> | www.meetdandy.com                             |
| <b>Complaints</b>        | complaints@meetdandy.com                      |

## DEVICE IDENTITY

### Device Identity

|                           |                                                                         |
|---------------------------|-------------------------------------------------------------------------|
| <b>Device Family</b>      | Bite Registration Rim Wax                                               |
| <b>GMDN / EMDN Code</b>   | 38602                                                                   |
| <b>Product Code(s)</b>    | WAXRIM, SWAXPH                                                          |
| <b>Device Type</b>        | Custom-made dental device (CMD) / Patient Matched Medical Device (PMMD) |
| <b>Sterility</b>          | Not supplied sterile                                                    |
| <b>Single-use / Reuse</b> | Single-use                                                              |
| <b>Intended User</b>      | Prescriber: Dental professional / dental practice                       |

## TABLE OF CONTENT

|                                                         |           |
|---------------------------------------------------------|-----------|
| <b>AUSTRALIA (ENGLISH)</b>                              | <b>6</b>  |
| 1) Intended Purpose (Intended Use)                      | 6         |
| 2) Device Description / Key Specifications              | 6         |
| 3) Contraindications                                    | 6         |
| 4) Warnings / Precautions / Potential Risks             | 7         |
| 5) Cleaning and Care                                    | 7         |
| 6) Storage                                              | 7         |
| 7) Expected Life / Service                              | 8         |
| 8) Incident / Complaint Reporting                       | 8         |
| 9) Disposal                                             | 8         |
| <b>NEW ZEALAND (ENGLISH)</b>                            | <b>9</b>  |
| 1) Intended Purpose (Intended Use)                      | 9         |
| 2) Device Description / Key Specifications              | 9         |
| 3) Contraindications                                    | 9         |
| 4) Warnings / Precautions / Potential Risks             | 10        |
| 5) Cleaning and Care                                    | 10        |
| 6) Storage                                              | 10        |
| 7) Expected Life / Service                              | 11        |
| 8) Incident / Complaint Reporting                       | 11        |
| 9) Disposal                                             | 11        |
| <b>UNITED KINGDOM (ENGLISH)</b>                         | <b>12</b> |
| 1) Intended Purpose (Intended Use)                      | 12        |
| 2) Device Description / Key Specifications              | 12        |
| 3) Contraindications                                    | 12        |
| 4) Warnings / Precautions / Potential Risks             | 13        |
| 5) Cleaning and Care                                    | 13        |
| 6) Storage                                              | 13        |
| 7) Expected Life / Service                              | 14        |
| 8) Incident / Complaint Reporting                       | 14        |
| 9) Disposal                                             | 14        |
| <b>EUROPEAN UNION (ENGLISH)</b>                         | <b>15</b> |
| 1) Intended Purpose (Intended Use)                      | 15        |
| 2) Device Description / Key Specifications              | 15        |
| 3) Contraindications                                    | 15        |
| 4) Warnings / Precautions / Potential Risks             | 16        |
| 5) Cleaning and Care                                    | 16        |
| 6) Storage                                              | 16        |
| 7) Expected Life / Service                              | 17        |
| 8) Incident / Complaint Reporting                       | 17        |
| 9) Disposal                                             | 17        |
| <b>UNIÓN EUROPEA (ESPAÑOL)</b>                          | <b>18</b> |
| 1) Finalidad prevista (uso previsto)                    | 18        |
| 2) Descripción del dispositivo / Especificaciones clave | 18        |
| 3) Contraindicaciones                                   | 18        |
| 4) Advertencias / Precauciones / Riesgos potenciales    | 19        |

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## AUSTRALIA (ENGLISH)

## INSTRUCTIONS FOR USE (IFU)

Australia — Patient Matched Medical Device (PMMD)

|                              |                                                                                                                                                                                                              |
|------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>GMDN / EMDN</b>           | 38602                                                                                                                                                                                                        |
| <b>Product Code(s)</b>       | WAXRIM, SWAXPH                                                                                                                                                                                               |
| <b>Manufacturer</b>          | <b>Zima International, Inc. Also DBA Dandy</b><br>1320 N. 300 W., Lehi, UT 84043<br>United States                                                                                                            |
| <b>Australian Sponsor</b>    | <b>Name: DANDY Labs AUS Pty Ltd</b><br><b>Company No.:</b> ABN: 33693955761<br><b>Address:</b> Tower One - International Towers, Level 46, 100<br>Barangaroo Avenue, Sydney, New South Wales 2000, Australia |
| <b>Intended User</b>         | Prescriber — Dental professional / dental practice                                                                                                                                                           |
| <b>Device Type</b>           | Patient Matched Medical Device (PMMD)                                                                                                                                                                        |
| <b>Sterility</b>             | Not supplied sterile                                                                                                                                                                                         |
| <b>Single-use / Reusable</b> | Single-use                                                                                                                                                                                                   |

### 1) Intended Purpose (Intended Use)

Dandy's Wax Rim is a dental appliance that creates an accurate and comfortable fit for restorations. It serves as a temporary framework to guide the bite registration process and establish a VDO. Dandy's Wax Bite Rim is a digitally designed and manufactured dental appliance. It is made from high-quality materials and crafted to provide a precise fit for optimal patient comfort.

### 2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: Acrylic base with wax blocks.
- Material:

[Full Dentures Material Safety Data Sheets](#)

[Partial Dentures Material Safety Data Sheets](#)

### 3) Contraindications

**The Wax Bite Rim should not be used in the following situations:**

- Severe gum disease or oral infections
- Uncontrolled bleeding disorders

- Allergy or hypersensitivity to the materials used in the device
- Patients who are unable to follow instructions or communicate effectively

## 4) Warnings / Precautions / Potential Risks

### ⚠ WARNING

- Using the Wax Bite Rim carries certain risks, including
- potential allergic reactions and temporary discomfort.
- To ensure safe and effective use of the Wax Bite Rim, please be aware of the following:
- This device is intended for use by dental professionals only.
- Follow the case-specific instructions provided by Dandy carefully.
- Do not modify or adjust the device without professional guidance.
- Discontinue use if any discomfort or adverse reactions occur.

### ⚠ CAUTION

NA

## 5) Cleaning and Care

### Directions for Use

- The Wax Bite Rim aims to provide a precise fit, resulting in improved functionality and aesthetics.
- The associated procedure involves taking a bite registration to record the patient's occlusal relationship.

### Setup instructions:

- Before using the device, thoroughly wash your hands with soap and water.
- Position the Denture Wax Bite Rim in the patient's mouth, evaluate the fit and dimensions of the wax rim intra-orally, and adjust accordingly.
- Border mold the periphery of the wax rim using compound, wax, or heavy body PVS material, having patients perform functional movements. (It is recommended to use PVS tray adhesive.)
- Apply a light body wash and have patients repeat functional movements.
- Trim any excess impression material that may be covering the occlusal or facial portion of the wax rim.
- Adjust the wax rim for proper VDO, bite, and lip support.
- Evaluate the "F, V," and "S" sounds to verify the dimensions of the wax rim.
- Clearly mark the wax bite rim with the patient's smile lines, canine lines, and midlines.
- Also, cut notches in the occlusal and mark up the rest of the unmarked wax rim to break up the monotony of the flat wax. This will greatly assist in capturing the 360-degree scan.
- Capture bite registration and trim excess material away from the buccal and facial sides of the wax rim.
- Perform a 360-degree extraoral scan of the wax bite rim(s).
- Follow any additional instructions provided regarding specific adjustments or tasks.

## 6) Storage

To store dental Wax Rims, keep them clean, cool, and away from direct sun/heat to prevent distortion; for short-term, a sealed container works, but for longer storage, especially if they're part of a denture base, keep them moist (in water/solution) in a secure, covered case to maintain shape and prevent drying, but avoid hot water which can warp them

## 7) Expected Life / Service

Wax rims are custom-made for each individual patient (will not fit other patients). They come in contact with the patient for only a few minutes (long enough to be adjusted and establish the VDO). They are intended for single use; however, they can be reused on the same patient when:

- VDO needs to be reattained
- They haven't been distorted or crushed
- High Heat
- Occlusal forces
- Material residue

## 8) Incident / Complaint Reporting

### Report suspected serious incidents associated with this device to:

- The Manufacturer and/or Dandy Australia entity using the contacts above, and
- Therapeutic Goods Administration (TGA)

In Australia, serious incidents involving medical devices must be reported to the TGA in accordance with the Therapeutic Goods Act 1989.

### INFO

*For full contact details of the local Authorized Representative, refer to: Authorized Representatives — end of this document.*

## 9) Disposal

- Wax rim disposal typically involves handling it as either regulated biohazardous waste or general waste, often through specialized handlers for incineration or secure landfill disposal.
- Always avoid putting them in regular trash without checking local rules.
- Biohazard/Infectious Waste: Treat used or unclaimed wax rims as biohazardous due to saliva/blood; use red bags and authorized carriers.
- Incineration: A common method for final destruction.
- Confidentiality: Redact or destroy patient info on lab boxes before disposal.

## New Zealand (ENGLISH)

## INSTRUCTIONS FOR USE (IFU)

New Zealand— Patient Matched Medical Device (PMMD)

|                       |                                                                                                                       |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------|
| GMDN / EMDN           | 38602                                                                                                                 |
| Product Code(s)       | WAXRIM, SWAXPH                                                                                                        |
| Manufacturer          | <b>Zima International, Inc. Also DBA Dandy</b><br>1320 N. 300 W., Lehi, UT 84043<br>United States                     |
| New Zealand Sponsor   | <b>Name: DANDY Labs NZ Limited</b><br><b>Address:</b> Quigg Partners, Level 7, 36 Brandon Street, Wellington, NZ 6011 |
| Intended User         | Prescriber — Dental professional / dental practice                                                                    |
| Device Type           | Patient Matched Medical Device (PMMD)                                                                                 |
| Sterility             | Not supplied sterile                                                                                                  |
| Single-use / Reusable | Single-use                                                                                                            |

## 1) Intended Purpose (Intended Use)

Dandy's Wax Rim is a dental appliance that creates an accurate and comfortable fit for restorations. It serves as a temporary framework to guide the bite registration process and establish a VDO. Dandy's Wax Bite Rim is a digitally designed and manufactured dental appliance. It is made from high-quality materials and crafted to provide a precise fit for optimal patient comfort.

## 2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: Acrylic base with wax blocks.
- Material:

[Full Dentures Material Safety Data Sheets](#)[Partial Dentures Material Safety Data Sheets](#)

## 3) Contraindications

**The Wax Bite Rim should not be used in the following situations:**

- Severe gum disease or oral infections
- Uncontrolled bleeding disorders

- Allergy or hypersensitivity to the materials used in the device
- Patients who are unable to follow instructions or communicate effectively

#### 4) Warnings / Precautions / Potential Risks

##### ⚠ WARNING

- Using the Wax Bite Rim carries certain risks, including
- potential allergic reactions and temporary discomfort.
- To ensure safe and effective use of the Wax Bite Rim, please be aware of the following:
- This device is intended for use by dental professionals only.
- Follow the case-specific instructions provided by Dandy carefully.
- Do not modify or adjust the device without professional guidance.
- Discontinue use if any discomfort or adverse reactions occur.

##### ⚠ CAUTION

NA

#### 5) Cleaning and Care

##### Directions for Use

- The Wax Bite Rim aims to provide a precise fit, resulting in improved functionality and aesthetics.
- The associated procedure involves taking a bite registration to record the patient's occlusal relationship.

##### Setup instructions:

- Before using the device, thoroughly wash your hands with soap and water.
- Position the Denture Wax Bite Rim in the patient's mouth, evaluate the fit and dimensions of the wax rim intra-orally, and adjust accordingly.
- Border mold the periphery of the wax rim using compound, wax, or heavy body PVS material, having patients perform functional movements. (It is recommended to use PVS tray adhesive.)
- Apply a light body wash and have patients repeat functional movements.
- Trim any excess impression material that may be covering the occlusal or facial portion of the wax rim.
- Adjust the wax rim for proper VDO, bite, and lip support.
- Evaluate the "F, V," and "S" sounds to verify the dimensions of the wax rim.
- Clearly mark the wax bite rim with the patient's smile lines, canine lines, and midlines.
- Also, cut notches in the occlusal and mark up the rest of the unmarked wax rim to break up the monotony of the flat wax. This will greatly assist in capturing the 360-degree scan.
- Capture bite registration and trim excess material away from the buccal and facial sides of the wax rim.
- Perform a 360-degree extraoral scan of the wax bite rim(s).
- Follow any additional instructions provided regarding specific adjustments or tasks.

#### 6) Storage

To store dental Wax Rims, keep them clean, cool, and away from direct sun/heat to prevent distortion; for short-term, a sealed container works, but for longer storage, especially if they're part of a denture base, keep them moist (in water/solution) in a secure, covered case to maintain shape and prevent drying, but avoid hot water which can warp them

## 7) Expected Life / Service

Wax rims are custom-made for each individual patient (will not fit other patients). They come in contact with the patient for only a few minutes (long enough to be adjusted and establish the VDO). They are intended for single use; however, they can be reused on the same patient when:

- VDO needs to be reattained
- They haven't been distorted or crushed
- High Heat
- Occlusal forces
- Material residue

## 8) Incident / Complaint Reporting

### Report suspected serious incidents associated with this device to:

- The Manufacturer and/or Dandy New Zealand entity using the contacts above, and
- Medsafe

In New Zealand, serious incidents and quality issues involving medical devices should be reported to Medsafe in accordance with New Zealand's medical device adverse event reporting requirements.

### INFO

*For full contact details of the local Authorized Representative, refer to: Authorized Representatives — end of this document.*

## 9) Disposal

- Wax rim disposal typically involves handling it as either regulated biohazardous waste or general waste, often through specialized handlers for incineration or secure landfill disposal.
- Always avoid putting them in regular trash without checking local rules.
- Biohazard/Infectious Waste: Treat used or unclaimed wax rims as biohazardous due to saliva/blood; use red bags and authorized carriers.
- Incineration: A common method for final destruction.
- Confidentiality: Redact or destroy patient info on lab boxes before disposal.

## UNITED KINGDOM (ENGLISH)

## INSTRUCTIONS FOR USE (IFU)

United Kingdom — Custom-Made Device (CMD)

|                                     |                                                                                                                                                         |
|-------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>GMDN / EMDN</b>                  | 38602                                                                                                                                                   |
| <b>Product Code(s)</b>              | WAXRIM, SWAXPH                                                                                                                                          |
| <b>Manufacturer</b>                 | <b>Zima International, Inc. Also DBA Dandy</b><br>1320 N. 300 W., Lehi, UT 84043<br>United States                                                       |
| <b>UK Responsible Person (UKRP)</b> | <b>Name:</b> DANDY LABS GB, LTD<br><b>Company No.:</b> Company Number: 16873608<br><b>Address:</b> 5 New Street Square, London EC4A 3TW, United Kingdom |
| <b>Intended User</b>                | Prescriber — Dental professional / dental practice                                                                                                      |
| <b>Device Type</b>                  | Custom-made dental device (CMD)                                                                                                                         |
| <b>Sterility</b>                    | Not supplied sterile (if applicable)                                                                                                                    |
| <b>Single-use / Reusable</b>        | Single-use                                                                                                                                              |

## 1) Intended Purpose (Intended Use)

Dandy's Wax Rim is a dental appliance that creates an accurate and comfortable fit for restorations. It serves as a temporary framework to guide the bite registration process and establish a VDO. Dandy's Wax Bite Rim is a digitally designed and manufactured dental appliance. It is made from high-quality materials and crafted to provide a precise fit for optimal patient comfort.

## 2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: Acrylic base with wax blocks.
- Material:

[Full Dentures Material Safety Data Sheets](#)[Partial Dentures Material Safety Data Sheets](#)

## 3) Contraindications

**The Wax Bite Rim should not be used in the following situations:**

- Severe gum disease or oral infections
- Uncontrolled bleeding disorders

- Allergy or hypersensitivity to the materials used in the device

Patients who are unable to follow instructions or communicate effectively

#### 4) Warnings / Precautions / Potential Risks

##### **⚠ WARNING**

Using the Wax Bite Rim carries certain risks, including potential allergic reactions and temporary discomfort.

To ensure safe and effective use of the Wax Bite Rim, please be aware of the following:

This device is intended for use by dental professionals only.

Follow the case-specific instructions provided by Dandy carefully.

Do not modify or adjust the device without professional guidance.

Discontinue use if any discomfort or adverse reactions occur.

##### **⚠ CAUTION**

NA

#### 5) Cleaning and Care

##### Directions for Use

- The Wax Bite Rim aims to provide a precise fit, resulting in improved functionality and aesthetics.
- The associated procedure involves taking a bite registration to record the patient's occlusal relationship.

##### Setup instructions:

- Before using the device, thoroughly wash your hands with soap and water.
- Position the Denture Wax Bite Rim in the patient's mouth, evaluate the fit and dimensions of the wax rim intra-orally, and adjust accordingly.
- Border mold the periphery of the wax rim using compound, wax, or heavy body PVS material, having patients perform functional movements. (It is recommended to use PVS tray adhesive.)
- Apply a light body wash and have patients repeat functional movements.
- Trim any excess impression material that may be covering the occlusal or facial portion of the wax rim.
- Adjust the wax rim for proper VDO, bite, and lip support.
- Evaluate the "F, V," and "S" sounds to verify the dimensions of the wax rim.
- Clearly mark the wax bite rim with the patient's smile lines, canine lines, and midlines.
- Also, cut notches in the occlusal and mark up the rest of the unmarked wax rim to break up the monotony of the flat wax. This will greatly assist in capturing the 360-degree scan.
- Capture bite registration and trim excess material away from the buccal and facial sides of the wax rim.
- Perform a 360-degree extraoral scan of the wax bite rim(s).
- Follow any additional instructions provided regarding specific adjustments or tasks.

#### 6) Storage

To store dental Wax Rims, keep them clean, cool, and away from direct sun/heat to prevent distortion; for short-term, a sealed container works, but for longer storage, especially if they're part of a denture base, keep them moist (in water/solution) in a secure, covered case to maintain shape and prevent drying, but avoid hot water which can warp them

## 7) Expected Life / Service

The Wax Bite Rim is designed to withstand normal usage. However, like any dental appliance, it may eventually require replacement due to wear and tear.

The expected lifespan of the device can vary depending on individual factors

## 8) Incident / Complaint Reporting

Report suspected serious incidents associated with this device to:

- The Manufacturer and/or UK Responsible Person (UKRP) using the contacts above, and
- The UK Medicines and Healthcare products Regulatory Agency (MHRA) — [www.gov.uk/mhra](http://www.gov.uk/mhra)

In accordance with The Medical Devices Regulations 2002 (as amended).

### INFO

For full contact details of the local Authorized Representative, refer to: Authorized Representatives — end of this document.

## 9) Disposal

- Wax rim disposal typically involves handling it as either regulated biohazardous waste or general waste, often through specialized handlers for incineration or secure landfill disposal.
- Always avoid putting them in regular trash without checking local rules.
- Biohazard/Infectious Waste: Treat used or unclaimed wax rims as biohazardous due to saliva/blood; use red bags and authorized carriers.
- Incineration: A common method for final destruction.
- Confidentiality: Redact or destroy patient info on lab boxes before disposal.

## EUROPEAN UNION (ENGLISH)

## INSTRUCTIONS FOR USE (IFU)

European Union — Custom-Made Device (CMD)

|                                            |                                                                                                                                                      |
|--------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>GMDN / EMDN</b>                         | Q01020203: Occlusal Waxes                                                                                                                            |
| <b>Product Code(s)</b>                     | WAXRIM, SWAXPH                                                                                                                                       |
| <b>Manufacturer</b>                        | <b>Zima International, Inc. Also DBA Dandy</b><br>1320 N. 300 W., Lehi, UT 84043<br>United States                                                    |
| <b>EU Authorized representative (EURP)</b> | <b>Name:</b> Dandy Labs Europe<br><b>Company No.:</b> Registration No.: 989857065<br><b>Address:</b> 3, boulevard de Sebastopol, 75001 Paris, France |
| <b>Intended User</b>                       | Prescriber — Dental professional / dental practice                                                                                                   |
| <b>Device Type</b>                         | Custom-made dental device (CMD)                                                                                                                      |
| <b>Sterility</b>                           | Not supplied sterile (if applicable)                                                                                                                 |
| <b>Single-use / Reusable</b>               | Single-use                                                                                                                                           |

## 1) Intended Purpose (Intended Use)

Dandy's Wax Rim is a dental appliance that creates an accurate and comfortable fit for restorations. It serves as a temporary framework to guide the bite registration process and establish a VDO. Dandy's Wax Bite Rim is a digitally designed and manufactured dental appliance. It is made from high-quality materials and crafted to provide a precise fit for optimal patient comfort.

## 2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: Acrylic base with wax blocks.
- Material:

[Full Dentures Material Safety Data Sheets](#)[Partial Dentures Material Safety Data Sheets](#)

## 3) Contraindications

**The Wax Bite Rim should not be used in the following situations:**

- Severe gum disease or oral infections
- Uncontrolled bleeding disorders

- Allergy or hypersensitivity to the materials used in the device

Patients who are unable to follow instructions or communicate effectively

## 4) Warnings / Precautions / Potential Risks

### ⚠ WARNING

- Using the Wax Bite Rim carries certain risks, including
- potential allergic reactions and temporary discomfort.
- To ensure safe and effective use of the Wax Bite Rim, please be aware of the following:
- This device is intended for use by dental professionals only.
- Follow the case-specific instructions provided by Dandy carefully.
- Do not modify or adjust the device without professional guidance.
- Discontinue use if any discomfort or adverse reactions occur.

### ⚠ CAUTION

NA

## 5) Cleaning and Care

### Directions for Use

- The Wax Bite Rim aims to provide a precise fit, resulting in improved functionality and aesthetics.
- The associated procedure involves taking a bite registration to record the patient's occlusal relationship.

### Setup instructions:

- Before using the device, thoroughly wash your hands with soap and water.
- Position the Denture Wax Bite Rim in the patient's mouth, evaluate the fit and dimensions of the wax rim intra-orally, and adjust accordingly.
- Border mold the periphery of the wax rim using compound, wax, or heavy body PVS material, having patients perform functional movements. (It is recommended to use PVS tray adhesive.)
- Apply a light body wash and have patients repeat functional movements.
- Trim any excess impression material that may be covering the occlusal or facial portion of the wax rim.
- Adjust the wax rim for proper VDO, bite, and lip support.
- Evaluate the "F, V," and "S" sounds to verify the dimensions of the wax rim.
- Clearly mark the wax bite rim with the patient's smile lines, canine lines, and midlines.
- Also, cut notches in the occlusal and mark up the rest of the unmarked wax rim to break up the monotony of the flat wax. This will greatly assist in capturing the 360-degree scan.
- Capture bite registration and trim excess material away from the buccal and facial sides of the wax rim.
- Perform a 360-degree extraoral scan of the wax bite rim(s).
- Follow any additional instructions provided regarding specific adjustments or tasks.

## 6) Storage

To store dental Wax Rims, keep them clean, cool, and away from direct sun/heat to prevent distortion; for short-term, a sealed container works, but for longer storage, especially if they're part of a denture base, keep them moist (in water/solution) in a secure, covered case to maintain shape and prevent drying, but avoid hot water which can warp them

## 7) Expected Life / Service

The Wax Bite Rim is designed to withstand normal usage. However, like any dental appliance, it may eventually require replacement due to wear and tear.

The expected lifespan of the device can vary depending on individual factors

## 8) Incident / Complaint Reporting

### Report suspected serious incidents associated with this device to:

- The Manufacturer and/or Dandy EU entity using the contacts above, and corresponding local requirements
- The Competent Authority of your Member State, in accordance with Regulation (EU) 2017/745 (MDR).

### Competent Authorities:

- France: Agence Nationale de Sécurité du Médicament et des Produits de Santé (ANSM)
- Spain : Spanish Agency for Medicines and Health Products (AEMPS)
- Germany : Bundesinstitut für Arzneimittel und Medizinprodukte (BfArM) (BfArM)
- Italy : Ministero della Salute / Istituto Superiore di Sanità (ISS)
- Spain: Agencia Española de Medicamentos y Productos Sanitarios (AEMPS)
- Netherlands: Inspectie Gezondheidszorg en Jeugd (IGJ)
- Austria: Bundesamt für Sicherheit im Gesundheitswesen (BASG)
- Sweden: Läkemedelsverket
- Denmark: Lægemiddelstyrelsen
- Portugal: Autoridade Nacional do Medicamento e Produtos de Saúde (INFARMED)
- Poland: Urząd Rejestracji Produktów Leczniczych, Wyrobów Medycznych i Produktów Biobójczych (URPL)
- Norway: Direktoratet for medisinske produkter (DMP)

### INFO

For the full contact information of the local Authorized Representative, see: Authorized Representatives — at the end of the document.

## 9) Disposal

- Wax rim disposal typically involves handling it as either regulated biohazardous waste or general waste, often through specialized handlers for incineration or secure landfill disposal.
- Always avoid putting them in regular trash without checking local rules.
- Biohazard/Infectious Waste: Treat used or unclaimed wax rims as biohazardous due to saliva/blood; use red bags and authorized carriers.
- Incineration: A common method for final destruction.
- Confidentiality: Redact or destroy patient info on lab boxes before disposal.

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The table below lists the local Authorized Representative or Sponsor for each country / jurisdiction where this device is marketed. For incident reporting, contact the representative for your country.

*Le tableau ci-dessous liste le Représentant Autorisé local pour chaque pays. | Die folgende Tabelle enthält die bevollmächtigten Vertreter je Land. | La tabella seguente elenca il rappresentante autorizzato o lo sponsor locale per ciascun paese. | La tabla siguiente enumera el Representante Autorizado o Patrocinador local para cada país. | Følgende tabel viser den autoriserede repræsentant eller lokale sponsor for hvert land. | De volgende tabel geeft een overzicht van de gemachtigde vertegenwoordiger of lokale sponsor voor elk land. | W poniższej tabeli wymieniono upoważnionych przedstawicieli lub lokalnych sponsorów dla każdego kraju. | Följande tabell listar den auktoriserade representanten eller den lokala sponsorn för varje land. | A tabela seguinte lista o Representante Autorizado ou Patrocinador local para cada país. | Tabellen nedenfor viser den autoriserte representanten eller den lokale sponsoren for hvert land.*

| Jurisdiction     | Authorised Rep / Sponsor | Registered Address / Contact Email                                                   | Reg. Number      | Regulatory Body                                 |
|------------------|--------------------------|--------------------------------------------------------------------------------------|------------------|-------------------------------------------------|
| Australia (AUS)  | DANDY Labs AUS Pty Ltd   | Tower One, Level 46, 100 Barangaroo Ave, Sydney NSW 2000<br>Compliance@meetdandy.com | ABN: 33693955761 | TGA (Therapeutic Goods Administration)          |
| EU - English     | Dandy Labs Europe        | 3, boulevard de Sebastopol, 75001 Paris, France<br>Compliance@meetdandy.com          | Reg: 989857065   | Competent Authority / EUDAMED                   |
| EU - España      | Dandy Labs Europa        | 3, boulevard de Sebastopol, 75001 París, Francia<br>Compliance@meetdandy.com         | Reg: 989857065   | AEMPS / Autoridad Competente                    |
| EU - France      | Dandy Labs Europe        | 3, boulevard de Sebastopol, 75001 Paris, France<br>Compliance@meetdandy.com          | Reg: 989857065   | ANSM / Autorité Compétente                      |
| EU - Deutschland | Dandy Labs Europe        | 3, boulevard de Sebastopol, 75001 Paris, Frankreich<br>Compliance@meetdandy.com      | Reg: 989857065   | BfArM / Zuständige Behörde                      |
| EU - Italia      | Dandy Labs Europe        | 3, boulevard de Sebastopol, 75001 Parigi, Francia<br>Compliance@meetdandy.com        | Reg: 989857065   | Min. della Salute / ISS                         |
| EU - Denmark     | Dandy Labs Europe        | 3, boulevard de Sebastopol, 75001 Parigi, Francia<br>Compliance@meetdandy.com        | Reg: 989857065   | Lægemiddelstyre Isen                            |
| EU - Nederland   | Dandy Labs Europe        | 3, boulevard de Sebastopol, 75001 Parigi, Francia<br>Compliance@meetdandy.com        | Reg: 989857065   | Inspectie Gezondheidszorg en Jeugd (IGJ)        |
| EU - Polska      | Dandy Labs Europe        | 3, boulevard de Sebastopol, 75001 Parigi, Francia<br>Compliance@meetdandy.com        | Reg: 989857065   | Urząd Rejestracji Produktów Leczniczych (UR PL) |
| EU - Sverige     | Dandy Labs Europe        | 3, boulevard de Sebastopol, 75001 Parigi, Francia                                    | Reg: 989857065   | Läkemedelsverket                                |

|                     |                       |                                                                                  |                        |                                                                   |
|---------------------|-----------------------|----------------------------------------------------------------------------------|------------------------|-------------------------------------------------------------------|
|                     |                       | Compliance@meetdandy.com                                                         |                        |                                                                   |
| EU - Portugal       | Dandy Labs Europe     | 3, boulevard de Sebastopol,<br>75001 Parigi, Francia<br>Compliance@meetdandy.com | Reg:<br>989857065      | Autoridade Nacional do Medicamento e Produtos de Saúde (INFARMED) |
| EU - Norge          | Dandy Labs Europe     | 3, boulevard de Sebastopol,<br>75001 Parigi, Francia<br>Compliance@meetdandy.com | Reg:<br>989857065      | Direktoratet for medisinske produkter (DMP)                       |
| United Kingdom (UK) | DANDY LABS GB, LTD    | 5 New Street Square, London<br>EC4A 3TW<br>Compliance@meetdandy.com              | Co. No.:<br>16873608   | MHRA                                                              |
| New Zealand (NZ)    | DANDY Labs NZ Limited | Quigg Partners, Level 7, 36<br>Brandon Street, Wellington, NZ<br>6011            | NZBN:<br>9429053758446 | Medsafe                                                           |

#### NOTE

Contact details are subject to change. Always verify current information on the official product labelling or the manufacturer's website before reporting an incident. / Les coordonnées sont sujettes à modification. / Die Kontaktdaten können sich ändern. / I dati di contatto sono soggetti a modifiche. / Los datos de contacto están sujetos a cambios. / Kontaktopplysningerne kan ændres. / De contactgegevens zijn onder voorbehoud van wijzigingen. / Dane kontaktowe mogą ulec zmianie. / Kontaktoppgifterna kan komma att ändras. / Os dados de contacto estão sujeitos a alterações. / Kontaktopplysningene kan endres.

**Manufacturer / Fabricant /  
Hersteller / Produttore /  
Fabricante / Producent /  
Fabrikant / Producent /  
Tillverkare / Fabricante /  
Produsent**

**Zima International, Inc. Also DBA Dandy**  
1320 N. 300 W., Lehi, UT 84043, United States  
compliance@meetdandy.com

This is a representation of an electronic record that was signed electronically and this page is the manifestation of the electronic signatures.

Signatory Table

| Action Name | User Name      | Title                   | Signature Date    |
|-------------|----------------|-------------------------|-------------------|
| Approve     | Nicole Sheldon | Senior Quality Engineer | 07-Aug-2026 17:12 |

\* Datetime fields are displayed according to the system time zone: (GMT-04:00) Eastern Daylight Time (America/New\_York)