# INSTRUCTIONS FOR USE (IFU)

Dandy / Zima International, Inc. | AUS · NZ · UK · EU/EN · EU/ES · EU/FR · EU/DE · EU/IT · EU/DA ·  
EU/NL · EU/PL · EU/SV · EU/NB · EU/PT

## MANUFACTURER

|                          |                                               |
|--------------------------|-----------------------------------------------|
| <b>Manufacturer</b>      | Zima International, Inc. Also DBA Dandy       |
| <b>Address</b>           | 1320 N. 300 W., Lehi, UT 84043, United States |
| <b>Contact / Website</b> | www.meetdandy.com                             |
| <b>Complaints</b>        | complaints@meetdandy.com                      |

## DEVICE IDENTITY

### Device Identity

|                           |                                                                         |
|---------------------------|-------------------------------------------------------------------------|
| <b>Device Family</b>      | Dental Occlusal Splint                                                  |
| <b>GMDN / EMDN Code</b>   | 43025                                                                   |
| <b>Product Code(s)</b>    | TMJNTI, TMJAD, TMJDBS, TMJGA                                            |
| <b>Device Type</b>        | Custom-made dental device (CMD) / Patient Matched Medical Device (PMMD) |
| <b>Sterility</b>          | Not supplied sterile                                                    |
| <b>Single-use / Reuse</b> | Single-use                                                              |
| <b>Intended User</b>      | Prescriber: Dental professional / dental practice                       |

## TABLE OF CONTENT

|                                                         |           |
|---------------------------------------------------------|-----------|
| <b>AUSTRALIA (ENGLISH)</b>                              | <b>4</b>  |
| 1) Intended Purpose (Intended Use)                      | 4         |
| 2) Device Description / Key Specifications              | 4         |
| 3) Contraindications                                    | 4         |
| 4) Warnings / Precautions / Potential Risks             | 5         |
| 5) Cleaning and Care                                    | 5         |
| 6) Storage                                              | 5         |
| 7) Expected Life / Service                              | 5         |
| 8) Incident / Complaint Reporting                       | 6         |
| 9) Disposal                                             | 6         |
| <b>NEW ZEALAND (ENGLISH)</b>                            | <b>7</b>  |
| 1) Intended Purpose (Intended Use)                      | 7         |
| 2) Device Description / Key Specifications              | 7         |
| 3) Contraindications                                    | 7         |
| 4) Warnings / Precautions / Potential Risks             | 8         |
| 5) Cleaning and Care                                    | 8         |
| 6) Storage                                              | 8         |
| 7) Expected Life / Service                              | 8         |
| 8) Incident / Complaint Reporting                       | 9         |
| 9) Disposal                                             | 9         |
| <b>UNITED KINGDOM (ENGLISH)</b>                         | <b>10</b> |
| 1) Intended Purpose (Intended Use)                      | 10        |
| 2) Device Description / Key Specifications              | 10        |
| 3) Contraindications                                    | 10        |
| 4) Warnings / Precautions / Potential Risks             | 11        |
| 5) Cleaning and Care                                    | 11        |
| 6) Storage                                              | 11        |
| 7) Expected Life / Service                              | 11        |
| 8) Incident / Complaint Reporting                       | 11        |
| 9) Disposal                                             | 12        |
| <b>EUROPEAN UNION (ENGLISH)</b>                         | <b>13</b> |
| 1) Intended Purpose (Intended Use)                      | 13        |
| 2) Device Description / Key Specifications              | 13        |
| 3) Contraindications                                    | 13        |
| 4) Warnings / Precautions / Potential Risks             | 14        |
| 5) Cleaning and Care                                    | 14        |
| 6) Storage                                              | 14        |
| 7) Expected Life / Service                              | 14        |
| 8) Incident / Complaint Reporting                       | 14        |
| 9) Disposal                                             | 15        |
| <b>UNIÓN EUROPEA (ESPAÑOL)</b>                          | <b>16</b> |
| 1) Finalidad prevista (uso previsto)                    | 16        |
| 2) Descripción del dispositivo / Especificaciones clave | 16        |
| 3) Contraindicaciones                                   | 16        |

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## AUSTRALIA (ENGLISH)

## INSTRUCTIONS FOR USE (IFU)

Australia — Patient Matched Medical Device (PMMD)

|                              |                                                                                                                                                                                                              |
|------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>GMDN / EMDN</b>           | 43025                                                                                                                                                                                                        |
| <b>Product Code(s)</b>       | TMJNTI, TMJAD, TMJDBS, TMJGA                                                                                                                                                                                 |
| <b>Manufacturer</b>          | <b>Zima International, Inc. Also DBA Dandy</b><br>1320 N. 300 W., Lehi, UT 84043<br>United States                                                                                                            |
| <b>Australian Sponsor</b>    | <b>Name: DANDY Labs AUS Pty Ltd</b><br><b>Company No.:</b> ABN: 33693955761<br><b>Address:</b> Tower One - International Towers, Level 46, 100<br>Barangaroo Avenue, Sydney, New South Wales 2000, Australia |
| <b>Intended User</b>         | Prescriber — Dental professional / dental practice                                                                                                                                                           |
| <b>Device Type</b>           | Patient Matched Medical Device (PMMD)                                                                                                                                                                        |
| <b>Sterility</b>             | Not supplied sterile (if applicable)                                                                                                                                                                         |
| <b>Single-use / Reusable</b> | Single-use                                                                                                                                                                                                   |

### 1) Intended Purpose (Intended Use)

Dandy TMJ Appliances are intended to reduce pain and discomfort in the TMJ and sometimes help re-train muscles to avoid these movements.

### 2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: The devices are worn while sleeping to support patients' TMJ in the position prescribed by the dentists. The appliances are removable by the patient.
- Material: Erkodent Disc

[Night Guard Material Safety Data Sheets](#)

### 3) Contraindications

For the select subset of patients who have surgery to treat TMD or any TMJ issues, any further treatment with TMJ appliances or splints should be prescribed by the surgeon.

## 4) Warnings / Precautions / Potential Risks

### ⚠ WARNING

- Tooth movement or changes in dental occlusion
- Gingival or dental soreness
- Pain or soreness to the temporomandibular joint
- Obstruction of oral breathing
- Excessive salivation

### ⚠ CAUTION

Dentists should consider the medical history of their patients. Including history of arthritis, surgery, impact trauma, any previous TMJ treatment, respiratory disorders, or other relevant health problems, and refer the patient to the appropriate healthcare provider before prescribing the device.

## 5) Cleaning and Care

### Appliance Fit & Adjustments

- The dental staff or dentist should remove the appliance from its packaging and rinse it under water. Then, place the device in the patient's mouth, checking for the proper fit and function of the appliance. The prescribing dentist should always perform the appliance adjustments. Any adjustments in the fit of the device should easily be made chairside to ensure that the patient is pleased with the fit and function of their new appliance.

### Operating Instructions

#### CLEANING

- Brush with a soft-bristle toothbrush before each use.

#### FIT ADJUSTMENT

- The prescribing dentist or physician should make all adjustments.

#### REMOVAL INSTRUCTIONS

- To remove the appliance properly, use both index fingers to pull down the front area of the upper appliance.

#### For Dawson B-splint,

- Do this for the upper first; second, for the lower, pull up in the front area.
- To avoid breakage, do not pull down the back of the appliance.

## 6) Storage

Dry and store in the provided case when not in use..

## 7) Expected Life / Service

Expect a 2-3 year product life with the appliance.

Appliance replacement signs are cracks in the acrylic and failure of the device to maintain retention in the patient's mouth.

## 8) Incident / Complaint Reporting

---

### **Report suspected serious incidents associated with this device to:**

- The Manufacturer and/or Dandy Australia entity using the contacts above, and
- Therapeutic Goods Administration (TGA)

In Australia, serious incidents involving medical devices must be reported to the TGA in accordance with the Therapeutic Goods Act 1989.

#### **INFO**

*For full contact details of the local Authorized Representative, refer to: Authorized Representatives — end of this document.*

## 9) Disposal

---

Disposing of a TMJ appliance depends on its material and your location, but generally involves treating it as a small electronic/medical item, focusing on local rules for special waste.

If the item requires refabrication, Dandy requests that the item be returned so the hardware can be retrieved.

## NEW ZEALAND (ENGLISH)

## INSTRUCTIONS FOR USE (IFU)

New Zealand— Patient Matched Medical Device (PMMD)

|                              |                                                                                                                       |
|------------------------------|-----------------------------------------------------------------------------------------------------------------------|
| <b>GMDN / EMDN</b>           | 43025                                                                                                                 |
| <b>Product Code(s)</b>       | TMJNTI, TMJAD, TMJDBS, TMJGA                                                                                          |
| <b>Manufacturer</b>          | <b>Zima International, Inc. Also DBA Dandy</b><br>1320 N. 300 W., Lehi, UT 84043<br>United States                     |
| <b>New Zealand Sponsor</b>   | <b>Name: DANDY Labs NZ Limited</b><br><b>Address:</b> Quigg Partners, Level 7, 36 Brandon Street, Wellington, NZ 6011 |
| <b>Intended User</b>         | Prescriber — Dental professional / dental practice                                                                    |
| <b>Device Type</b>           | Patient Matched Medical Device (PMMD)                                                                                 |
| <b>Sterility</b>             | Not supplied sterile (if applicable)                                                                                  |
| <b>Single-use / Reusable</b> | Single-use                                                                                                            |

## 1) Intended Purpose (Intended Use)

Dandy TMJ Appliances are intended to reduce pain and discomfort in the TMJ and sometimes help re-train muscles to avoid these movements.

## 2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: The devices are worn while sleeping to support patients' TMJ in the position prescribed by the dentists. The appliances are removable by the patient.
- Material: Erkodent Disc

[Night Guard Material Safety Data Sheets](#)

## 3) Contraindications

For the select subset of patients who have surgery to treat TMD or any TMJ issues, any further treatment with TMJ appliances or splints should be prescribed by the surgeon.

## 4) Warnings / Precautions / Potential Risks

### ⚠ WARNING

- Tooth movement or changes in dental occlusion
- Gingival or dental soreness
- Pain or soreness to the temporomandibular joint
- Obstruction of oral breathing
- Excessive salivation

### ⚠ CAUTION

Dentists should consider the medical history of their patients. Including history of arthritis, surgery, impact trauma, any previous TMJ treatment, respiratory disorders, or other relevant health problems, and refer the patient to the appropriate healthcare provider before prescribing the device.

## 5) Cleaning and Care

### Appliance Fit & Adjustments

- The dental staff or dentist should remove the appliance from its packaging and rinse it under water. Then, place the device in the patient's mouth, checking for the proper fit and function of the appliance. The prescribing dentist should always perform the appliance adjustments. Any adjustments in the fit of the device should easily be made chairside to ensure that the patient is pleased with the fit and function of their new appliance.

### Operating Instructions

#### CLEANING

- Brush with a soft-bristle toothbrush before each use.

#### FIT ADJUSTMENT

- The prescribing dentist or physician should make all adjustments.

#### REMOVAL INSTRUCTIONS

- To remove the appliance properly, use both index fingers to pull down the front area of the upper appliance.

#### For Dawson B-splint,

- Do this for the upper first; second, for the lower, pull up in the front area.
- To avoid breakage, do not pull down the back of the appliance.

## 6) Storage

Dry and store in the provided case when not in use..

## 7) Expected Life / Service

Expect a 2-3 year product life with the appliance.

Appliance replacement signs are cracks in the acrylic and failure of the device to maintain retention in the patient's mouth.

## 8) Incident / Complaint Reporting

---

### Report suspected serious incidents associated with this device to:

- The Manufacturer and/or Dandy New Zealand entity using the contacts above, and
- Medsafe

In New Zealand, serious incidents and quality issues involving medical devices should be reported to Medsafe in accordance with New Zealand's medical device adverse event reporting requirements.

#### **INFO**

*For full contact details of the local Authorized Representative, refer to: Authorized Representatives — end of this document.*

## 9) Disposal

---

Disposing of a TMJ appliance depends on its material and your location, but generally involves treating it as a small electronic/medical item, focusing on local rules for special waste.

If the item requires refabrication, Dandy requests that the item be returned so the hardware can be retrieved.

## UNITED KINGDOM (ENGLISH)

## INSTRUCTIONS FOR USE (IFU)

United Kingdom — Custom-Made Device (CMD)

|                                     |                                                                                                                                                         |
|-------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>GMDN / EMDN</b>                  | 43025                                                                                                                                                   |
| <b>Product Code(s)</b>              | TMJNTI, TMJAD, TMJDBS, TMJGA                                                                                                                            |
| <b>Manufacturer</b>                 | <b>Zima International, Inc. Also DBA Dandy</b><br>1320 N. 300 W., Lehi, UT 84043<br>United States                                                       |
| <b>UK Responsible Person (UKRP)</b> | <b>Name:</b> DANDY LABS GB, LTD<br><b>Company No.:</b> Company Number: 16873608<br><b>Address:</b> 5 New Street Square, London EC4A 3TW, United Kingdom |
| <b>Intended User</b>                | Prescriber — Dental professional / dental practice                                                                                                      |
| <b>Device Type</b>                  | Custom-made dental device (CMD)                                                                                                                         |
| <b>Sterility</b>                    | Not supplied sterile                                                                                                                                    |
| <b>Single-use / Reusable</b>        | Single-use                                                                                                                                              |

### 1) Intended Purpose (Intended Use)

Dandy TMJ Appliances are intended to reduce pain and discomfort in the TMJ and sometimes help re-train muscles to avoid these movements.

### 2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: The devices are worn while sleeping to support patients' TMJ in the position prescribed by the dentists. The appliances are removable by the patient.
- Material: Erkodent Disc

[Night Guard Material Safety Data Sheets](#)

### 3) Contraindications

For the select subset of patients who have surgery to treat TMD or any TMJ issues, any further treatment with TMJ appliances or splints should be prescribed by the surgeon.

## 4) Warnings / Precautions / Potential Risks

### ⚠ WARNING

- Tooth movement or changes in dental occlusion
- Gingival or dental soreness
- Pain or soreness to the temporomandibular joint
- Obstruction of oral breathing
- Excessive salivation

### ⚠ CAUTION

Dentists should consider the medical history of their patients. Including history of arthritis, surgery, impact trauma, any previous TMJ treatment, respiratory disorders, or other relevant health problems, and refer the patient to the appropriate healthcare provider before prescribing the device.

## 5) Cleaning and Care

### Appliance Fit & Adjustments

- The dental staff or dentist should remove the appliance from its packaging and rinse it under water. Then, place the device in the patient's mouth, checking for the proper fit and function of the appliance. The prescribing dentist should always perform the appliance adjustments. Any adjustments in the fit of the device should easily be made chairside to ensure that the patient is pleased with the fit and function of their new appliance.

### Operating Instructions

#### CLEANING

- Brush with a soft-bristle toothbrush before each use.

#### FIT ADJUSTMENT

- The prescribing dentist or physician should make all adjustments.

#### REMOVAL INSTRUCTIONS

- To remove the appliance properly, use both index fingers to pull down the front area of the upper appliance.

#### For Dawson B-splint,

- Do this for the upper first; second, for the lower, pull up in the front area.
- To avoid breakage, do not pull down the back of the appliance.

## 6) Storage

Dry and store in the provided case when not in use. .

## 7) Expected Life / Service

Expect a 2-3 year product life with the appliance.

Appliance replacement signs are cracks in the acrylic and failure of the device to maintain retention in the patient's mouth.

## 8) Incident / Complaint Reporting

Report suspected serious incidents associated with this device to:

- The Manufacturer and/or UK Responsible Person (UKRP) using the contacts above, and
- The UK Medicines and Healthcare products Regulatory Agency (MHRA) — [www.gov.uk/mhra](http://www.gov.uk/mhra)

In accordance with The Medical Devices Regulations 2002 (as amended).

## INFO

For full contact details of the local Authorized Representative, refer to: Authorized Representatives — end of this document.

## 9) Disposal

---

- Disposing of a TMJ appliance depends on its material and your location, but generally involves treating it as a small electronic/medical item, focusing on local rules for special waste.
- If the item requires refabrication, Dandy requests that the item be returned so the hardware can be retrieved.

## EUROPEAN UNION (ENGLISH)

## INSTRUCTIONS FOR USE (IFU)

European Union — Custom-Made Device (CMD)

|                                            |                                                                                                                                                      |
|--------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>GMDN / EMDN</b>                         | Q0199: ODONTOLOGY DEVICES - OTHER                                                                                                                    |
| <b>Product Code(s)</b>                     | TMJNTI, TMJAD, TMJDBS, TMJGA                                                                                                                         |
| <b>Manufacturer</b>                        | <b>Zima International, Inc. Also DBA Dandy</b><br>1320 N. 300 W., Lehi, UT 84043<br>United States                                                    |
| <b>EU authorized representative (EURP)</b> | <b>Name:</b> Dandy Labs Europe<br><b>Company No.:</b> Registration No.: 989857065<br><b>Address:</b> 3, boulevard de Sebastopol, 75001 Paris, France |
| <b>Intended User</b>                       | Prescriber — Dental professional / dental practice                                                                                                   |
| <b>Device Type</b>                         | Custom-made dental device (CMD)                                                                                                                      |
| <b>Sterility</b>                           | Not supplied sterile                                                                                                                                 |
| <b>Single-use / Reusable</b>               | Single-use                                                                                                                                           |

## 1) Intended Purpose (Intended Use)

Dandy TMJ Appliances are intended to reduce pain and discomfort in the TMJ and sometimes help re-train muscles to avoid these movements.

## 2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: The devices are worn while sleeping to support patients' TMJ in the position prescribed by the dentists. The appliances are removable by the patient.
- Material: Erkodent Disc

[Night Guard Material Safety Data Sheets](#)

## 3) Contraindications

For the select subset of patients who have surgery to treat TMD or any TMJ issues, any further treatment with TMJ appliances or splints should be prescribed by the surgeon.

## 4) Warnings / Precautions / Potential Risks

### ⚠ WARNING

- Tooth movement or changes in dental occlusion
- Gingival or dental soreness
- Pain or soreness to the temporomandibular joint
- Obstruction of oral breathing
- Excessive salivation

### ⚠ CAUTION

Dentists should consider the medical history of their patients. Including history of arthritis, surgery, impact trauma, any previous TMJ treatment, respiratory disorders, or other relevant health problems, and refer the patient to the appropriate healthcare provider before prescribing the device.

## 5) Cleaning and Care

### Appliance Fit & Adjustments

- The dental staff or dentist should remove the appliance from its packaging and rinse it under water. Then, place the device in the patient's mouth, checking for the proper fit and function of the appliance. The prescribing dentist should always perform the appliance adjustments. Any adjustments in the fit of the device should easily be made chairside to ensure that the patient is pleased with the fit and function of their new appliance.

### Operating Instructions

#### CLEANING

- Brush with a soft-bristle toothbrush before each use.

#### FIT ADJUSTMENT

- The prescribing dentist or physician should make all adjustments.

#### REMOVAL INSTRUCTIONS

- To remove the appliance properly, use both index fingers to pull down the front area of the upper appliance.
- For Dawson B-splint,
- Do this for the upper first; second, for the lower, pull up in the front area.
- To avoid breakage, do not pull down the back of the appliance.

## 6) Storage

Dry and store in the provided case when not in use..

## 7) Expected Life / Service

Expect a 2-3 year product life with the appliance.

Appliance replacement signs are cracks in the acrylic and failure of the device to maintain retention in the patient's mouth.

## 8) Incident / Complaint Reporting

**Report suspected serious incidents associated with this device to:**

- The Manufacturer and/or Dandy EU entity using the contacts above, and corresponding local requirements
- The Competent Authority of your Member State, in accordance with Regulation (EU) 2017/745 (MDR).

**Competent Authorities:**

- France: Agence Nationale de Sécurité du Médicament et des Produits de Santé (ANSM)
- Spain : Spanish Agency for Medicines and Health Products (AEMPS)
- Germany : Bundesinstitut für Arzneimittel und Medizinprodukte (BfArM) (BfArM)
- Italy : Ministero della Salute / Istituto Superiore di Sanità (ISS)
- Spain: Agencia Española de Medicamentos y Productos Sanitarios (AEMPS)
- Netherlands: Inspectie Gezondheidszorg en Jeugd (IGJ)
- Austria: Bundesamt für Sicherheit im Gesundheitswesen (BASG)
- Sweden: Läkemedelsverket
- Denmark: Lægemiddelstyrelsen
- Portugal: Autoridade Nacional do Medicamento e Produtos de Saúde (INFARMED)
- Poland: Urząd Rejestracji Produktów Leczniczych, Wyrobów Medycznych i Produktów Biobójczych (URPL)
- Norway: Direktoratet for medisinske produkter (DMP)

**INFO**

For the full contact information of the local Authorized Representative, see: Authorized Representatives — at the end of the document.

## 9) Disposal

- Disposing of a TMJ appliance depends on its material and your location, but generally involves treating it as a small electronic/medical item, focusing on local rules for special waste.
- If the item requires refabrication, Dandy requests that the item be returned so the hardware can be retrieved.

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The table below lists the local Authorized Representative or Sponsor for each country / jurisdiction where this device is marketed. For incident reporting, contact the representative for your country.

*Le tableau ci-dessous liste le Représentant Autorisé local pour chaque pays. | Die folgende Tabelle enthält die bevollmächtigten Vertreter je Land. | La tabella seguente elenca il rappresentante autorizzato o lo sponsor locale per ciascun paese. | La tabla siguiente enumera el Representante Autorizado o Patrocinador local para cada país. | Følgende tabel viser den autoriserede repræsentant eller lokale sponsor for hvert land. | De volgende tabel geeft een overzicht van de gemachtigde vertegenwoordiger of lokale sponsor voor elk land. | W poniższej tabeli wymieniono upoważnionych przedstawicieli lub lokalnych sponsorów dla każdego kraju. | Följande tabell listar den auktoriserade representanten eller den lokala sponsorn för varje land. | A tabela seguinte lista o Representante Autorizado ou Patrocinador local para cada país. | Tabellen nedenfor viser den autoriserede representanten eller den lokale sponsoren for hvert land.*

| Jurisdiction     | Authorised Rep / Sponsor | Registered Address / Contact Email                                                   | Reg. Number      | Regulatory Body                                 |
|------------------|--------------------------|--------------------------------------------------------------------------------------|------------------|-------------------------------------------------|
| Australia (AUS)  | DANDY Labs AUS Pty Ltd   | Tower One, Level 46, 100 Barangaroo Ave, Sydney NSW 2000<br>Compliance@meetdandy.com | ABN: 33693955761 | TGA (Therapeutic Goods Administration)          |
| EU - English     | Dandy Labs Europe        | 3, boulevard de Sebastopol, 75001 Paris, France<br>Compliance@meetdandy.com          | Reg: 989857065   | Competent Authority / EUDAMED                   |
| EU - España      | Dandy Labs Europa        | 3, boulevard de Sebastopol, 75001 Paris, Francia<br>Compliance@meetdandy.com         | Reg: 989857065   | AEMPS / Autoridad Competente                    |
| EU - France      | Dandy Labs Europe        | 3, boulevard de Sebastopol, 75001 Paris, France<br>Compliance@meetdandy.com          | Reg: 989857065   | ANSM / Autorité Compétente                      |
| EU - Deutschland | Dandy Labs Europe        | 3, boulevard de Sebastopol, 75001 Paris, Frankreich<br>Compliance@meetdandy.com      | Reg: 989857065   | BfArM / Zuständige Behörde                      |
| EU - Italia      | Dandy Labs Europe        | 3, boulevard de Sebastopol, 75001 Parigi, Francia<br>Compliance@meetdandy.com        | Reg: 989857065   | Min. della Salute / ISS                         |
| EU - Denmark     | Dandy Labs Europe        | 3, boulevard de Sebastopol, 75001 Parigi, Francia<br>Compliance@meetdandy.com        | Reg: 989857065   | Lægemiddelstyrelsen                             |
| EU - Nederland   | Dandy Labs Europe        | 3, boulevard de Sebastopol, 75001 Parigi, Francia<br>Compliance@meetdandy.com        | Reg: 989857065   | Inspectie Gezondheidszorg en Jeugd (IGJ)        |
| EU - Polska      | Dandy Labs Europe        | 3, boulevard de Sebastopol, 75001 Parigi, Francia<br>Compliance@meetdandy.com        | Reg: 989857065   | Urząd Rejestracji Produktów Leczniczych (UR PL) |

|                     |                       |                                                                               |                     |                                                                   |
|---------------------|-----------------------|-------------------------------------------------------------------------------|---------------------|-------------------------------------------------------------------|
| EU - Sverige        | Dandy Labs Europe     | 3, boulevard de Sebastopol, 75001 Parigi, Francia<br>Compliance@meetdandy.com | Reg: 989857065      | Läkemedelsverket                                                  |
| EU - Portugal       | Dandy Labs Europe     | 3, boulevard de Sebastopol, 75001 Parigi, Francia<br>Compliance@meetdandy.com | Reg: 989857065      | Autoridade Nacional do Medicamento e Produtos de Saúde (INFARMED) |
| EU - Norge          | Dandy Labs Europe     | 3, boulevard de Sebastopol, 75001 Parigi, Francia<br>Compliance@meetdandy.com | Reg: 989857065      | Direktoratet for medisinske produkter (DMP)                       |
| United Kingdom (UK) | DANDY LABS GB, LTD    | 5 New Street Square, London EC4A 3TW<br>Compliance@meetdandy.com              | Co. No.: 16873608   | MHRA                                                              |
| New Zealand (NZ)    | DANDY Labs NZ Limited | Quigg Partners, Level 7, 36 Brandon Street, Wellington, NZ 6011               | NZBN: 9429053758446 | Medsafe                                                           |

#### NOTE

Contact details are subject to change. Always verify current information on the official product labelling or the manufacturer's website before reporting an incident. / Les coordonnées sont susceptibles de changer. / Angaben können sich ändern. / I recapiti sono soggetti a modifiche.

**Manufacturer / Hersteller /  
Fabbricante / Fabricant /  
Fabricante**

**Zima International, Inc. Also DBA Dandy**  
1320 N. 300 W., Lehi, UT 84043, United States  
compliance@meetdandy.com

This is a representation of an electronic record that was signed electronically and this page is the manifestation of the electronic signatures.

Signatory Table

| Action Name | User Name      | Title                   | Signature Date    |
|-------------|----------------|-------------------------|-------------------|
| Approve     | Nicole Sheldon | Senior Quality Engineer | 07-Aug-2026 17:13 |

\* Datetime fields are displayed according to the system time zone: (GMT-04:00) Eastern Daylight Time (America/New\_York)