INSTRUCTIONS FOR USE (IFU)

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MANUFACTURER

Manufacturer	Zima International, Inc. Also DBA Dandy
Address	1320 N. 300 W., Lehi, UT 84043, United States
Contact / Website	www.meetdandy.com
Complaints	complaints@meetdandy.com

DEVICE IDENTITY

Device Identity

Device Family	Fixed Restoration
GMDN / EMDN Code	Q019005
Product Code(s)	INLCRCC, INLCRFB, INLCRZR
Device Type	Custom-made dental device (CMD) / Patient Matched Medical Device (PMMD)
Sterility	Not supplied sterile
Single-use / Reuse	Single-Use
Intended User	Prescriber: Dental professional / dental practice

TABLE OF CONTENT	
EUROPEAN UNION (ENGLISH)	3
1) Intended Purpose (Intended Use)	3
2) Device Description / Key Specifications	3
3) Contraindications	4
4) Warnings / Precautions / Potential Risks	4
5) Cleaning and Care	4
6) Storage	5
7) Expected Life / Service	5
8) Incident / Complaint Reporting	5
9) Disposal	6
UNION EUROPÉENNE (FRANÇAIS)	7
1) Usage prévu (utilisation prévue)	7
2) Description du dispositif / Caractéristiques principales	7
3) Contre-indications	8
4) Avertissements / Précautions / Risques potentiels	8
5) Nettoyage et entretien	9
6) Stockage	9
7) Durée de vie / Service attendu	9
8) Signalement des incidents / Réclamations	9
9) Élimination	10
AUTHORIZED REPRESENTATIVES / REPRÉSENTANTS AUTORISÉS / BEVOLLMÄCHTIGTE VERTRETER / RAPPRESENTANTI AUTORIZZATI / REPRESENTANTES AUTORIZADOS	11

EUROPEAN UNION (ENGLISH)

INSTRUCTIONS FOR USE (IFU)

European Union — Custom-Made Device (CMD)

GMDN / EMDN	Q019005
Product Code(s)	INLCRCC, INLCRFB, INLCRZR
Manufacturer	Zima International, Inc. Also DBA Dandy 1320 N. 300 W., Lehi, UT 84043 United States
EU Authorized Rep	Name: Dandy Labs Europe Company No.: Registration No.: 989857065 Address: 3, boulevard de Sebastopol, 75001 Paris, France
Intended User	Prescriber — Dental professional / dental practice
Device Type	Custom made Medical Device (CMD)
Sterility	Not supplied sterile (if applicable)
Single-use / Reusable	Single-Use

1) Intended Purpose (Intended Use)

An inlay core (often called a post and core) is a custom-made dental restoration used to rebuild a severely damaged or devitalized tooth. It consists of a post that anchors inside the tooth's root and a "false stump" that protrudes above the gumline to support a dental crown.

2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.

- Configuration:

The Post (Pivot): Fits directly into the prepared root canal to provide mechanical strength and anchoring.

The Core: Reconstructs the visible coronal portion of the tooth, acting as the stump over which a final prosthetic crown is fitted.

Both procedures are performed by a dental professional.

- Material: Metal (non-precious), Fiber, or Zirconia

[Metal](#) / [Fiber](#) / [Zirconia](#)

3) Contraindications

An Inlay Core is contraindicated if the tooth has insufficient or severely curved root structure, active endodontic pathology, or adequate remaining dentin to retain a build-up. Placing a post in these scenarios risks root fracture, perforation, or premature failure.

4) Warnings / Precautions / Potential Risks

WARNING

Surgery in the oral cavity and oral rehabilitation include general risks.

These include (but are not limited to):

- Small components accidentally dropped in the patient's mouth can be swallowed or aspirated.
- Allergy or hypersensitivity to the chemical ingredients of the material used.

CAUTION

A post and core restores a heavily compromised, root-canal-treated tooth.

Leaving 4–5 mm of gutta-percha apical seal

- Ensures a 2 mm "ferrule" (collar of natural tooth structure)

The post preparation should be smaller than one-third of the root's diameter to keep the root structurally sound.

Avoid Rigid/Threaded Posts: Threaded (screwed-in) metal posts exert intense pressure on the brittle root structure.

- Use a more passive and flexible fiber posts which flex with the tooth to reduce stress.

Root Perforation

- If the drill strays from the canal's natural path or the dentist applies too much force, the root can be punctured. This significantly reduces the tooth's lifespan and may require extraction.

Root Fracture

- Teeth are much more brittle after root canals. Poor post placement or an inadequate ferrule height can cause the root to crack under chewing forces.

Microleakage and Decay

- If the temporary restoration is loose or if the crown fails to seal completely, bacteria can seep in and cause decay within the remaining natural tooth structure.

5) Cleaning and Care

The Inlay Core is protected by either a temporary restoration or the final restoration, so home care is not required; however.

1. **Protect the Temporary:** Post-and-core buildup teeth are inherently weak until the permanent crown is placed. Avoid chewing hard, sticky, or crunchy foods on the affected side until the permanent crown is securely bonded.
2. **Avoid Delays:** Get your permanent crown placed as soon as your dentist advises. Temporary crowns can loosen over time, exposing the post to contamination.

6) Storage

Inlay Cores should be stored in their original packaging.

- Temporary (Provisional) Storage
 - If the patient must wait for the permanent crown, the prepared post-and-core space must be temporarily sealed:
 - The most secure method involves filling the post space with a material like Teflon, covered with petroleum jelly, and fabricating a provisional crown that seals the coronal portion.
 - Leaving a post-and-core without a crown risks rapid bacterial decay, root fracture, and bond failure.

7) Expected Life / Service

The average life expectancy of an Inlay Core ranges from 10 to 15 years. While some can last a lifetime, their longevity depends heavily on the type of restoration, the remaining tooth structure (ferrule), and oral hygiene.

8) Incident / Complaint Reporting

Report suspected serious incidents associated with this device to:

- The Manufacturer and/or Dandy EU entity using the contacts above, and corresponding local requirements
- The Competent Authority of your Member State, in accordance with Regulation (EU) 2017/745 (MDR).

Competent Authorities:

- France: Agence Nationale de Sécurité du Médicament et des Produits de Santé (ANSM)
- Spain : Spanish Agency for Medicines and Health Products (AEMPS)
- Germany : Bundesinstitut für Arzneimittel und Medizinprodukte (BfArM) (BfArM)
- Italy : Ministero della Salute / Istituto Superiore di Sanità (ISS)
- Spain: Agencia Española de Medicamentos y Productos Sanitarios (AEMPS) for Spain
- Netherlands: Inspectie Gezondheidszorg en Jeugd (IGJ)
- Austria: Bundesamt für Sicherheit im Gesundheitswesen (BASG)
- Sweden: Läkemedelsverket
- Denmark: Lægemiddelstyrelsen
- Portugal: Autoridade Nacional do Medicamento e Produtos de Saúde (INFARMED)

INFO

For the full contact information of the local Authorized Representative, see: Authorized Representatives — at the end of the document.

9) Disposal

Dispose of extracted teeth with any posts as regulated medical waste, but never in containers destined for incineration. Separate the metal, zirconia, or fiber using a licensed medical waste transporter or recycler.

Comply with the current applicable national waste disposal regulations.

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The table below lists the local Authorized Representative or Sponsor for each country / jurisdiction where this device is marketed. For incident reporting, contact the representative for your country.

La tabla siguiente enumera el Representante Autorizado o Patrocinador local para cada país. | Le tableau ci-dessous liste le Représentant Autorisé local pour chaque pays. | Die folgende Tabelle enthält die bevollmächtigten Vertreter je Land.

Jurisdiction	Authorised Rep / Sponsor	Registered Address / Contact Email	Reg. Number	Regulatory Body
Australia (AUS)	DANDY Labs AUS Pty Ltd	Tower One, Level 46, 100 Barangaroo Ave, Sydney NSW 2000 Compliance@meetdandy.com	ABN: 33693955761	TGA (Therapeutic Goods Administration)
EU — English	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, France Compliance@meetdandy.com	Reg: 989857065	Competent Authority / EUDAMED
EU — España	Dandy Labs Europa	3, boulevard de Sebastopol, 75001 París, Francia Compliance@meetdandy.com	Reg: 989857065	AEMPS / Autoridad Competente
EU — France	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, France Compliance@meetdandy.com	Reg: 989857065	ANSM / Autorité Compétente
EU — Deutschland	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Paris, Frankreich Compliance@meetdandy.com	Reg: 989857065	BfArM / Zuständige Behörde
EU — Italia	Dandy Labs Europe	3, boulevard de Sebastopol, 75001 Parigi, Francia Compliance@meetdandy.com	Reg: 989857065	Min. della Salute / ISS
United Kingdom (UK)	DANDY LABS GB, LTD	5 New Street Square, London EC4A 3TW Compliance@meetdandy.com	Co. No.: 16873608	MHRA

NOTE

Contact details are subject to change. Always verify current information on the official product labelling or the manufacturer's website before reporting an incident. / Les coordonnées sont susceptibles de changer. / Angaben können sich ändern. / I recapiti sono soggetti a modifiche.

**Manufacturer / Hersteller /
Fabbricante / Fabricant /
Fabricante**

Zima International, Inc. Also DBA Dandy
1320 N. 300 W., Lehi, UT 84043, United States
compliance@meetdandy.com

This is a representation of an electronic record that was signed electronically and this page is the manifestation of the electronic signatures.

Signatory Table

Action Name	User Name	Title	Signature Date
Approve	Kasha Lim	Supplier Quality Engineer	29-Jun-2026 14:07

* Datetime fields are displayed according to the system time zone: (GMT-04:00) Eastern Daylight Time (America/New_York)