Scan Bodies

Instructions for use

dandy

Scan Bodies

Instructions for use

UK	Instructions for use	3	NO	Bruksanvisning	51
FR	Mode d'emploi	7	DK	Brugsvejledning	55
DE	Gebrauchsanweisung	11	BG	Инструкции за употреба	59
ES	Instrucciones de uso	15	HR	Upute za uporabu	63
IT	Istruzioni per l'uso	19	CZ	Návod k použití	67
PT	Instruções de utilização	23	ET	Kasutusjuhend	71
PL	Instrukcje użytkowania	27	HU	Használati útmutató	75
RO	Инструкция по эксплуатации	31	LV	Lietošanas pamācība	79
GR	Οδηγίες χρήσης	35	LT	Naudojimo instrukcija	83
FI	Käyttöohje	39	SK	Návod na použitie	87
NL	Gebruiksaanwijzing	43	TR	Kullanım talimatları	91
SE	Bruksanvisning	47			

Legend of labels and symbols

	WARNING Safety information regarding risk
	CAUTION Important notice that must be strictly observed
	IMPORTANT NOTICE Important notice that must be strictly observed
	Date of manufacture
	Manufacturer
	CE marking
	Keep away from sunlight
	Use by date

	Consult Instructions for Use
	Catalogue or reference number
	LOT number
	WEEE Product disposal symbol
	Prescription use only
	Temperature limit
	The product is delivered non-sterile

	Instructions for use A printed Instructions for Use will be provided at no additional cost within 7 calendar days, if requested. Contact your reseller or customer support at compliance@meetdandy.com .
---	--

Intended use

Scan bodies are used to capture the location, angulation and connection orientation of dental implants in a patient's mouth. For intra-oral use, the scan bodies are mounted directly on the patient's implant and scanned (captured) with an intraoral scanner.

The scan bodies can also be used in gypsum models. For model use, the scan bodies are mounted on an implant analog and scanned with a laboratory scanner.

In patient use – step by step

Step 1: Selection of the scan body

Select the scan body that matches the implant system and platform.

	IMPORTANT NOTICE Each scan body contains human readable markings with a code indicating the corresponding implant platform and system.
	CAUTION To avoid damage to the implant and ensure that correct scan data is obtained, it is recommended that an in-depth knowledge of preferred systems is obtained enabling a quick confirmation of both the connection and the screwdriver.
	WARNING Replace the scan body if any damage, i.e. scratches, deformation, excessive wear, etc., is detected on either the scan body or the system screwdriver.

Step 2: Preparation for INTRA-ORAL use

To avoid cross-contamination between patients or contamination of operators, the scan bodies and the associated screwdriver must be cleaned and sterilized prior to each use.

Prior to use, the cleaned scan bodies and associated screwdrivers must be autoclaved at 121°C for 15 min, or at 134°C for 3 min.

After autoclaving, let the scan bodies and associated screwdrivers cool to room temperature. If still hot, skin contact with either the scan bodies or the screwdrivers may cause discomfort to the patient. The scan bodies must not be warmer than body temperature.

	WARNING Replace the scan body if any damage (scratches, deformation, excessive wear, etc.) is detected on either the scan body or the system screwdriver.
	WARNING The sterile barrier must be maintained until use!

Step 3: Mounting the scan body on the implant in the patient's mouth

Place the scan body on a clean implant and rotate until the interface locks and no further rotational movement is possible. Use the screwdriver to fix the screw finger-tight.

	IMPORTANT NOTICE To avoid damage and unnecessary wear, the recommended mounting torque is 5 N·cm.
	CAUTION The screwdriver does not have torque limitation!
	CAUTION To avoid damage to the scan body connection and screw, avoid exceeding the recommended mounting torque. In case of unintended damage to the scan body, screw or screwdriver making it impossible to remove the items, follow standard procedures and contact the nearest medical facility.
	WARNING To avoid aspiration or swallowing of the scan body during intra-oral scanning, ensure that the scan body matches the implant system and platform and that the screw is tightened. The scan body's ID (identifying the implant system and platform) is labeled on the head section. The standard recommendation for handling small accessories intra-orally is to tie a 20 cm long cord (Dental Floss) around the scan body. Once the scan body is in the correct position and tightened, the cord can be removed. Reverse the order of handling for safe removal. Furthermore, position the patient in such a way that, should the scan body be dropped or come loose nevertheless, it will not fall into the patient's respiratory or digestive tract. In the case of the scan body being dropped, it can be securely retrieved using the cord to pull it out.
	WARNING The standard recommendation for handling small instruments intra-orally is to tie a 20 cm long cord (Dental Floss) around the screwdriver. To avoid aspiration or swallowing of the screwdriver during mounting/unmounting the scan body, ensure to hold it firmly and position the patient in such a way that, should the screwdriver be dropped, it will not fall into the patient's respiratory or digestive tract. In the case of the screwdriver being dropped, it can be securely retrieved using the cord to pull it out.
	WARNING To avoid aspiration or swallowing despite all above preventive measures, always have an assistant in place with suction equipment ready for use. In case of swallowing of the scan body or screwdriver, follow standard procedures and contact the nearest medical facility.

Step 4: Scanning with the scan body in place

Scan the scan body following the scanner's instructions for use. Ensure that the entire surface of the scan body is scanned fully and without surface gaps, especially the head of the scan body, both top and sides.

	WARNING Position the patient in such a way that, should the scan body be dropped or come loose nevertheless, it will not fall into the patient's respiratory or digestive tract.
	WARNING To avoid aspiration or swallowing, always have an assistant in place with suction equipment ready for use. In case of swallowing of the scan body or screwdriver, follow standard procedures and contact the nearest medical facility.

Step 5: Unmounting the scan body

Use the screwdriver to unmount the scan body.

	WARNING Refer to the IMPORTANT NOTICE, CAUTIONS and WARNINGS in STEP 3 above.
---	---

Step 6: Cleaning and storage

Clean the scan body and screwdriver, if needed, and place them back in the storage box.

	WARNING The scan bodies and associated screwdrivers must be thoroughly cleaned immediately after use.
---	---

In model use – step by step**Step 1: Selection of the scan body**

Select the scan body that matches the implant system and platform.

	IMPORTANT NOTICE Each scan body and screwdriver contain human-readable markings with a code indicating the corresponding implant platform and system. The screwdrivers contain a 2-character system-specific alphanumeric code, and the scan bodies have an additional 3rd character referring to the connection platform.
	CAUTION To avoid damage to the implant analog and ensure that correct scan data is obtained, it is recommended that an in-depth knowledge of preferred systems is obtained, enabling a quick confirmation of both the connection and the screwdriver.

**CAUTION**

Replace the scan body if any damage, i.e. scratches, deformation, excessive wear, etc., is detected on either the scan body or the system screwdriver.

Step 2: Preparation for IN-MODEL use

Clean the scan body and the associated screwdriver for use.

**CAUTION**

The scan body and the analog connection must be undamaged and free from cast material.

**CAUTION**

Replace the scan body if any damage, i.e. scratches, deformation, excessive wear, etc., is detected on either scan body or system screwdriver.

Step 3: Mounting the scan body on the implant analog in the gypsum model

Place the scan body on the clean implant analog and rotate until the interface locks and no further rotational movement is possible. Use the supplied screwdriver to fix the screw finger-tight.

**IMPORTANT NOTICE**

To avoid damage and unnecessary wear, the recommended mounting torque is 5 Ncm.

**CAUTION**

The screwdriver does not have torque limitation!

Step 4: Scanning with the scan body in place

Scan the scan body following the scanner's instructions for use. Ensure that the entire surface of the scan body is scanned fully and without surface gaps, especially the head of the scan body, both top and sides.

Step 5: Unmounting the scan body

Use the supplied screwdriver to unmount the scan body.

Step 6: Cleaning and storage

Clean the scan body and screwdriver, if needed, and place them back in the storage box.

**CAUTION**

The scan bodies and associated screwdrivers must be thoroughly cleaned immediately after use.

Fading

	IMPORTANT NOTICE Towards the end of the scan bodies' lifespan, a slight change in the surface appearance can be detected. The fading is purely visual and has no impact on product functionality or scan quality.
	CAUTION Replace the scan body if any damage, i.e. scratches, deformation, excessive wear etc. is detected on either the scan body or the system screwdriver.

Storage

	The scan body kit should be stored away from direct sunlight.
	Store at room temperature.

It is recommended to use the system storage box that allows for sorting of scan bodies by implant system and platform.

Lifetime

	Discard the scan body after 100 autoclave / use cycles
	CAUTION The scan bodies may change dimensions beyond those limits and no longer allow for accurate scanning results. It is recommended that the scan bodies be inspected prior to use and discarded at the first sign of wear or deterioration.

Disposal

	Discarded items should be cleaned and disinfected prior to disposal to avoid cross-contamination. It is recommended to recycle them according to local regulations.
---	---

Contraindications

Dandy Scan Bodies should not be used for elastomeric impressions or bite registrations.

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