# INSTRUCTIONS FOR USE (IFU)

Dandy / Zima International, Inc. | AUS · NZ · UK · EU/EN · EU/ES · EU/FR · EU/DE · EU/IT · EU/DA · EU/NL · EU/PL · EU/SV · EU/NB · EU/PT

## MANUFACTURER

|                   |                                                      |
|-------------------|------------------------------------------------------|
| Manufacturer      | Zima International, Inc. Also DBA Dandy              |
| Address           | 22 Cortlandt St, 30th Floor, New York, NY, 10007, US |
| Contact / Website | www.meetdandy.com                                    |
| Complaints        | complaints@meetdandy.com                             |

## DEVICE IDENTITY

### Device Identity

|                    |                                                                                                                                                                                                                                                              |
|--------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Device Family      | Crown - Ceramic                                                                                                                                                                                                                                              |
| GMDN / EMDN Code   | 38594 / Q010206                                                                                                                                                                                                                                              |
| Product Code(s)    | EMXLPSCRA, EMXLPSCR, EMXSCRA, EMXSCR, IMPANTEM, IMPANTHTZ, IMPANTLEMM, IMPANTLVD, IMPANTZIR, IMPPOSEMM, IMPPOSHTZ, IMPPOSLAY, IMPPOSLEMM, IMPPOSZIR, PFZSCRA, PFZSCR, SHDGLXBASE, SHDGLXDARK, SHDGLXLGHT, ZIRDSSCR, ZIRHTSCRA, ZIRHTSCR, ZIRMOSCRA, ZIRMOSCR |
| Device Type        | Custom-made dental device (CMD) / Patient Matched Medical Device (PMMD)                                                                                                                                                                                      |
| Sterility          | Not supplied sterile                                                                                                                                                                                                                                         |
| Single-use / Reuse | Single-use                                                                                                                                                                                                                                                   |
| Intended User      | Prescriber: Dental professional / dental practice                                                                                                                                                                                                            |

**TABLE OF CONTENT**

|                                             |           |
|---------------------------------------------|-----------|
| <b>AUSTRALIA (ENGLISH)</b>                  | <b>6</b>  |
| 1) Intended Purpose (Intended Use)          | 6         |
| 2) Device Description / Key Specifications  | 7         |
| 3) Contraindications                        | 7         |
| 4) Warnings / Precautions / Potential Risks | 7         |
| 5) Cleaning and Care                        | 7         |
| 6) Storage                                  | 8         |
| 7) Expected Life / Service                  | 8         |
| 8) Incident / Complaint Reporting           | 8         |
| 9) Disposal                                 | 8         |
| <b>NEW ZEALAND (ENGLISH)</b>                | <b>10</b> |
| 1) Intended Purpose (Intended Use)          | 10        |
| 2) Device Description / Key Specifications  | 11        |
| 3) Contraindications                        | 11        |
| 4) Warnings / Precautions / Potential Risks | 11        |
| 5) Cleaning and Care                        | 11        |
| 6) Storage                                  | 12        |
| 7) Expected Life / Service                  | 12        |
| 8) Incident / Complaint Reporting           | 12        |
| 9) Disposal                                 | 12        |
| <b>UNITED KINGDOM (ENGLISH)</b>             | <b>14</b> |
| 1) Intended Purpose (Intended Use)          | 14        |
| 2) Device Description / Key Specifications  | 14        |
| 3) Contraindications                        | 15        |
| 4) Warnings / Precautions / Potential Risks | 15        |
| 5) Cleaning and Care                        | 15        |
| 6) Storage                                  | 16        |
| 7) Expected Life / Service                  | 16        |
| 8) Incident / Complaint Reporting           | 16        |
| 9) Disposal                                 | 16        |
| <b>EUROPEAN UNION (ENGLISH)</b>             | <b>17</b> |
| 1) Intended Purpose (Intended Use)          | 17        |
| 2) Device Description / Key Specifications  | 17        |
| 3) Contraindications                        | 18        |
| 4) Warnings / Precautions / Potential Risks | 18        |
| 5) Cleaning and Care                        | 18        |
| 6) Storage                                  | 19        |
| 7) Expected Life / Service                  | 19        |

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**AUSTRALIA (ENGLISH)****INSTRUCTIONS FOR USE (IFU)**

Australia — Patient Matched Medical Device (PMMD)

|                              |                                                                                                                                                                                                                                                         |
|------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>GMDN / EMDN</b>           | 38594 / Q010206                                                                                                                                                                                                                                         |
| <b>Product Code(s)</b>       | EMXLPSCRA, EMXLPSCR, EMXSCRA, EMXSCR, IMPANTEM, IMPANTHTZ, IMPANTLEMM, IMPANTLVD, IMPANTZIR, IMPOSEMM, IMPOSHTZ, IMPOSLAY, IMPOSLEMM, IMPOSZIR, PFZSCRA, PFZSCR, SHDGLXBASE, SHDGLXDARK, SHDGLXLGHT, ZIRDSSCR, ZIRHTSCRA, ZIRHTSCR, ZIRMOSCRA, ZIRMOSCR |
| <b>Manufacturer</b>          | <b>Zima International, Inc. Also DBA Dandy</b><br>22 Cortlandt St, 30th Floor, New York, NY, 10007<br>US                                                                                                                                                |
| <b>Australian Sponsor</b>    | <b>Name: DANDY Labs AUS Pty Ltd</b><br><b>Company No.:</b> ABN: 33693955761<br><b>Address:</b> Tower One - International Towers, Level 46, 100 Barangaroo Avenue, Sydney, New South Wales 2000, Australia                                               |
| <b>Intended User</b>         | Prescriber — Dental professional / dental practice                                                                                                                                                                                                      |
| <b>Device Type</b>           | Patient Matched Medical Device (PMMD)                                                                                                                                                                                                                   |
| <b>Sterility</b>             | Not supplied sterile (if applicable)                                                                                                                                                                                                                    |
| <b>Single-use / Reusable</b> | Single use                                                                                                                                                                                                                                              |

**1) Intended Purpose (Intended Use)**

Dandy's ceramic crowns are made to restore worn, carious, broken, and aesthetically displeasing dentition.

**They:**

- Restore single, permanent teeth (crowns) that are fractured, carious, discolored, or severely worn.
- Replace one or more missing permanent teeth by spanning the edentulous space and supporting them with dental implants.
- Restore the form, function, and esthetics of natural dentition.
- Provide occlusal stability and maintain proper interarch relationships.

**If used to restore an implant:**

- Connected to the implant fixture with an abutment
- Cases that are screw-retained are permanently fixed to the abutment in the lab

## 2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: Dandy's ceramic crowns are fabricated using digital technology, allowing for a precise fit and customization to the patient's specific needs. They are made of a highly durable and biocompatible material, ensuring long-term comfort and function.
- Material: Zirconia Monolithic (3Y & 4Y), Porcelain-Fused-To-Zirconia (PFZ), E.Max (full contour), E.Max (layered)

[Zirconia Material Safety Data Sheet](#)

[PFZ Material Safety Data Sheet](#)

[Full Contour E.MAX Material Safety Data Sheet](#)

## 3) Contraindications

Crowns should be used with caution in patients with uncontrolled periodontal disease, severe bone resorption, or oral lesions.

## 4) Warnings / Precautions / Potential Risks

### **WARNING**

The use of fixed restorations may have potential risks, including failure, misfit, breakage, and recurrent caries.

#### **Warning – Sensitivity:**

This device may contain trace levels of metals used as colorants in the ceramic material. The treating clinician should inform patients with sensitivity or allergy of the presence of metals in the device before treatment. In rare cases, allergic reactions can occur in individuals who are sensitized.

#### **Warning – Irritation:**

Localized irritation of surrounding oral tissues may occur in certain individuals due to individual variability in tissue response. Patients should be monitored post-placement, and the treating clinician should evaluate any signs of irritation, inflammation, or discomfort.

### **CAUTION**

NA

## 5) Cleaning and Care

Daily Cleaning

### **Oral Hygiene:**

#### **Brushing:**

- Recommended brushing teeth twice a day, using a soft-bristled toothbrush and a fluoride toothpaste.

#### **Flossing:**

- Recommend daily flossing to remove plaque and debris from between the teeth and around the fixed restoration. A floss threader may be needed for hard-to-reach areas (for bridges).

**Mouthwash:**

- Alcohol-based mouthwashes tend to dry out the mouth, leading to increased plaque and tartar buildup.

**Dietary and Habit Recommendations**

- Avoid abrasive foods: Limit hard, crunchy, or sticky foods that can damage the restoration or cause it to loosen.
- Avoid sugary foods: where new cavities can form or cause premature decay of the restorative margins.
- Stay hydrated: Drinking plenty of water will help keep your patient's mouth clean and prevent gum disease.
- Avoid bad habits: Refrain from using teeth to open packages or bite on hard objects.

## 6) Storage

Use dedicated dental crown/veneer boxes with soft foam or silicone inserts to prevent chipping, scratching, or movement, ensuring they are kept clean, dry, and protected during transport or while awaiting placement in a clinic or lab setting..

## 7) Expected Life / Service

Ceramic dental restorations (crowns, fillings, veneers) generally last 10 to 20 years, with modern materials like zirconia potentially exceeding this, but lifespans vary significantly based on care, location (back teeth wear faster), and type (implants are more permanent than crowns on top). Excellent oral hygiene (brushing, flossing) and avoiding hard foods or teeth grinding (bruxism) are key to maximizing longevity, though some chipping or wear is possible.

## 8) Incident / Complaint Reporting

**Report suspected serious incidents associated with this device to:**

- The Manufacturer and/or Dandy Australia entity using the contacts above, and
- Therapeutic Goods Administration (TGA)

In Australia, serious incidents involving medical devices must be reported to the TGA in accordance with the Therapeutic Goods Act 1989.

**INFO**

*For full contact details of the local Authorized Representative, refer to: Authorized Representatives — end of this document.*

## 9) Disposal

- All-ceramic (porcelain) restorations are often disposed of as regular solid waste or taken to specific recycling centers that accept mixed ceramics/construction debris, as they're challenging to recycle like typical items.

- Household Trash: Place in your regular garbage bin, as they don't break down in standard recycling
- Always check with your local municipality or waste hauler for specific rules; however, ceramic items are generally not accepted in standard curbside recycling bins.
- Recycling Centers: Take them to a local recycling center that accepts mixed ceramics, concrete, or construction debris, as they can sometimes be repurposed.
- If the restoration has been in the patient's mouth or has come into contact with saliva, discard it as biomedical waste.

**NEW ZEALAND (ENGLISH)****INSTRUCTIONS FOR USE (IFU)**

New Zealand — Patient Matched Medical Device (PMMD)

|                              |                                                                                                                                                                                                                                                         |
|------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>GMDN / EMDN</b>           | 38594 / Q010206                                                                                                                                                                                                                                         |
| <b>Product Code(s)</b>       | EMXLPSCRA, EMXLPSCR, EMXSCRA, EMXSCR, IMPANTEM, IMPANTHTZ, IMPANTLEMM, IMPANTLVD, IMPANTZIR, IMPOSEMM, IMPOSHTZ, IMPOSLAY, IMPOSLEMM, IMPOSZIR, PFZSCRA, PFZSCR, SHDGLXBASE, SHDGLXDARK, SHDGLXLGHT, ZIRDSSCR, ZIRHTSCRA, ZIRHTSCR, ZIRMOSCRA, ZIRMOSCR |
| <b>Manufacturer</b>          | <b>Zima International, Inc. Also DBA Dandy</b><br>22 Cortlandt St, 30th Floor, New York, NY, 10007<br>US                                                                                                                                                |
| <b>New Zealand Sponsor</b>   | <b>Name: DANDY Labs NZ Limited</b><br><b>Address:</b> Quigg Partners, Level 7, 36 Brandon Street, Wellington, NZ 6011                                                                                                                                   |
| <b>Intended User</b>         | Prescriber — Dental professional / dental practice                                                                                                                                                                                                      |
| <b>Device Type</b>           | Patient Matched Medical Device (PMMD)                                                                                                                                                                                                                   |
| <b>Sterility</b>             | Not supplied sterile (if applicable)                                                                                                                                                                                                                    |
| <b>Single-use / Reusable</b> | Single use                                                                                                                                                                                                                                              |

**1) Intended Purpose (Intended Use)**

Dandy's ceramic crowns are made to restore worn, carious, broken, and aesthetically displeasing dentition.

**They:**

- Restore single, permanent teeth (crowns) that are fractured, carious, discolored, or severely worn.
- Replace one or more missing permanent teeth by spanning the edentulous space and supporting them with dental implants.
- Restore the form, function, and esthetics of natural dentition.
- Provide occlusal stability and maintain proper interarch relationships.

**If used to restore an implant:**

- Connected to the implant fixture with an abutment
- Cases that are screw-retained are permanently fixed to the abutment in the lab

## 2) Device Description / Key Specifications

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: Dandy's ceramic crowns are fabricated using digital technology, allowing for a precise fit and customization to the patient's specific needs. They are made of a highly durable and biocompatible material, ensuring long-term comfort and function.
- Material: Zirconia Monolithic (3Y & 4Y), Porcelain-Fused-To-Zirconia (PFZ), E.Max (full contour), E.Max (layered)

[Zirconia Material Safety Data Sheet](#)

[PFZ Material Safety Data Sheet](#)

[Full Contour E.MAX Material Safety Data Sheet](#)

## 3) Contraindications

Crowns should be used with caution in patients with uncontrolled periodontal disease, severe bone resorption, or oral lesions.

## 4) Warnings / Precautions / Potential Risks

### **WARNING**

The use of fixed restorations may have potential risks, including failure, misfit, breakage, and recurrent caries.

#### **Warning – Sensitivity:**

This device may contain trace levels of metals used as colorants in the ceramic material. The treating clinician should inform patients with sensitivity or allergy of the presence of metals in the device before treatment. In rare cases, allergic reactions can occur in individuals who are sensitized.

#### **Warning – Irritation:**

Localized irritation of surrounding oral tissues may occur in certain individuals due to individual variability in tissue response. Patients should be monitored post-placement, and the treating clinician should evaluate any signs of irritation, inflammation, or discomfort.

### **CAUTION**

NA

## 5) Cleaning and Care

Daily Cleaning

### **Oral Hygiene:**

#### **Brushing:**

- Recommended brushing teeth twice a day, using a soft-bristled toothbrush and a fluoride toothpaste.

#### **Flossing:**

- Recommend daily flossing to remove plaque and debris from between the teeth and around the fixed restoration. A floss threader may be needed for hard-to-reach areas (for bridges).

#### Mouthwash:

- Alcohol-based mouthwashes tend to dry out the mouth, leading to increased plaque and tartar buildup.

#### Dietary and Habit Recommendations

- Avoid abrasive foods: Limit hard, crunchy, or sticky foods that can damage the restoration or cause it to loosen.
- Avoid sugary foods: where new cavities can form or cause premature decay of the restorative margins.
- Stay hydrated: Drinking plenty of water will help keep your patient's mouth clean and prevent gum disease.
- Avoid bad habits: Refrain from using teeth to open packages or bite on hard objects.

## 6) Storage

Use dedicated dental crown/veneer boxes with soft foam or silicone inserts to prevent chipping, scratching, or movement, ensuring they are kept clean, dry, and protected during transport or while awaiting placement in a clinic or lab setting..

## 7) Expected Life / Service

Ceramic dental restorations (crowns, fillings, veneers) generally last 10 to 20 years, with modern materials like zirconia potentially exceeding this, but lifespans vary significantly based on care, location (back teeth wear faster), and type (implants are more permanent than crowns on top). Excellent oral hygiene (brushing, flossing) and avoiding hard foods or teeth grinding (bruxism) are key to maximizing longevity, though some chipping or wear is possible.

## 8) Incident / Complaint Reporting

#### Report suspected serious incidents associated with this device to:

- The Manufacturer and/or Dandy New Zealand entity using the contacts above, and
- Medsafe

In New Zealand, serious incidents and quality issues involving medical devices should be reported to Medsafe in accordance with New Zealand's medical device adverse event reporting requirements.

#### INFO

*For full contact details of the local Authorized Representative, refer to: Authorized Representatives — end of this document.*

## 9) Disposal

- All-ceramic (porcelain) restorations are often disposed of as regular solid waste or taken to specific recycling centers that accept mixed ceramics/construction debris, as they're challenging to recycle like typical items.

- Household Trash: Place in your regular garbage bin, as they don't break down in standard recycling
- Always check with your local municipality or waste hauler for specific rules; however, ceramic items are generally not accepted in standard curbside recycling bins.
- Recycling Centers: Take them to a local recycling center that accepts mixed ceramics, concrete, or construction debris, as they can sometimes be repurposed.
- If the restoration has been in the patient's mouth or has come into contact with saliva, discard it as biomedical waste.

**UNITED KINGDOM (ENGLISH)****INSTRUCTIONS FOR USE (IFU)**

United Kingdom — Custom-Made Device (CMD)

|                                     |                                                                                                                                                                                                                                                         |
|-------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>GMDN / EMDN</b>                  | 38594 / Q010206                                                                                                                                                                                                                                         |
| <b>Product Code(s)</b>              | EMXLPSCRA, EMXLPSCR, EMXSCRA, EMXSCR, IMPANTEM, IMPANTHTZ, IMPANTLEMM, IMPANTLVD, IMPANTZIR, IMPOSEMM, IMPOSHTZ, IMPOSLAY, IMPOSLEMM, IMPOSZIR, PFZSCRA, PFZSCR, SHDGLXBASE, SHDGLXDARK, SHDGLXLGHT, ZIRDSSCR, ZIRHTSCRA, ZIRHTSCR, ZIRMOSCRA, ZIRMOSCR |
| <b>Manufacturer</b>                 | <b>Zima International, Inc. Also DBA Dandy</b><br>22 Cortlandt St, 30th Floor, New York, NY, 10007<br>US                                                                                                                                                |
| <b>UK Responsible Person (UKRP)</b> | <b>Name:</b> DANDY LABS GB, LTD<br><b>Company No.:</b> Company Number: 16873608<br><b>Address:</b> 5 New Street Square, London EC4A 3TW, United Kingdom                                                                                                 |
| <b>Intended User</b>                | Prescriber — Dental professional / dental practice                                                                                                                                                                                                      |
| <b>Device Type</b>                  | Custom-made dental device (CMD)                                                                                                                                                                                                                         |
| <b>Sterility</b>                    | Not supplied sterile (if applicable)                                                                                                                                                                                                                    |
| <b>Single-use / Reusable</b>        | Single use                                                                                                                                                                                                                                              |

**1) Intended Purpose (Intended Use)**

Dandy's ceramic crowns are made to restore worn, carious, broken, and aesthetically displeasing dentition.

**They:**

- Restore single, permanent teeth (crowns) that are fractured, carious, discolored, or severely worn.
- Replace one or more missing permanent teeth by spanning the edentulous space and supporting them with dental implants.
- Restore the form, function, and esthetics of natural dentition.
- Provide occlusal stability and maintain proper interarch relationships.

**If used to restore an implant:**

- Connected to the implant fixture with an abutment
- Cases that are screw-retained are permanently fixed to the abutment in the lab

**2) Device Description / Key Specifications**

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.
- Configuration: Dandy's ceramic crowns are fabricated using digital technology, allowing for a precise fit and customization to the patient's specific needs. They are made of a highly durable and biocompatible material, ensuring long-term comfort and function.
- Material: Zirconia Monolithic (3Y & 4Y), Porcelain-Fused-To-Zirconia (PFZ), E.Max (full contour), E.Max (layered)

[Zirconia Material Safety Data Sheet](#)

[PFZ Material Safety Data Sheet](#)

[Full Contour E.MAX Material Safety Data Sheet](#)

### 3) Contraindications

Crowns should be used with caution in patients with uncontrolled periodontal disease, severe bone resorption, or oral lesions.

### 4) Warnings / Precautions / Potential Risks

#### **WARNING**

The use of fixed restorations may have potential risks, including failure, misfit, breakage, and recurrent caries.

#### **Warning – Sensitivity:**

This device may contain trace levels of metals used as colorants in the ceramic material. The treating clinician should inform patients with sensitivity or allergy of the presence of metals in the device before treatment. In rare cases, allergic reactions can occur in individuals who are sensitized.

#### **Warning – Irritation:**

Localized irritation of surrounding oral tissues may occur in certain individuals due to individual variability in tissue response. Patients should be monitored post-placement, and the treating clinician should evaluate any signs of irritation, inflammation, or discomfort.

#### **CAUTION**

NA

### 5) Cleaning and Care

Daily Cleaning

#### **Oral Hygiene:**

##### **Brushing:**

- Recommended brushing teeth twice a day, using a soft-bristled toothbrush and a fluoride toothpaste.

##### **Flossing:**

- Recommend daily flossing to remove plaque and debris from between the teeth and around the fixed restoration. A floss threader may be needed for hard-to-reach areas (for bridges).

##### **Mouthwash:**

- Alcohol-based mouthwashes tend to dry out the mouth, leading to increased plaque and tartar buildup.

#### Dietary and Habit Recommendations

- Avoid abrasive foods: Limit hard, crunchy, or sticky foods that can damage the restoration or cause it to loosen.
- Avoid sugary foods: where new cavities can form or cause premature decay of the restorative margins.
- Stay hydrated: Drinking plenty of water will help keep your patient's mouth clean and prevent gum disease.
- Avoid bad habits: Refrain from using teeth to open packages or bite on hard objects.

## 6) Storage

Use dedicated dental crown/veneer boxes with soft foam or silicone inserts to prevent chipping, scratching, or movement, ensuring they are kept clean, dry, and protected during transport or while awaiting placement in a clinic or lab setting..

## 7) Expected Life / Service

Ceramic dental restorations (crowns, fillings, veneers) generally last 10 to 20 years, with modern materials like zirconia potentially exceeding this, but lifespans vary significantly based on care, location (back teeth wear faster), and type (implants are more permanent than crowns on top). Excellent oral hygiene (brushing, flossing) and avoiding hard foods or teeth grinding (bruxism) are key to maximizing longevity, though some chipping or wear is possible.

## 8) Incident / Complaint Reporting

#### Report suspected serious incidents associated with this device to:

- The Manufacturer and/or UK Responsible Person (UKRP) using the contacts above, and
- The UK Medicines and Healthcare products Regulatory Agency (MHRA) — [www.gov.uk/mhra](http://www.gov.uk/mhra)

In accordance with The Medical Devices Regulations 2002 (as amended).

#### INFO

*For full contact details of the local Authorized Representative, refer to: Authorized Representatives — end of this document.*

## 9) Disposal

- All-ceramic (porcelain) restorations are often disposed of as regular solid waste or taken to specific recycling centers that accept mixed ceramics/construction debris, as they're challenging to recycle like typical items.
- Household Trash: Place in your regular garbage bin, as they don't break down in standard recycling
- Always check with your local municipality or waste hauler for specific rules; however, ceramic items are generally not accepted in standard curbside recycling bins.
- Recycling Centers: Take them to a local recycling center that accepts mixed ceramics, concrete, or construction debris, as they can sometimes be repurposed.
- If the restoration has been in the patient's mouth or has come into contact with saliva, discard it as biomedical waste.

**EUROPEAN UNION (ENGLISH)****INSTRUCTIONS FOR USE (IFU)**

European Union — Custom-Made Device (CMD)

|                                            |                                                                                                                                                                                                                                                        |
|--------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>GMDN / EMDN</b>                         | 38594 / Q010206- Dental Prostheses                                                                                                                                                                                                                     |
| <b>Product Code(s)</b>                     | EMXLPSCRA, EMXLPSCR, EMXSCRA, EMXSCR, IMPANTEM, IMPANTHTZ, IMPANTLEMM, IMPANTLVD, IMPANTZIR, IMPOSEMM, IMPOSHTZ, IMPOSLEY, IMPOSEMM, IMPOSZIR, PFZSCRA, PFZSCR, SHDGLXBASE, SHDGLXDARK, SHDGLXLGHT, ZIRDSSCR, ZIRHTSCRA, ZIRHTSCR, ZIRMOSCRA, ZIRMOSCR |
| <b>Manufacturer</b>                        | <b>Zima International, Inc. Also DBA Dandy</b><br>22 Cortlandt St, 30th Floor, New York, NY, 10007<br>US                                                                                                                                               |
| <b>EU authorized representative (EURP)</b> | <b>Name:</b> Dandy Labs Europe<br><b>Company No.:</b> Registration No.: 989857065<br><b>Address:</b> 3, boulevard de Sebastopol, 75001 Paris, France                                                                                                   |
| <b>Intended User</b>                       | Prescriber — Dental professional / dental practice                                                                                                                                                                                                     |
| <b>Device Type</b>                         | Custom Made Device (CMD)                                                                                                                                                                                                                               |
| <b>Sterility</b>                           | Not supplied sterile (if applicable)                                                                                                                                                                                                                   |
| <b>Single-use / Reusable</b>               | Single-use                                                                                                                                                                                                                                             |

**1) Intended Purpose (Intended Use)**

Dandy's ceramic crowns are made to restore worn, carious, broken, and aesthetically displeasing dentition.

**They:**

- Restore single, permanent teeth (crowns) that are fractured, carious, discolored, or severely worn.
- Replace one or more missing permanent teeth by spanning the edentulous space and supporting them with dental implants.
- Restore the form, function, and esthetics of natural dentition.
- Provide occlusal stability and maintain proper interarch relationships.

**If used to restore an implant:**

- Connected to the implant fixture with an abutment
- Cases that are screw-retained are permanently fixed to the abutment in the lab

**2) Device Description / Key Specifications**

- Custom-made, patient-specific device manufactured per dental professional prescription and digital input data.

- Configuration: Dandy's ceramic crowns are fabricated using digital technology, allowing for a precise fit and customization to the patient's specific needs. They are made of a highly durable and biocompatible material, ensuring long-term comfort and function.
- Material: Zirconia Monolithic (3Y & 4Y), Porcelain-Fused-To-Zirconia (PFZ), E.Max (full contour), E.Max (layered)

[Zirconia Material Safety Data Sheet](#)

[PFZ Material Safety Data Sheet](#)

[Full Contour E.MAX Material Safety Data Sheet](#)

### 3) Contraindications

Crowns should be used with caution in patients with uncontrolled periodontal disease, severe bone resorption, or oral lesions.

### 4) Warnings / Precautions / Potential Risks

#### **WARNING**

The use of fixed restorations may have potential risks, including failure, misfit, breakage, and recurrent caries.

#### **Warning – Sensitivity:**

This device may contain trace levels of metals used as colorants in the ceramic material. The treating clinician should inform patients with sensitivity or allergy of the presence of metals in the device before treatment. In rare cases, allergic reactions can occur in individuals who are sensitized.

#### **Warning – Irritation:**

Localized irritation of surrounding oral tissues may occur in certain individuals due to individual variability in tissue response. Patients should be monitored post-placement, and the treating clinician should evaluate any signs of irritation, inflammation, or discomfort.

#### **CAUTION**

NA

### 5) Cleaning and Care

Daily Cleaning

#### **Oral Hygiene:**

##### **Brushing:**

- Recommended brushing teeth twice a day, using a soft-bristled toothbrush and a fluoride toothpaste.

##### **Flossing:**

- Recommend daily flossing to remove plaque and debris from between the teeth and around the fixed restoration. A floss threader may be needed for hard-to-reach areas (for bridges).

##### **Mouthwash:**

- Alcohol-based mouthwashes tend to dry out the mouth, leading to increased plaque and tartar buildup.

#### **Dietary and Habit Recommendations**

- Avoid abrasive foods: Limit hard, crunchy, or sticky foods that can damage the restoration or cause it to loosen.
- Avoid sugary foods: where new cavities can form or cause premature decay of the restorative margins.
- Stay hydrated: Drinking plenty of water will help keep your patient's mouth clean and prevent gum disease.
- Avoid bad habits: Refrain from using teeth to open packages or bite on hard objects.

## 6) Storage

Use dedicated dental crown/veneer boxes with soft foam or silicone inserts to prevent chipping, scratching, or movement, ensuring they are kept clean, dry, and protected during transport or while awaiting placement in a clinic or lab setting..

## 7) Expected Life / Service

Ceramic dental restorations (crowns, fillings, veneers) generally last 10 to 20 years, with modern materials like zirconia potentially exceeding this, but lifespans vary significantly based on care, location (back teeth wear faster), and type (implants are more permanent than crowns on top). Excellent oral hygiene (brushing, flossing) and avoiding hard foods or teeth grinding (bruxism) are key to maximizing longevity, though some chipping or wear is possible.

## 8) Incident / Complaint Reporting

### **Report suspected serious incidents associated with this device to:**

- The Manufacturer and/or Dandy EU entity using the contacts above, and corresponding local requirements
- The Competent Authority of your Member State, in accordance with Regulation (EU) 2017/745 (MDR).

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**INFO**

For the full contact information of the local Authorized Representative, see: Authorized Representatives — at the end of the document.

## 9) Disposal

- All-ceramic (porcelain) restorations are often disposed of as regular solid waste or taken to specific recycling centers that accept mixed ceramics/construction debris, as they're challenging to recycle like typical items.
- Household Trash: Place in your regular garbage bin, as they don't break down in standard recycling
- Always check with your local municipality or waste hauler for specific rules; however, ceramic items are generally not accepted in standard curbside recycling bins.
- Recycling Centers: Take them to a local recycling center that accepts mixed ceramics, concrete, or construction debris, as they can sometimes be repurposed.
- If the restoration has been in the patient's mouth or has come into contact with saliva, discard it as biomedical waste.

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The table below lists the local Authorized Representative or Sponsor for each country / jurisdiction where this device is marketed. For incident reporting, contact the representative for your country.

*Le tableau ci-dessous liste le Représentant Autorisé local pour chaque pays. | Die folgende Tabelle enthält die bevollmächtigten Vertreter je Land. | La tabella seguente elenca il rappresentante autorizzato o lo sponsor locale per ciascun paese. | La tabla siguiente enumera el Representante Autorizado o Patrocinador local para cada país. | Følgende tabel viser den autoriserede repræsentant eller lokale sponsor for hvert land. | De volgende tabel geeft een overzicht van de gemachtigde vertegenwoordiger of lokale sponsor voor elk land. | W poniższej tabeli wymieniono upoważnionych przedstawicieli lub lokalnych sponsorów dla każdego kraju. | Följande tabell listar den auktoriserade representanten eller den lokala sponsorn för varje land. | A tabela seguinte lista o Representante Autorizado ou Patrocinador local para cada país. | Tabellen nedenfor viser den autoriserte representanten eller den lokale sponsoren for hvert land.*

| Jurisdiction     | Authorised Rep / Sponsor | Registered Address / Contact Email                                                   | Reg. Number      | Regulatory Body                                     |
|------------------|--------------------------|--------------------------------------------------------------------------------------|------------------|-----------------------------------------------------|
| Australia (AUS)  | DANDY Labs AUS Pty Ltd   | Tower One, Level 46, 100 Barangaroo Ave, Sydney NSW 2000<br>Compliance@meetdandy.com | ABN: 33693955761 | TGA (Therapeutic Goods Administration)              |
| EU - English     | Dandy Labs Europe        | 3, boulevard de Sebastopol, 75001 Paris, France<br>Compliance@meetdandy.com          | Reg: 989857065   | Competent Authority / EUDAMED                       |
| EU - España      | Dandy Labs Europa        | 3, boulevard de Sebastopol, 75001 Paris, Francia<br>Compliance@meetdandy.com         | Reg: 989857065   | AEMPS / Autoridad Competente                        |
| EU - France      | Dandy Labs Europe        | 3, boulevard de Sebastopol, 75001 Paris, France<br>Compliance@meetdandy.com          | Reg: 989857065   | ANSM / Autorité Compétente                          |
| EU - Deutschland | Dandy Labs Europe        | 3, boulevard de Sebastopol, 75001 Paris, Frankreich<br>Compliance@meetdandy.com      | Reg: 989857065   | BfArM / Zuständige Behörde                          |
| EU - Italia      | Dandy Labs Europe        | 3, boulevard de Sebastopol, 75001 Parigi, Francia<br>Compliance@meetdandy.com        | Reg: 989857065   | Min. della Salute / ISS                             |
| EU - Austria     | Dandy Labs Europe        | 3, boulevard de Sebastopol, 75001 Paris, France<br>Compliance@meetdandy.com          | Reg: 989857065   | Bundesamt für Sicherheit im Gesundheitswesen (BASG) |
| EU - Ireland     | Dandy Labs Europe        | 3, boulevard de Sebastopol, 75001 Paris, France<br>Compliance@meetdandy.com          | Reg: 989857065   | Health Products Regulatory Authority (HPRA)         |
| EU - Denmark     | Dandy Labs Europe        | 3, boulevard de Sebastopol, 75001 Paris, France<br>Compliance@meetdandy.com          | Reg: 989857065   | Lægemiddelstyrelsen                                 |

|                      |                       |                                                                             |                     |                                                                   |
|----------------------|-----------------------|-----------------------------------------------------------------------------|---------------------|-------------------------------------------------------------------|
| EU - Nederland       | Dandy Labs Europe     | 3, boulevard de Sebastopol, 75001 Paris, France<br>Compliance@meetdandy.com | Reg: 989857065      | Inspectie Gezondheidszorg en Jeugd (IGJ)                          |
| EU - Polska          | Dandy Labs Europe     | 3, boulevard de Sebastopol, 75001 Paris, France<br>Compliance@meetdandy.com | Reg: 989857065      | Urząd Rejestracji Produktów Leczniczych (UR PL)                   |
| EU - Belgique/België | Dandy Labs Europe     | 3, boulevard de Sebastopol, 75001 Paris, France<br>Compliance@meetdandy.com | Reg: 989857065      | Federal Agency for Medicines and Health Products (FAMHP)          |
| EU - Sverige         | Dandy Labs Europe     | 3, boulevard de Sebastopol, 75001 Paris, France<br>Compliance@meetdandy.com | Reg: 989857065      | Läkemedelsverket                                                  |
| EU - Portugal        | Dandy Labs Europe     | 3, boulevard de Sebastopol, 75001 Paris, France<br>Compliance@meetdandy.com | Reg: 989857065      | Autoridade Nacional do Medicamento e Produtos de Saúde (INFARMED) |
| EU - Norge           | Dandy Labs Europe     | 3, boulevard de Sebastopol, 75001 Paris, France<br>Compliance@meetdandy.com | Reg: 989857065      | Direktoratet for medisinske produkter (DMP)                       |
| United Kingdom (UK)  | DANDY LABS GB, LTD    | 5 New Street Square, London EC4A 3TW<br>Compliance@meetdandy.com            | Co. No.: 16873608   | MHRA                                                              |
| New Zealand (NZ)     | DANDY Labs NZ Limited | Quigg Partners, Level 7, 36 Brandon Street, Wellington, NZ 6011             | NZBN: 9429053758446 | Medsafe                                                           |

**NOTE**

Contact details are subject to change. Always verify current information on the official product labelling or the manufacturer's website before reporting an incident. / Les coordonnées sont sujettes à modification. / Die Kontaktdaten können sich ändern. / I dati di contatto sono soggetti a modifiche. / Los datos de contacto están sujetos a cambios. / Kontaktopplysningerne kan ændres. / De contactgegevens zijn onder voorbehoud van wijzigingen. / Dane kontaktowe mogą ulec zmianie. / Kontaktoppgifterna kan komma att ändras. / Os dados de contacto estão sujeitos a alterações. / Kontaktopplysningene kan endres.

**Manufacturer / Fabricant /  
Hersteller / Produttore /  
Fabricante / Producent /  
Fabrikant / Producent /  
Tillverkare / Fabricante /  
Produsent**

**Zima International, Inc. Also DBA Dandy**  
22 Cortlandt St, 30th Floor, New York, NY, 10007, US  
compliance@meetdandy.com

This is a representation of an electronic record that was signed electronically and this page is the manifestation of the electronic signatures.

Signatory Table

| Action Name | User Name | Title                     | Signature Date    |
|-------------|-----------|---------------------------|-------------------|
| Approve     | Kasha Lim | Supplier Quality Engineer | 24-Sep-2026 14:01 |

\* Datetime fields are displayed according to the system time zone: (GMT-04:00) Eastern Daylight Time (America/New\_York)