



Lucitone Digital Print Denture Implant Overdenture


Implant Overdenture

Dentsply Sirona recommends the following pick-up materials for placing attachments for removable overdentures.

EZ PickUp™ by Sterngold Dental†

Quick Up® by Voco†

Chairside® Attachment Processing Material by Zest†



Dentsply Sirona continues to test and validate materials compatible with Lucitone Digital Print Dentures. Visit the [Lucitone Digital Print Denture Learning Portal](#) to get the most up-to-date material and workflow information.

Dentsply Sirona | 570 West College Avenue | York PA. 17401 | Tel. 800-243-1942 | dentsplysirona.com

DP-00000496 Rev. 0 (12/2019)

*Carbon® is a registered mark owned by Carbon, Inc.

†EZ PickUp™ by Sterngold Dental, Quick Up® by Voco and Chairside® Attachment Processing Material by Zest are not registered trademarks of Dentsply Sirona.